**WICHTIGE KUNDENINFORMATION zu ETHICON Morcellations Geräten
Produktcodes MX0100, MX0200, MX0100R, MX0200R, MD0100, MD0100C,
MD0200, MD0200C, P-MX20A und P-MX20B**

06. Mai, 2014

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Ethicon, Inc.

Referenzen:

- Bulun SE. Uterine fibroids. *N Engl J Med.* 2013 Oct 3;369(14):1344-55.
- Wu JM, Wechter ME, Geller EJ, Nguyen TV, Visco AG. Hysterectomy rates in the United States, 2003. *Obstet Gynecol* 2007;110:1091-5.
- Perri T, Korach J, Sadetzki S, Oberman B, Fridman E, Ben-Baruch G. Uterine leiomyosarcoma: does the primary surgical procedure matter? *Int J Gynecol Cancer.* 2009 Feb;19(2):257-60.
- Park JY, Kim DY, Kim JH, Kim YM, Kim YT, Nam JH. The impact of tumor morcellation during surgery on the outcomes of patients with apparently early low-grade endometrial stromal sarcoma of the uterus. *Ann Surg Oncol.* 2011 Nov;18(12):3453-61.
- Oduyebo T, Rauh-Hain AJ, Meserve EE, Seidman MA, Hinchcliff E, George S, Quade B, Nucci MR, Del Carmen MG, Muto MG. The value of re-exploration in patients with inadvertently morcellated uterine sarcoma. *Gynecol Oncol.* 2014 Feb;132(2):360-5.
- Morice P, Rodriguez A, Rey A, Pautier P, Atallah D, Genestie C, Pomel C, Lhommé C, Haie-Meder C, Du villard P, Castaigne D. Prognostic value of initial surgical procedure for patients with uterine sarcoma: analysis of 123 patients. *Eur J Gynaecol Oncol.* 2003;24(3-4):237-40.
- Kho KA, Nezhat CH. Evaluating the risks of electric uterine morcellation. *JAMA.* 2014 Mar 5;311(9):905-6.
- Leibsohn S, d'Ablaing G, Mishell DR, Schlaerth JB. Leiomyosarcoma in a series of hysterectomies performed for presumed uterine leiomyomas. *Am J Obstet Gynecol.* 1990;162(4):968-974.
- Reiter RC, Wagner PL, Gambone JC. Routine hysterectomy for larger asymptomatic uterine leiomyomata – a reappraisal. *Obstet Gynecol.* 1992;79(4):481-4.
- Parker WH, Fu YS, Berek JS. Uterine sarcoma in patients operated on for presumed leiomyoma and rapidly growing leiomyoma. *Obstet Gynecol.* 1994;83(3):414-8. 4
- Takamizawa S, Minakami H, Usui R, Noguchi S, Ohwada M, Suzuki M, et al. Risk of complications and uterine malignancies in women undergoing hysterectomy for presumed benign leiomyomas. *Gynecol Obstet Invest.* 1999;48(3):193-6.
- Sinha R, Hegde A, Mahajan C, Dubey N, Sundaram M. Laparoscopic myomectomy: do size, number, and location of the myomas form limiting factors for laparoscopic myomectomy? *J Minim Invasive Gynecol.* 2008;15(3):292-300.
- Kamikabeya TS, Etchebehere RM, Nomellini RS, Murta EF. Gynecological malignant neoplasias diagnosed after hysterectomy performed for leiomyoma in a university hospital. *European journal of gynaecological oncology.* 2010;31(6):651-3.
- Rowland M, Lesnock J, Edwards R, Richard S, Zorn K, Sukumvanich P, et al. Occult uterine cancer in patients undergoing laparoscopic hysterectomy with morcellation. *Gynecol Oncol.* 2012;127(1):S29.
- Leung F, Terzibackian JJ. Re “The impact of tumor morcellation during surgery on the prognosis of patients with apparently early uterine leiomyosarcoma”. *Gynecol Oncol.* 2012;124(1):172-173
- Seidman MA, Oduyebo T, Muto MG, Crum CP, Nucci MR, Quade BJ. Peritoneal dissemination complicating morcellation of uterine mesenchymal neoplasms. *PLoS One.* 2012;7(11):e50058.
- Hampton T. Critics of fibroid removal procedure question risks it may pose for women with undetected uterine cancer. *JAMA.* 2014 Mar 5;311(9):891-3.
- Hagemann IS, Hagemann AR, LiVolsi VA, Montone KT, Chu CS. Risk of occult malignancy in morcellated hysterectomy: a case series. *Int J Gynecol Pathol.* 2011 Sep;30(5):476-83.
- Goto A, Takeuchi S, Sugimura K, Maruo T. Usefulness of Gd-DTPA contrast-enhanced dynamic MRI and serum determination of LDH and its isozymes in the differential diagnosis of leiomyosarcoma from degenerated leiomyoma of the uterus. *Int J Gynecol Cancer.* 2002 Jul-Aug;12(4):354-61.