[Firma]  
[Abteilung]  
[Strasse Nummer]  
[PLZ Ort]

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Durchwahl

Datum

QIL 153-010

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11.09.2020

## SICHERHEITSINFORMATION

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Link: <https://www.olympusprofed.com/gi/powerspiral/met>

- Navigation: Gastroenterology > PowerSpiral Enteroscopy > Procedure information > Antegrade Approach

**1-1. Patients with altered GI anatomy**

## Patients with altered GI anatomy

The safety of this endoscope has not been established in patients with altered GI anatomy. Perform endoscopy and endoscopic treatment in these patients only when its potential benefits are greater than its risks.

### Recommendation of PowerSpiral experts

In the absence of studies confirming the safety of PowerSpiral in patients with altered anatomy, when considering its use, one should consider alternative techniques and conclude that the potential benefit of PowerSpiral outweighs the potential risks.

At least 10 cases experiences (per physician) and at least 20 cases experiences (per center) of deep enteroscopy with PowerSpiral is recommended before treating patients with altered GI anatomy.

This website is a mere summary of some important procedure steps and/or product information. It is not meant to replace the instructions for use. Any user of this product must at all times observe all mandatory information for the product, found, in particular, on the labels and the instructions for use.

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Link: <https://www.olympusprofed.com/gi/powerspiral/met>

- Navigation: Gastroenterology > PowerSpiral Enteroscopy > Procedure information > Reported Adverse Events

### 2-2. Reported Adverse Events

|   | Adverse Event   | Location           | Age, gender, and body measurements         | Patient's condition   | Insertion direction | Anesthesia    | Experts' comments  | Relevant information link  |
|---|---|--------------------|--|---|---------------------|---------------|--|--|
| 1 | Bowel perforation   | NA                 | NA   | Crohn's disease (Diagnosed via exam)  | Retrograde          | NA            | · Safety unknown for known stricturing diseases such as Crohn's disease: Perform exam only when benefits outweigh risks.   | <a href="#">2-4</a> Precautions  |
| 2 | · Esophageal tear<br>· Esophageal perforation                   | NA                 | NA   | Dysphagia   | Antegrade           | NA            | · Safety unknown for history of dysphagia and esophageal stricture: Perform exam only when benefits outweigh risks.<br>· All patients with a history of dysphagia should undergo a formal evaluation which includes an EGD before performing PowerSpiral.  | <a href="#">2-4</a> Precautions  |
| 3 | Ileal perforation   | Terminal ileum     | 81-year-old, Female, 158cm, 58kg, BMI:23.2 | · Small intestine bleeding<br>· Diverticulosis  | Retrograde          | Deep sedation | · Continuous rotation in the same location may result in bowel injury with mucosal disruption, bleeding, and deeper injury with perforation, especially within the terminal ileum.<br>· Repositioning the patient and the scope or adding manual abdominal compression if no advancement is achieved. If the rotating power segment is in the same position, this should be a warning to assess continuing efforts or terminating the procedure. | <a href="#">4-1</a> General Techniques<br><a href="#">4-2</a> Insertion Tips   |
| 4 | · Superficial esophageal tear<br>· Small esophageal perforation | Upper esophagus    | 50-year old, Female, 173cm, 51kg, BMI:17.0 | · Tight upper esophageal sphincter<br>· No past medical history of dysphagia                  | Antegrade           | GA            | · Patient's thin body habitus (BMI 17.0) may be problematic for scope intubation.<br>· Confirmation of PowerSpiral insertion accommodation using a bougie (54-60Fr) is highly recommended, especially in early experience (first 10 cases).<br>· If passage into the esophagus remains difficult, repositioning of the patient's position and use of head and neck repositioning by the anesthesiologist or anesthetist may be useful.           | <a href="#">4-2</a> Insertion Tips<br><a href="#">4-3</a> Withdrawal Tips<br><a href="#">4-4</a> Antegrade Approach Note<br><a href="#">5-1</a> Emergency Withdrawal |
| 5 | Severe mucosal injury in proximal esophagus                     | Proximal esophagus | 78-year-old, Female, 160cm, 71kg, BMI:27.7 | · The proximal esophagus may have been less relaxed<br>· No past medical history of dysphagia | Antegrade           | Deep sedation | · Pre-intubation bougie dilation useful if there is any resistance to intubation<br>· General anesthesia may provide more complete muscle relaxation<br>· Neck positioning is important  | <a href="#">4-2</a> Insertion Tips<br><a href="#">4-3</a> Withdrawal Tips<br><a href="#">4-4</a> Antegrade Approach Note   |

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