


Hørsholm, 12 September 2011

URGENT FIELD SAFETY NOTICE CONCERNING THE XO OPERATION LAMP ON XO 4 UNIT

Unfortunately, there have been instances of the reflector becoming damaged and subsequently coming loose from the lamp. To prevent accidents, we have developed a safety foil that we would like you to affix to the back of the reflector as shown in the photo below. See also the attached instructions.



When you have attached the safety foil, please fill in and return the attached "Annex 1" to XO CARE A/S to document that this has been done.

You are welcome to contact your XO dealer / service technician or write to [REDACTED] [REDACTED]@xo-care.com or [REDACTED]@xo-care.com should you have any questions.

The relevant authorities have been notified of this preventive action.

Best regards,

[REDACTED], CEO
XO CARE A/S



Installation instructions Installations of Safety Foil MG-979

The XO 4 must be turned off and the reflector must be cooled down before starting the installation.

1. Clean the central area of the reflectors rear side with alcohol. Wait for about a minute until the alcohol has evaporated.
2. Remove the centre portion of the protective paper.



3. Place the foil under the arc so that it is centred and the middle cut-out clears the lamp's cooling vents.



4. Apply pressure to the foil where the protective paper has been removed.



5. Pull back the one tab to partially remove the protective paper a little over halfway; position the foil under the arc again and press it into place while removing the rest of the protective paper from the opposite side.



6. Do the same on the other side.
7. Finish by thoroughly rubbing the entire adhesive portion, so that the foil sits smoothly over the whole surface.
8. The installation is now complete.



Annex 1

The undersigned hereby confirms mounting of the safety foil on the XO Operation Lamp.

Date: _____

Signature: _____

Filled in by DDS/technician

Serial number (XO 4)	
Manufacturing date	
Distributor	
End-user (name)	
Street address	
City	
Postcode	
E-mail	
Installed the foil on (date)	

Kindly return the completed attachment within five workdays.

You may either fax it or e-mail a scanned copy to:

FAX: +45 70 20 55 10

E-mail: safetyfoil@xo-care.com

XO CARE A/S
Håndværkersvinget 6
2970 Hørsholm

