17. Dezember 2007

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Important Product Information

December 10th, 2007

PRODUCT: BD Plastipak™ Luer Slip Syringe products, sizes 2, 5 and 10ml, supplied in bulk and non sterile

CATALOGUE NUMBERS: #300202, # 301294, and #301295

LOT NUMBERS:

REF	LOT
301294	0705021, 0706016, 0707028
301295	0705028, 0706023
300202	0707029

Dear Mrs. [REDACTED]:

BD has received from you a report of device disconnection involving a BD bulk non-sterile syringe product included in a medical device kit manufactured by your company.

The report is similar to those previously received from some customers in relation to specific lots of BD sterile Plastipak™ Luer Slip Syringes. BD's preliminary investigation concluded that bulk non-sterile syringe products were not affected. However, BD's subsequent investigation concluded that bulk non-sterile syringe products were indeed affected, and action has already been taken to prevent similar effects in current production.

Disconnections that have been reported occurred particularly when syringes were connected to needleless infusion sites that incorporate spring actuated sealing mechanisms. However, disconnections from other devices have also been reported.

We confirm that the dimensional and physical characteristics of the luer tip for these syringes have at all times been compliant with current International Standards.

BD recommends that you provide your customers, as appropriate, with the following information.

To reduce the potential for disconnection and to minimize clinical risk, we advise you to **firmly** attach all luer connections to these syringes, and manually secure syringe/luer fitting connections during infusion procedures, particularly when needleless connection sites are being used. In addition, do not leave these connections unattended until the injection/infusion is complete. If a syringe is used in a critical application, when a secure connection is required, we advise you to use luer lock syringes and luer lock connections only.

Please do not hesitate to contact your contact manager, [REDACTED] (0473/55.43.43), if you have any further questions in this regard.

Please accept our sincere apologies for the inconvenience this issue may have caused. We will be in contact with you in the next few days to arrange a meeting to discuss the issue in greater detail.

Sincerely,

[REDACTED]

Platform Leader – Core Injection