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10F656
10G175
11A900
11A901
11A902
11B107
11B262
11D054
11D539
11F040
11F411

11F412
12A704
12C239
12C795
12D697
13C353
13C354
13C355

# Humeral Components & Glenoid Sphere Impaction Techniques

Aequalis™ Reversed II, Aequalis™ Reversed Fracture  
Surgical Technique Addendum

## Surgical technique

## Humeral Components & Glenoid Sphere Impaction Techniques

Aequalis™ Reversed II, Aequalis™ Reversed Fracture  
Surgical Technique Addendum

# table of contents

## Humeral Components & Glenoid Sphere Impaction Techniques

Aequalis™ Reversed II, Aequalis™ Reversed Fracture & Surgical Technique Addendum

- Humeral Poly Impaction Technique  
p. 4
- Glenoid Sphere Impaction Technique  
p. 7
- Aequalis™ Reversed II Spacer Impaction Technique  
p. 8
- Aequalis™ Hemi-Adaptor Impaction Technique  
p. 9

## —● Humeral Poly Impaction Technique

Following implantation of the final humeral prosthesis, the final poly insert is placed into the stem. An assembly composed of the drill guide MWD074, drill guide handles MWB260, and tilted sphere trials (MWD191 for Dia. 36mm or MWD194 for Dia. 42mm) should be used for impaction of the final poly into the final humeral prosthesis. These items are located in the Aequalis Reversed 2 glenoid tray (YKAD98).

FIGURE 1: Drill Guide, Central Guide Handle, and Sphere Trials

# Surgical Technique

First, the 29mm diameter drill guide (MWD074) is assembled to the central drill guide handle (MWB260). (Figure 2)

Next, the tilted sphere trial is chosen that matches the diameter of the final poly implant to be impacted into the humeral prosthesis (The centered sphere trial may also be chosen). For impaction of a 36mm poly, use the 36mm diameter titled sphere trial for the 29mm diameter baseplate (MWD191). For impaction of a 42mm poly, use the 42mm diameter titled sphere trial for the 29mm diameter baseplate (MWD194).

The trial is then assembled to the drill guide with the pilot of the sphere trial positioned inside the drill guide assembly. (Figure 3)

**NOTE:** This technique should use 29mm diameter instruments only. Do not use the 25mm drill guide or 25mm sphere trials.

**NOTE:** This technique should use tilted, or centered, sphere trials only. Do not use the eccentric sphere trials.

**NOTE:** The connection between the sphere trial and the drill guide assembly is a loose connection. Care should be taken to ensure the parts do not dissociate.

Figure 2: Assembly of drill guide to Handle.

Figure 3: Assembly of trial sphere to Drill Guide assembly.

# Surgical Technique

The metaphyseal component should be thoroughly cleaned and dried.  
If a lateralized spacer is used, impact the spacer according to this addendum (see page 8) into the metaphyseal cup.  
The selected polyethylene insert is then positioned by aligning the insert, orientation notch with the metaphyseal tab.

The impaction assembly (Figure 4) is then positioned inside the final poly implant in the axis of the poly insert. Using a mallet, impact the end of the drill guide handle until the poly implant is fully seating into the metaphyseal component to achieve final fixation (Figure 5).

**NOTE:** Technique for impaction of the Aequalis Reversed Fracture and Adjustable Reversed poly is identical to the Aequalis Reversed technique shown.

Figure 4: 36mm and 42mm diameter impaction assemblies.

Figure 5: Impaction of the poly into the metaphyseal component.

## —● Glenoid Sphere Impaction Technique

Following implantation of the baseplate components, the final glenoid sphere is positioned onto the baseplate per the original surgical technique recommendations. Then, the following handle options are available for impaction of the sphere onto the baseplate (Figure 6):

- A. Handle (9000259) and Impaction Tip (9000359) from the Ascend™ trays (YKAD103)
- B. Handle (9722894) and Impaction Tip (9722902) from the Simpliciti™ trays (YKAD202)
- C. Impaction handle (MWA108 or MWA109) from:
  - The Aequalis™ Cemented trays (YKAD22)
  - The Aequalis™ Fracture trays (YKAD52)
  - The Aequalis™ Press-Fit trays (YKAD69)
- D. Impaction handle (MWB208) from the Aequalis™ Reversed trays (YKAD39/YKAD59)

Figure 6: (A) Simpliciti™ impaction assembly, (B) Ascend™ Impaction assembly, (C) Aequalis™ Cemented/Press-Fit/Fracture impactor and (D) Aequalis™ Reversed Impactor.

## Spacer Impaction Technique

The metaphyseal component should be thoroughly cleaned and dried.

To impact the spacer into the metaphyseal cup, the following impactor options can be used:

- A. Handle (9000259) and Impaction Tip (9000359) from the Ascend™ trays (YKAD103)
- B. Handle (9722894) and Impaction Tip (9722902) from the Simpliciti™ trays (YKAD202)
- C. Impaction handle (MWA108 or MWA109) from
  - The Aequalis Cemented trays (YKAD22)
  - The Aequalis Fracture trays (YKAD52)
  - The Aequalis Press-Fit trays (YKAD69)
  -
- D. Impaction handle (MWB208) from the Aequalis Reversed trays (YKAD39/YKAD59)

After impaction, the central screw is inserted and fully tightened with the 4.5mm screwdriver to secure the spacer onto the metaphysis.


Figure 8: (A) Spacer positioned onto metaphysis (B) Spacer impacted onto metaphysis, and (C) Spacer screw tightened onto spacer.

## Hemi-Prosthesis Adaptor Impaction Technique

Per the recommendation of the original surgical technique, fix the adaptor/metaphysis Union Screw (DWD054 or DWD990). Then, the internal cup of the metaphyseal component is thoroughly cleaned and dried.

For the Aequalis Reversed, the hemi-prosthesis adaptor size must be the same size as the metaphyseal component. For the Aequalis Reversed Fracture, the hemi-prosthesis adaptor is available in only one unique size.

The adaptor is then positioned over the union screw, into the metaphyseal component, and engaged by hand pressure into the metaphysis. Use a small bone tamp to fully seat the hemi-adaptor. The small bone tamp should be positioned around the union screw in several different places.

After impaction, hand check the adaptor to ensure it is well fixed into the metaphysis. A small gap will remain between the adaptor and metaphysis.


Figure 9: (A) Union screw fix onto metaphysis (B) Hemi-adaptor well engaged onto metaphysis by hand pressure and (C) Hemi-adaptor impacted with a small bone tamp.

**US HEADQUARTERS**

Tornier, Inc.  
10801 Nesbitt Avenue  
South  
Bloomington, MN  
55437  
USA  
+ 1 952 426 7600

**TORNIER INTERNATIONAL  
MANUFACTURER**

Tornier S.A.S  
161 rue Lavoisier, Montbonnot  
38330 Saint-Ismier Cedex  
France  
+ 33 (0)4 76 61 35 00  
+ 33 (0)4 76 61 35 33

[www.tornier.com](http://www.tornier.com)