P M H
P R O D U T O S
M E D I C O
H O S P I T A L A R E S

Sociedade Anónima

Capital Social: €6.924.000,00

NIPC: PT 502 376 899

Reg. na C.R.C de Benavente

nº 502 376 899

FSN Ref: 04/2023

FSCA Ref: FSCA03/2023

Date: 2023/10/27

Urgent Field Safety Notice Dringende Sicherheitsinformation

**PMH Extention set with 3 stopcock manifold (HAHNBANK 3-FACH BLAU
VERLÄNGERUNG 180cm)**
Ref.ª AMT 0043 - LOT 31 23 04G

For Attention of: **all costumers**

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nº 502 376 899

Probability of problem arising

No serious incidents have been reported as this problem can be detected during purge. However, as the product is intended for the administration of medication, and given the risk of the leak occurring during the administration of the chemical/drug, there is a high probability that the medication will not be administered with consequences for the patient.

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Background on Issue Hintergurd des Falles

No serious incidents have been reported.

The problem was detected by one customer during the system purging process and the device was discarded. There was no serious incident involving risk to patient or health care professional as the problem was detected before the drug was administered.

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- Review this Field Safety Notice in its entirety and ensure that all users of the above mentioned product in your organization and other concerned persons are informed about this Field Safety Notice.
 - Identify, quarantine, and return affected articles.
 - not use affected devices anymore.
 - Confirm receipt of this information by completing the attached form and return to PMH, SA
 - As distributor, please forward this correction notification to your customers, and collect the batch from customers
 - Return the devices collected from the batch to PMH

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Contact details of local representative

AMT Medica GmbH
Robert-Bosch-Str. 3
71088 Holzgerlingen
Germany
Email: rajkumar.daniel@amt-medica.de
Tel. +49 7031 209 412 27

Enclosure: Customer response form

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P M H
P R O D U T O S
M É D I C O
H O S P I T A L A R E S

Sociedade Anónima
Capital Social: €6.924.000,00
NIPC: PT 502 376 899
Reg. na C.R.C de Benavente
nº 502 376 899

Field Safety Notice distributor Reply Form Feldsicherheitshinweis Verteiler Antwortformular

1. Field Safety Notice (FSN) information Informationen zu Field Safety Notice (FSN)	
FSN Reference number* Referenznummer*	FSN 04/2023
FSN Date*	27/10/2023
Product/ Device name*	Extention Set
Product Code(s)	AMT 0043
Batch/Serial Number (s)	31 23 04G

2. Distributor/Importer Details	
Company Name*	AMT Medica GmbH
Address*	Robert-Bosch-Str. 3 71088 Holzgerlingen Germany
Shipping address	Robert-Bosch-Str. 3 71088 Holzgerlingen Germany
Contact Name*	Rajkumar.Daniel
Title or Function	Quality Control In-Charge
Telephone number*	rajkumar.daniel@amt-medica.de
Email*	Tel. +49 7031 209 412 27

3. Distributors/Importers (Tick all that apply)		
<input type="checkbox"/>	*I confirm the receipt, the reading and understanding of the Field Safety Notice.	Distributor/Importer to complete or enter N/A
<input type="checkbox"/>	I have checked my stock and quarantined inventory	Distributor/Importer to enter quantity and date
<input type="checkbox"/>	I have identified customers that received or may have received this device	
<input type="checkbox"/>	I have attached customer list	
<input type="checkbox"/>	I have informed the identified customers of this FSN	Date of communication:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers	
<input type="checkbox"/>	I have returned affected devices - enter number of devices returned and date complete.	Add quantity, Lot/Serial Number/Date Returned

P M H

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Sociedade Anónima

Capital Social: €6.924.000,00

NIPC: PT 502 376 899

Reg. na C.R.C de Beja

nº 502 376 899

<input type="checkbox"/>	I have destroyed affected devices – enter number destroyed and date complete.	Add quantity, Lot/Serial Number/Date Returned
<input type="checkbox"/>	Neither I nor any of my customers has any affected devices in inventory	
Print Name*		Distributor/Importer print name here
Signature*		Distributor/Importer sign Here
Date *		