

URGENT FIELD SAFETY NOTICE

Software Anomaly Using

VITROS[®] 5600 Integrated System Software Version 1.6 & Below

Dear Customer,

As part of a Field Safety Corrective Action, the purpose of this communication is to inform you that Ortho Clinical Diagnostics (OCD) has identified an anomaly with VITROS[®] 5600 Integrated System Software Version 1.6 and below. Internal testing determined that under very *specific* conditions (listed below), a surfactant addition (i.e., VITROS[®] Chemistry Products FS Diluent Pack 4) may inadvertently be added to a CuveTip or an empty CuveTip position for which it was not intended.

Investigation Summary & Resolution

Our internal investigation determined that while performing Drugs of Abuse Tests (DATs) on the VITROS[®] 5600 System, all VITROS[®] MicroTip assays could potentially be affected by the anomaly. In addition, samples that are scheduled for on-board dilutions for VITROS[®] MicroSlides could also be affected.

For the anomaly to occur on a VITROS[®] 5600 System, all of the following must occur in order:

1. Surfactant from Diluent Pack 4 for a DAT assay is not added to the CuveTip within the maximal allowable time limit, the test is cancelled (i.e., *No Result*) and the CuveTip is discarded. A condition code (i.e., PEW-026) is generated. **Note:** The rate of this occurring is infrequent; however, the system responds as intended.
2. Due to the software anomaly, the pending surfactant addition is not cancelled as expected.
3. The next available CuveTip Ring position is the same one that was previously cancelled. The system places a new CuveTip from a different sample into the same CuveTip position as the cancelled position.
4. Surfactant is then inadvertently added to the new (unintended) CuveTip.

This anomaly was identified during internal testing. Subsequently, we reviewed the last four months of e-Connectivity[™] data and found two occurrences out of 1.2 million DATs processed. **Refer to the enclosure that contains the list of potentially affected assays and the average biases that could be seen if this anomaly occurs on your VITROS[®] 5600 System.**

The resolution to this anomaly is contained in Software Version 2.0 that is scheduled to be released within a few weeks. Until Software Version 2.0 is installed on your VITROS[®] 5600 System, we recommend that you group all DAT assays into a single run as described below.

Actions Required

To help reduce the potential of the anomaly from occurring, OCD recommends:

- Allow all testing to complete prior to processing DAT assays.
- Process DAT assays in groups of no more than 8 patient samples per run.
NOTE: Do not load additional reagents or slides while DAT assays are processing.
- Allow the testing to complete prior to starting the next run of DAT assays.
- Once DAT processing is complete, it is acceptable to process other types of assays.
- Complete and return the Confirmation of Receipt form no later than **May x, 2012.**
- **Refer to the enclosure to determine the potential affect of this anomaly on your test results.**
- Post this notification by each VITROS[®] System in your facility or with the user documentation.

We apologize for the inconvenience this may cause your laboratory. We have anticipated some questions you may have in the following Questions and Answers section. If you have any additional questions, please contact Customer Technical Services at *insert appropriate number*.

Sincerely,

insert appropriate name

insert appropriate title

Enclosure:

Potential Impact to Results on VITROS® 5600 Systems

Questions and Answers

1. How often did this anomaly occur?

This anomaly was identified during internal testing. As a result, we reviewed the last four months of customer data using e-Connectivity™. We found two occurrences out of 1.2 million DATs processed whereby a surfactant addition was inadvertently added to a CuveTip.

We are proactively monitoring systems that are e-connected. OCD will notify your laboratory if surfactant is inadvertently added to an unintended CuveTip on your e-connected VITROS® 5600 System.

2. What VITROS® Systems are affected by this anomaly?

This anomaly only affects VITROS® 5600 Systems with Software Version 1.6 & below that are used to process Drugs of Abuse Tests (DATs). VITROS® 5,1 FS Systems and VITROS® 4600 Systems are not affected.

3. What are the Drugs of Abuse Tests (DATs)?

The following DATs are the only assays that require surfactant from VITROS® FS Diluent Pack 4:

| Product Code | Product Name |
|--------------|------------------------------------|
| 6801988 | VITROS® BARBITUATES Reagent |
| 6801989 | VITROS® BENZODIAZEPINES Reagent |
| 6801991 | VITROS® AMPHETAMINES Reagent |
| 6801994 | VITROS® THC Reagent |
| 6801995 | VITROS® COCAINE METABOLITE Reagent |
| 6801996 | VITROS® METHADONE Reagent |
| 6801997 | VITROS® OPIATES Reagent |
| 6801998 | VITROS® PHENCYCLIDINE Reagent |

4. What are the assays that could potentially be affected?

Our internal investigation determined that while processing DATs, all VITROS® MicroTip assays including User Defined Assays (UDAs) and Manufacturer Validated Applications (MVAs) could potentially be affected by this anomaly. In addition, samples that require on-board dilutions (out of range or standard dilutions) for VITROS® MicroSlides could also be affected. **Refer to enclosure that contains list of potentially affected assays and the average biases that could be seen if this anomaly occurs on your VITROS® 5600 System.**

VITROS® MicroWell assays are not affected.

5. Should I take any action on previously reported patient results using my VITROS® 5600 System?

We recommend that you discuss any concerns you may have regarding previously reported results with your Laboratory Medical Director to determine the appropriate course of action. Results should be reviewed and interpreted in the context of the overall clinical picture.

6. If this anomaly occurs, how will my test results be affected?

An erroneous result could potentially be generated by the VITROS[®] 5600 System under a very specific set of circumstances (defined on page one).

7. How can I determine if this anomaly has occurred on my VITROS[®] 5600 System?

If your VITROS[®] 5600 System is e-connected, OCD is monitoring your system. We have reviewed the last 4 months of data and we would have contacted you if surfactant was inadvertently added to a CuveTip. We will continue to do so until Software Version 2.0 is installed.

For systems that are not e-connected, upon request, OCD can perform a review of your electronic Datalogger files to understand the impact to your laboratory. Please indicate your preference for assistance on your Confirmation of Receipt form. A follow up notification will be sent to you with detailed information.

NOTE: To help reduce the potential of the anomaly from occurring, we recommend that you group all DAT assays into a single run as described on page one until Software Version 2.0 is installed on your VITROS[®] 5600 System.

8. What actions are required if my VITROS[®] 5600 System is configured with an enGen[™] Laboratory Automation System?

For VITROS[®] 5600 System configured with an enGen[™] Laboratory Automation System:

- Do not use the enGen[™] Laboratory Automation System to process samples for DAT assays.
- Prior to performing DATs, stop both Bypass Modules for the VITROS[®] 5600 System you wish to use.
- Load samples for DAT assays at the front of the VITROS[®] 5600 System, group all DAT assays as described on page one.

9. When will this anomaly be resolved?

This anomaly is resolved in Software Version 2.0 that is scheduled to be released within a few weeks.

Confirmation of Receipt – Response Required

URGENT FIELD SAFETY NOTICE

Software Anomaly Using VITROS[®] 5600 Integrated System Software Version 1.6 & Below

So that we can complete our records, please return this form no later than **May xx, 2012**.

FAX TO: *insert appropriate name*

FAX: *insert appropriate number*

Section I – Acknowledgement & Product Assessment

I received the Urgent Field Safety Notice (Ref. CL12-133_EU) regarding an anomaly using VITROS[®] 5600 Integrated System Software Version 1.6 and below. I understand that until Software Version 2.0 is installed on my system, OCD recommends that I group all Drugs of Abuse (DATs) tests into a single run as described on page one.

For VITROS[®] 5600 Systems that are e-connected, OCD is proactively monitoring your system.

For systems that are NOT e-connected, OCD can perform a review of your data upon request.

Please choose one of the following options:

I would like OCD to review my electronic Datalogger files. I understand that I will be contacted by OCD to initiate the process of data retrieval.

or

I do NOT wish to have OCD perform a review of my electronic Datalogger files.

Section II– Confirmation

Your Name: _____ Name of Facility: _____

*Your Signature/Date: _____ Telephone: _____

****Your signature provides confirmation that you have received and understood this notification***

Verify your name and mailing address or add it to Section III if none is present:

Section III – Please Complete If Your Contact Information Has Changed

Customer Number: _____

Customer Name: _____

Address: _____

Address: _____

City: _____

State/Province/Postal: _____

**Potential Impact to results if Anomaly Occurs Using
VITROS® 5600 Integrated System Software Version 1.6 & Below**

| VITROS® MicroSlide Assays Out-of-Range (OOR) on-board dilutions using ADD Default Dilution factors values | | | | |
|--|------------------------|---------------------------|----------------|-----------------------|
| Assay | Body Fluid Type | Test Concentration | Units | %Bias Observed |
| ACET | Serum/plasma | 167 | µg/mL | -16 |
| ALC | Serum/plasma | 221 | mg/dL | -17 |
| ALKP | Serum/plasma | 489 | U/L | -23 |
| ALT | Serum/plasma | 207 | U/L | -12 |
| AMON | Plasma | 168 | µmol/L | -28 |
| AMYL | Urine | 105 | U/L | -69 |
| | Serum/plasma | 424 | U/L | -25 |
| AST | Serum/plasma | 189 | U/L | -24 |
| Bc | Serum/plasma | 4.8 | mg/dL | -22 |
| Bu | Serum/plasma | 9.8 | mg/dL | -22 |
| BUN/UREA | Serum/plasma | 51 | mg/dL (urea N) | -20 |
| Ca | Urine | 10.5 | mg/dL | -26 |
| | Serum/plasma | 10.8 | mg/dL | -24 |
| CHE | Serum/plasma | 8.77 | U/mL | -19 |
| CHOL | Serum/plasma | 248 | mg/dL | -18 |
| CK | Serum/plasma | 1043 | U/L | -24 |
| CKMB | Serum | 47 | U/L | -20 |
| CREA | Serum/plasma | 5.30 | mg/dL | -23 |
| CRP | Serum/plasma | 87 | mg/L | -34 |
| DGXN | Serum/plasma | 2.4 | ng/mL | +5 |
| dHDL | Serum/plasma | 61.3 | mg/dL | -16 |
| ECO2 | Serum/plasma | 29 | mmol/L | -38 |
| Fe | Serum/plasma | 228 | µg/dL | -32 |
| GGT | Serum/plasma | 471 | U/L | -16 |
| GLU | Serum/plasma | 281 | mg/dL | -25 |
| | Urine | 301 | mg/dL | -22 |
| LAC | Plasma | 3.2 | mmol/L | -26 |
| LDH | Serum/plasma | 1804 | U/L | -17 |
| Li | Serum/plasma | 2.4 | mmol/L | -18 |
| LIPA | Serum/plasma | 637 | U/L | -28 |
| Mg | Serum/plasma | 4.4 | mg/dL | -22 |
| PHBR | Serum/plasma | 52.3 | µg/mL | -21 |
| PHOS | Serum/plasma | 6.9 | mg/dL | -18 |
| PHYT | Serum/plasma | 26.9 | µg/mL | -28 |
| PROT | Cerebrospinal Fluid | 167 | mg/dL | -19 |
| SALI | Serum/plasma | 30 | mg/dL | -23 |
| TBIL | Serum/plasma | 13.4 | mg/dL | -25 |

**Potential Impact to results if Anomaly Occurs Using
VITROS® 5600 Integrated System Software Version 1.6 & Below**

| VITROS® MicroSlide Assays (Continued) | | | | |
|--|------------------------|---------------------------|--------------|-----------------------|
| Out-of-Range (OOR) On-Board Dilutions using ADD Default Dilution Factors Values | | | | |
| Assay Name | Body Fluid Type | Test Concentration | Units | %Bias Observed |
| THEO | Serum/plasma | 23.9 | µg/mL | -13 |
| TP | Serum/plasma | 7.1 | g/dL | -17 |
| TRIG | Serum/plasma | 255 | mg/dL | -24 |
| UPRO | Urine | 126 | mg/dL | -12 |
| URIC | Serum/plasma | 10.6 | mg/dL | -25 |

| VITROS® MicroSlide Assays | | | | |
|--|------------------------|---------------------------|--------------|-----------------------|
| Standard On-Board Dilutions using ADD Default Dilution Factors Values | | | | |
| Assay Name | Body Fluid Type | Test Concentration | Units | %Bias Observed |
| CREA | Urine | 129 | mg/dL | -28 |
| BUN/UREA | Urine | 712 | mg/dL urea N | -30 |
| K+ | Urine | 9.8 | mmol/L | -25 |
| | | 50.4 | mmol/L | -28 |
| Mg | Urine | 12.1 | mg/dL | -29 |
| Na+ | Urine | 29 | mmol/L | -1 |
| | | 52 | mmol/L | +1 |
| PHOS | Urine | 52.9 | mg/dL | -30 |
| URIC | Urine | 23.3 | mg/dL | -34 |

| VITROS® MicroTip | | | | |
|---|------------------------|---------------------------|--------------|-----------------------|
| Standard Processing - NO Dilutions | | | | |
| Assay Name | Body Fluid Type | Test Concentration | Units | %Bias Observed |
| AAT | Serum | 274 | mg/dL | -27 |
| AMPH | Urine | 369 | ng/mL | -9 |
| | | 613 | ng/mL | -11 |
| ApoA1 | Serum/Plasma | 114 | mg/dL | -29 |
| ApoB | Serum/Plasma | 121 | mg/dL | -32 |
| ASO | Serum/Plasma | 282 | IU/mL | -61 |
| BARB | Urine | 145 | ng/mL | -10 |
| | | 259 | ng/mL | -10 |
| BENZ | Urine | 130 | ng/mL | +1 |
| | | 288 | ng/mL | -9 |
| C3 | Serum/Plasma | 264 | mg/dL | -22 |
| C4 | Serum/Plasma | 46 | mg/dL | -26 |
| CAFFN | Serum/Plasma | 20 | µg/mL | -16 |

URGENT PRODUCT CORRECTION NOTIFICATION

Software Anomaly Using VITROS® 5600 Integrated System Software Version 1.6 & Below

**Potential Impact to results if Anomaly Occurs Using
VITROS® 5600 Integrated System Software Version 1.6 & Below**

| VITROS® MicroTip (Continued) Standard Processing - NO Dilutions | | | | |
|--|------------------------|---------------------------|-----------------|-----------------------|
| Assay Name | Body Fluid Type | Test Concentration | Units | %Bias Observed |
| COCM | Urine | 93 | ng/mL | -14 |
| | | 169 | ng/mL | -6 |
| d%A1c | Whole Blood | 11.5 | %A1c NGSP Units | -2 |
| dHA1c | Whole Blood | 1.5 | g/dL | -19 |
| dHb | Whole Blood | 14.2 | g/dL | -17 |
| dLDL | Serum/Plasma | 158 | mg/dL | -35 |
| dTIBC | Serum | 335 | µg/dL | -27 |
| GENT | Serum/Plasma | 7.13 | µg/mL | -22 |
| HCY | Serum/Plasma | 38.5 | µmol/L | -32 |
| HPT | Serum/Plasma | 225.27 | mg/dL | -31 |
| hsCRP | Serum/Plasma | 10.8 | mg/L | -26 |
| IgA | Serum/Plasma | 379 | mg/dL | -32 |
| IgG | Serum/Plasma | 1731 | mg/dL | -18 |
| IgM | Serum/Plasma | 208 | mg/dL | -29 |
| METD | Urine | 103 | ng/mL | -6 |
| | | 182 | ng/mL | -12 |
| OP-HI | Urine | 1537 | ng/mL | -10 |
| | | 2527 | ng/mL | -5 |
| OP-LO | Urine | 260 | ng/mL | -9 |
| | | 400 | ng/mL | -10 |
| PALB | Serum/Plasma | 23.7 | mg/dL | -28 |
| PCP | Urine | 17.4 | ng/mL | -6 |
| | | 28.8 | ng/mL | -7 |
| RF | Serum/Plasma | 29 | IU/mL | -22 |
| THC | Urine | 15.9 | ng/mL | -8 |
| | | 28.1 | ng/mL | -12 |
| TOBRA | Serum/Plasma | 7.62 | µg/mL | -31 |
| TRFRN | Serum/Plasma | 530 | mg/dL | -18 |
| VALP | Serum/Plasma | 125.9 | µg/mL | -26 |
| VANC | Serum/Plasma | 35.34 | µg/mL | -27 |

URGENT PRODUCT CORRECTION NOTIFICATION

Software Anomaly Using VITROS® 5600 Integrated System Software Version 1.6 & Below