Urgent Field Safety Notice

Fujirebio Europe N.V.
Technologiepark 6
B-9052 Gent, Belgium

May 30, 2016

Product name: Lumipulse G E2-III Immunoreaction Cartridges Cross Reactivity of Fulvestrant in Estradiol Assay

Dear Customers;

Fujirebio Inc. would like to bring to your immediate attention a situation with the Lumipulse G E2-III Immunoreaction Cartridges (Lumipulse G E2-III) require actions for your laboratory to take.

Description of the problem

We have recognized that there is a cross reactivity of fulvestrant with Lumipulse G E2-III (Estradiol II assay) as below. When the Lumipulse G E2-III assay is used to assess such patients, it might be lead to falsely elevated estradiol test results.

Table 1 – Affected product details

Product	Code number	Lot number
Lumipulse G E2-III Immunoreaction Cartridges	296011	All lots

Risk to Health

As a result of investigation, Fujirebio Inc. determined that there is a cross reactivity of fulvestrant with Lumipulse G E2-III; see table 2 below for details. Falsely elevated Estradiol levels could lead to misinterpreting the patient as pre-menopausal.

Table 2: Cross reaction and results of Lumipulse G E2-III caused by fulvestrant:

Product	Reference sample result [pg/mL] (a)	Spiked sample result [pg/mL] (b)	(b) - (a) [pg/mL]	Percent cross reaction [%]
Lumipulse G E2-III	77.0	140.0	63.0	0.25
	546.0	595.5	49.5	0.20

Reference sample was spiked with Lumipulse G Specimen Diluent 2

Sample was spiked with 25000 pg/mL fulvestrant.

Percent of cross reaction: (spiked sample-reference sample)/concentration of fulvestrant (25000 pg/mL)

Actions to be taken by the Customer

1. Lumipulse G E2-III should not be used for patients being treated with fulvestrant.
2. Lumipulse G E2-III may continue to be used to report results for patients not on fulvestrant therapy.

Transmission of this Field Safety notice

Please bring this notice to the attention of all those who need to be aware within your laboratory or to other who need to be aware of this safety issue. If you have forwarded the affected product to another laboratory, please provide a copy of this letter to them.

We sincerely apologies for any inconvenience this may have caused you. If you or any of the health care providers you serve have any questions regarding this information, please contact our Customer Support department on:



TEL: +32 9 329 1611

e-mail: *Customer.Support@fujirebio-europe.com*

Additional Information

The Lumipulse G E2-III is based on the principle of competitive binding between E2 in the test specimen and E2-ALP conjugate for an anti-monoclonal Estradiol.

Please complete and return the enclosed Response Form as below to [\[add local contact information: fax and email\]](#) by the 15th of June 2016 at the latest, so we are assured you have received this important communication.

//

I confirmed receipt of FIELD SAFETY NOTICE dated May 30, 2016.

(Confirmation)

Your Name: _____

Name of Facility: _____

Your Signature/Date: _____

Telephone: _____

*Your signature provides confirmation that you have received and understood this notification.

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