



**Abbott**

**CORE DIAGNOSTICS**  
 Abbott Laboratories  
 Diagnostics Division  
 Abbott Park, IL, 60064 USA

Single Registration Number (SRN):  
 IE-MF-000001822

## Urgent Product Correction

Immediate Action Required

Date Issued August xx, 2023

**Product**

Product Description	List Number	Lot Number	UDI
Multichem S Plus Level 1	08P88-10	0108032101	(01)05391523440776 (17)230930 (10)0108032101
Multichem S Plus Level 2	08P88-11	0108032102	(01)05391523440783(17)230930 (10)0108032102
Multichem S Plus Level 3	08P88-12	0108032103	(01)05391523440790 (17) 230930(10)0108032103
Multichem S Plus Level 1	08P88-10	0110072101	(01)05391523440776(17)231231 (10)0110072101
Multichem S Plus Level 2	08P88-11	0110072102	(01)05391523440783(17)231231 (10)0110072102
Multichem S Plus Level 3	08P88-12	0110072103	(01)05391523440790(17)231231 (10)0110072103
Multichem S Plus Level 1	08P88-10	0119112101	(01)05391523440776(17)240531 (10)0119112101
Multichem S Plus Level 2	08P88-11	0119112102	(01)05391523440783(17)240531 (10)0119112102
Multichem S Plus Level 3	08P88-12	0119112103	(01)05391523440790(17)240531 (10)0119112103
Multichem S Plus Level 1	08P88-10	0126052201	(01)05391523440776(17)241031 (10)0126052201
Multichem S Plus Level 2	08P88-11	0126052202	(01)05391523440783(17)241031 (10)0126052202
Multichem S Plus Level 3	08P88-12	0126052203	(01)05391523440790(17)241031 (10)0126052203
Multichem S Plus Level 1	08P88-10	0130112201	(01)05391523440776(17)250531 (10)0130112201
Multichem S Plus Level 2	08P88-11	0130112202	(01)05391523440783(17)250531 (10)0130112202
Multichem S Plus Level 3	08P88-12	0130112203	(01)05391523440790(17)250531 (10)0130112203
Multichem S Plus Level 1	08P88-10	0131112201	(01)05391523440776(17)250531 (10)0131112201
Multichem S Plus Level 2	08P88-11	0131112202	(01)05391523440783(17)250531 (10)0131112202

Multichem S Plus Level 3	08P88-12	0131112203	(01)05391523440790(17)250531 (10)0131112203
Multichem IA Plus	08P86-10	032004230	(01)05391523440745(17)250731 (10)032004230
Multichem IA Plus	08P86-10	035106210	(01)05391523440745(17)230930 (10)035106210
Multichem IA Plus	08P86-10	035810210	(01)05391523440745(17)240131 (10)035810210
Multichem IA Plus	08P86-10	036404220	(01)05391523440745 (17)240731 (10)036404220
Multichem IA Plus	08P86-10	037209220	(01)05391523440745 (17)241231 (10)037209220
Multichem IA Plus	08P86-19	035506210	(01)05391523441636 (17)240531 (10)035506210
Multichem IA Plus	08P86-19	036010210	(01)05391523441636 (17)240930 (10)036010210
Multichem IA Plus	08P86-19	036604220	(01)05391523441636 (17)250331 (10)036604220
Multichem IA Plus	08P86-19	037409220	(01)0539152344163(17)250831 (10)037409220
Multichem P	08P90-10	04530621P	(01)05391523440769 (17)231130 (10)04530621P
Multichem P	08P90-10	04561021P	(01)05391523440769 (17)240430 (10)04561021P
Multichem P	08P90-10	04610523P	(01)05391523440769 (17)251130 (10)04610523P

**Explanation**

The purpose of this letter is to inform you that Technopath Manufacturing Ltd, the manufacturer of Multichem IA Plus controls and Multichem S Plus (Assayed) controls, has identified an issue with the 8mL Alinity glass vial commodity used to fill the products listed above. Some customers have experienced broken vials when receiving or thawing Technopath Multichem IA Plus controls, LN 08P86 or Multichem S Plus (Assayed) controls, LN 08P88.

**Impact on  
 Patient Results/  
 User Safety**

Please refer to the attached Technopath Urgent Field Safety Notice.

**Necessary  
Actions to be  
Taken by  
Customer**

- Please review the attached Technopath Urgent Field Safety Notice and follow the required actions.
  - Complete and return the Customer Reply Form.
  - Contact Abbott Customer Service for assistance if needed.
  - If you have forwarded the product listed above to other laboratories, please inform them of this Product Correction and provide to them a copy of this letter.
  - Please retain this letter for your laboratory records.
- 

**Contact  
Information**

If you or any of the health care providers you serve have questions regarding this information, U.S. Customers please contact Customer Service at 1-877-4ABBOTT (available 24 hours a day, 7 days a week). Customers outside the U.S., please contact your local area Customer Service.

Adverse reactions or quality problems experienced with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program online (<http://www.fda.gov/MedWatch/report.htm>), by phone (1-800-332-1088), or by fax (1-800-FDA0178).

If you have experienced any patient or user injury associated with this Field Action, please immediately report the event to your local area Customer Service.

---

**URGENT FIELD SAFETY NOTICE – FSN-RE23001**  
 For attention of: Healthcare Professionals

**Date Issued:** 31/08/2023

**Type of Action:** Field Safety Corrective Action

**Products:**

Product	Part Code	Lot Number	Expiry Date	UDI
Multichem S Plus Level 1	08P88-10	0108032101	30/09/2023	(01)05391523440776 (17)230930(10)0108032101
Multichem S Plus Level 2	08P88-11	0108032102	30/09/2023	(01)05391523440783 (17)230930(10)0108032102
Multichem S Plus Level 3	08P88-12	0108032103	30/09/2023	(01)05391523440790 (17)230930(10)0108032103
Multichem S Plus Level 1	08P88-10	0110072101	31/12/2023	(01)05391523440776 (17)231231(10)0110072101
Multichem S Plus Level 2	08P88-11	0110072102	31/12/2023	(01)05391523440783 (17)231231(10)0110072102
Multichem S Plus Level 3	08P88-12	0110072103	31/12/2023	(01)05391523440790 (17)231231(10)0110072103
Multichem S Plus Level 1	08P88-10	0119112101	31/05/2024	(01)05391523440776 (17)240531(10)0119112101
Multichem S Plus Level 2	08P88-11	0119112102	31/05/2024	(01)05391523440783 (17)240531(10)0119112102
Multichem S Plus Level 3	08P88-12	0119112103	31/05/2024	(01)05391523440790 (17)240531(10)0119112103
Multichem S Plus Level 1	08P88-10	0126052201	31/10/2024	(01)05391523440776 (17)241031(10)0126052201
Multichem S Plus Level 2	08P88-11	0126052202	31/10/2024	(01)05391523440783 (17)241031(10)0126052202
Multichem S Plus Level 3	08P88-12	0126052203	31/10/2024	(01)05391523440790 (17)241031(10)0126052203
Multichem S Plus Level 1	08P88-10	0130112201	31/05/2025	(01)05391523440776 (17)250531(10)0130112201
Multichem S Plus Level 2	08P88-11	0130112202	31/05/2025	(01)05391523440783 (17)250531(10)0130112202
Multichem S Plus Level 3	08P88-12	0130112203	31/05/2025	(01)05391523440790 (17)250531(10)0130112203



<b>Product</b>	<b>Part Code</b>	<b>Lot Number</b>	<b>Expiry Date</b>	<b>UDI</b>
Multichem S Plus Level 1	08P88-10	0131112201	31/05/2025	(01)05391523440776 (17)250531(10)0131112201
Multichem S Plus Level 2	08P88-11	0131112202	31/05/2025	(01)05391523440783 (17)250531(10)0131112202
Multichem S Plus Level 3	08P88-12	0131112203	31/05/2025	(01)05391523440790 (17)250531(10)0131112203
Multichem IA Plus	08P86-10	032004230	31/07/2025	(01)05391523440745 (17)250731(10)032004230
Multichem IA Plus	08P86-10	035106210	30/09/2023	(01)05391523440745 (17)230930(10)035106210
Multichem IA Plus	08P86-10	035810210	31/01/2024	(01)05391523440745 (17)240131(10)035810210
Multichem IA Plus	08P86-10	036404220	31/07/2024	(01)05391523440745 (17)240731(10)036404220
Multichem IA Plus	08P86-10	037209220	31/12/2024	(01)05391523440745 (17)241231(10)037209220
Multichem IA Plus	08P86-19	035506210	31/05/2024	(01)05391523441636 (17)240531(10)035506210
Multichem IA Plus	08P86-19	036010210	30/09/2024	(01)05391523441636 (17)240930(10)036010210
Multichem IA Plus	08P86-19	036604220	31/03/2025	(01)05391523441636 (17)250331(10)036604220
Multichem IA Plus	08P86-19	037409220	31/08/2025	(01)05391523441636 (17)250831(10)037409220
Multichem P	08P90-10	04530621P	30/11/2023	(01)05391523440769 (17)231130(10)04530621P
Multichem P	08P90-10	04561021P	30/04/2024	(01)05391523440769 (17)240430(10)04561021P
Multichem P	08P90-10	04610523P	30/11/2025	(01)05391523440769 (17)251130(10)04610523P



Technopath is part of LGC Clinical Diagnostics

Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, Ireland. t: +353 61 525 700 | e: [technopath.info@lgcgroup.com](mailto:technopath.info@lgcgroup.com) | w: [www.technopathcd.com](http://www.technopathcd.com)

**Reason for Notice:**

Dear Valued Customer,

Technopath Manufacturing Ltd. is initiating a field action for the products listed above.  
 This letter contains important information that needs your immediate attention.


<p><b>ISSUE:</b></p>	<p>Technopath Manufacturing Ltd. has identified an issue with the glass vial commodity used to fill the products as listed above.          Some customers have experienced broken vials when receiving or thawing Technopath Multichem IA Plus controls, LN 08P86 or Multichem S Plus (Assayed) controls, LN 08P88.</p> <p>Breakage has been observed upon delivery, while thawing, or while opening the vial.</p> <p>The breakage can result in the bottom of the vial breaking off, breakage in the neck thread area, or cracking on the sides of the vial.</p> <p>Technopath, Abbott and the vial supplier evaluated the vial manufacturing process for potential enhancements. Extensive investigation and verification studies were performed, and a strengthened glass vial design has now been qualified for use.</p> <p>Therefore, all future lots will be manufactured at Technopath Manufacturing Ltd. using the strengthened vials and will become available through the course of 2024.</p>
<p><b>IMPACT ON PATIENT RESULTS / USER SAFETY:</b></p>	<p>No patient results are impacted by this issue as there is no change to the product formulation.</p> <p>There is the potential for minor injuries or biohazardous exposure due to the possibility of vials breaking during thawing or handling.</p>
<p><b>NECESSARY ACTION TO BE TAKEN BY CUSTOMER:</b></p>	<p>Continue to adhere to the following points from the Procedures for Use in the respective Instructions for Use:</p> <ul style="list-style-type: none"> <li>➤ Allow the control vial to stand at room temperature (18 to 25°C) for 30 minutes or until completely thawed. The control must be completely solubilized before continuing to the next step.</li> <li>➤ Place the thawed control vial on a suitable mixing device, e.g. roller mixer for 15 minutes, to achieve a homogenous solution.</li> </ul>



	<p>If you identify that a vial has broken or cracks are visible, safely dispose of the vial in line with Safety Data Sheet (SDS).</p> <p>Continue to use caution to prevent operator exposure and wear gloves when handling the control, as it contains material of human and animal origin. In addition, refer to the respective Safety Data Sheet (SDS) for more information regarding exposure controls and personal protection.</p> <p>Please share this information with your laboratory staff who need to be aware of this safety notice and retain this notification as part of your laboratory Quality System documentation. If you have forwarded any of the impacted product(s) listed above to another laboratory, please provide them a copy of this letter.</p>
<b>CONTACT</b>	<p>If you or any of the health care providers you serve have questions regarding this notice, U.S. Customers please contact Customer Service at 1-877-4ABBOTT (available 24 hours a day, 7 days a week). Customers outside the U.S., please contact your local area Customer Service.</p>

We sincerely apologize for any inconvenience this may cause your laboratory and hope for your understanding and cooperation.

Yours Sincerely



Technopath is part of LGC Clinical Diagnostics

Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, Ireland. t: +353 61 525 700 | e: [technopath.info@lgcgroup.com](mailto:technopath.info@lgcgroup.com) | w: [www.technopathcd.com](http://www.technopathcd.com)



**CORE DIAGNOSTICS**  
 Abbott Laboratories  
 Diagnostics Division  
 Abbott Park, IL, 60064 USA

Single Registration Number (SRN):  
 IE-MF-000001822

# Customer Reply

Immediate Action Required

Core Diagnostics Product Correction letter dated **Month/Day/Year** - FADDMMMYYYY

Product	Product Description			List Number			Lot Number			UDI		
	Product Description	List Number	Lot Number	Product Description	List Number	Lot Number	Product Description	List Number	Lot Number	Product Description	List Number	Lot Number
Product	Multichem S Plus Level 1	08P88-10	0108032101	Multichem S Plus Level 1	08P88-10	0110072101	Multichem S Plus Level 1	08P88-10	0126052201	Multichem S Plus Level 1	08P88-10	0130112201
	Multichem S Plus Level 2	08P88-11	0108032102	Multichem S Plus Level 2	08P88-11	0110072102	Multichem S Plus Level 2	08P88-11	0126052202	Multichem S Plus Level 2	08P88-11	0130112202
	Multichem S Plus Level 3	08P88-12	0108032103	Multichem S Plus Level 3	08P88-12	0110072103	Multichem S Plus Level 3	08P88-12	0126052203	Multichem S Plus Level 3	08P88-12	0130112203
	Multichem S Plus Level 1	08P88-10	0119112101	Multichem S Plus Level 1	08P88-10	0119112102	Multichem S Plus Level 1	08P88-10	0131112201	Multichem S Plus Level 1	08P88-10	0131112202
	Multichem S Plus Level 2	08P88-11	0119112102	Multichem S Plus Level 2	08P88-11	0119112103	Multichem S Plus Level 2	08P88-11	0131112203	Multichem S Plus Level 2	08P88-11	0131112204
	Multichem S Plus Level 3	08P88-12	0119112103	Multichem S Plus Level 3	08P88-12	0119112104	Multichem S Plus Level 3	08P88-12	0131112205	Multichem S Plus Level 3	08P88-12	0131112206
	Multichem S Plus Level 1	08P88-10	01240531	Multichem S Plus Level 1	08P88-10	01240532	Multichem S Plus Level 1	08P88-10	01240533	Multichem S Plus Level 1	08P88-10	01240534
	Multichem S Plus Level 2	08P88-11	01240532	Multichem S Plus Level 2	08P88-11	01240533	Multichem S Plus Level 2	08P88-11	01240534	Multichem S Plus Level 2	08P88-11	01240535
	Multichem S Plus Level 3	08P88-12	01240533	Multichem S Plus Level 3	08P88-12	01240534	Multichem S Plus Level 3	08P88-12	01240535	Multichem S Plus Level 3	08P88-12	01240536
	Multichem S Plus Level 1	08P88-10	01241031	Multichem S Plus Level 1	08P88-10	01241032	Multichem S Plus Level 1	08P88-10	01241033	Multichem S Plus Level 1	08P88-10	01241034
	Multichem S Plus Level 2	08P88-11	01241032	Multichem S Plus Level 2	08P88-11	01241033	Multichem S Plus Level 2	08P88-11	01241034	Multichem S Plus Level 2	08P88-11	01241035
	Multichem S Plus Level 3	08P88-12	01241033	Multichem S Plus Level 3	08P88-12	01241034	Multichem S Plus Level 3	08P88-12	01241035	Multichem S Plus Level 3	08P88-12	01241036
	Multichem S Plus Level 1	08P88-10	01250531	Multichem S Plus Level 1	08P88-10	01250532	Multichem S Plus Level 1	08P88-10	01250533	Multichem S Plus Level 1	08P88-10	01250534
	Multichem S Plus Level 2	08P88-11	01250532	Multichem S Plus Level 2	08P88-11	01250533	Multichem S Plus Level 2	08P88-11	01250534	Multichem S Plus Level 2	08P88-11	01250535
	Multichem S Plus Level 3	08P88-12	01250533	Multichem S Plus Level 3	08P88-12	01250534	Multichem S Plus Level 3	08P88-12	01250535	Multichem S Plus Level 3	08P88-12	01250536
	Multichem S Plus Level 1	08P88-10	01250534	Multichem S Plus Level 1	08P88-10	01250535	Multichem S Plus Level 1	08P88-10	01250536	Multichem S Plus Level 1	08P88-10	01250537
	Multichem S Plus Level 2	08P88-11	01250535	Multichem S Plus Level 2	08P88-11	01250536	Multichem S Plus Level 2	08P88-11	01250537	Multichem S Plus Level 2	08P88-11	01250538
	Multichem S Plus Level 3	08P88-12	01250536	Multichem S Plus Level 3	08P88-12	01250537	Multichem S Plus Level 3	08P88-12	01250538	Multichem S Plus Level 3	08P88-12	01250539
	Multichem IA Plus	08P86-10	032004230	Multichem IA Plus	08P86-10	035106210	Multichem IA Plus	08P86-10	035810210	Multichem IA Plus	08P86-10	035810210



	Multichem IA Plus	08P86-10	036404220	(01)05391523440745 (17)240731 (10)036404220
	Multichem IA Plus	08P86-10	037209220	(01)05391523440745 (17)241231 (10)037209220
	Multichem IA Plus	08P86-19	035506210	(01)05391523441636 (17)240531 (10)035506210
	Multichem IA Plus	08P86-19	036010210	(01)05391523441636 (17)240930 (10)036010210
	Multichem IA Plus	08P86-19	036604220	(01)05391523441636 (17)250331 (10)036604220
	Multichem IA Plus	08P86-19	037409220	(01)0539152344163(17)250831 (10)037409220
	Multichem P	08P90-10	04530621P	(01)05391523440769 (17)231130 (10)04530621P
	Multichem P	08P90-10	04561021P	(01)05391523440769 (17)240430 (10)04561021P
	Multichem P	08P90-10	04610523P	(01)05391523440769 (17)251130 (10)04610523P
<b>Instructions</b>	<ol style="list-style-type: none"> <li>Please provide a copy of the accompanying Product <b>Correction</b> letter to the laboratory manager, supervisor or health professional responsible for the impacted product.</li> <li>Please <b>complete</b> all sections and return this Customer Reply Form to the Abbott contact prior to <b>DDMMYYYY</b>. Even if you no longer have the instrument(s)/reagent(s), this form is required for the reconciliation of our records.</li> </ol>			
<b>Abbott contact</b>	<ul style="list-style-type: none"> <li>E-mail: PMS@abbott.com (insert local e-mail address here)</li> <li>Fax: 1-800-777-0051 (insert local number here)</li> </ul>			
<b>Acknowledgement</b>	<p>By completing and signing this document I confirm that the Product Correction Letter was understood and that the necessary actions for the customer were completed. If not, please choose one of the options below.</p> <p> <input type="checkbox"/> No, I would like to be contacted by an Abbott Representative.  <input type="checkbox"/> Not Applicable, Please Explain (e.g. no longer have the instrument):            _____         </p>			
<b>Customer number</b>		<b>Serial Number(s)</b>		
<b>Facility Name(s)</b>				
<b>Address</b>				
<b>City</b>		<b>State</b>		
<b>Phone Number</b>		<b>E-mail</b>		
<b>Name (print)</b>		<b>Title/Position</b>		
<b>Signature</b>				

**URGENT FIELD SAFETY NOTICE – FSN-RE23001**

For attention of: Healthcare Professionals

**Date Issued:** 31/08/2023

**Type of Action:** Field Safety Corrective Action

**Products:**

Product	Part Code	Lot No	Expiry Data	UDI from Kit Label
Multichem IA Plus	IA321A	0354062101	31.05.2024	(01)05391523442145(240)IA321A(17)230930(10)0354062101
Multichem IA Plus	IA322A	0354062102	31.05.2024	(01)05391523442152(240)IA322A(17)230930(10)0354062102
Multichem IA Plus	IA323A	0354062103	31.05.2024	(01)05391523442169(240)IA323A(17)230930(10)0354062103
Multichem IA Plus	IA321A	0368042201	31.03.2025	(01)05391523442145(240)IA321A(17)240731(10)0368042201
Multichem IA Plus	IA322A	0368042202	31.03.2025	(01)05391523442152(240)IA322A(17)240731(10)0368042202
Multichem IA Plus	IA323A	0368042203	31.03.2025	(01)05391523442169(240)IA323A(17)240731(10)0368042203
Multichem S Plus	CH114A.05	0116072101	31.12.2023	(01)05391523442107(240)CH114A.05(17)231231(10)0116072101
Multichem S Plus	CH115A.05	0116072102	31.12.2023	(01)05391523442114(240)CH115A.05(17)231231(10)0116072102
Multichem S Plus	CH116A.05	0116072103	31.12.2023	(01)08391523442138(240)CH116A.05(17)231231(10)0116072103
Multichem S Plus	CH114A.05	0123032201	28.02.2025	(01)05391523442107(240)CH114A.05(17)240831(10)0123032201
Multichem S Plus	CH115A.05	0123032202	28.02.2025	(01)05391523442114(240)CH115A.05(17)240831(10)0123032202
Multichem S Plus	CH116A.05	0123032203	28.02.2025	(01)05391523442138(240)CH116A.05(17)240831(10)0123032203

**Reason for Notice:**

Dear Valued Customer,

Technopath Manufacturing Ltd. is initiating a field action for the products listed above. This letter contains important information that needs your immediate attention.

<p><b>ISSUE:</b></p>	<p>Technopath Manufacturing Ltd. has identified an issue with the glass vial commodity used to fill the products as listed above. Some customers have experienced broken vials when receiving or thawing Technopath Multichem IA Plus controls (LN IA321A, IA322A and IA323A) or Multichem S Plus controls (LN CH114A.05, CH115A.05 and CH116A.05).</p> <p>Breakage has been observed upon delivery, while thawing, or while opening the vial.</p> <p>The breakage can result in the bottom of the vial breaking off, breakage in the neck thread area, or cracking on the sides of the vial.</p> <p>Technopath, Abbott and the vial supplier evaluated the vial manufacturing process for potential enhancements. Extensive investigation and verification studies were performed, and a strengthened glass vial design has now been qualified for use.</p> <p>Therefore, all future lots will be manufactured at Technopath Manufacturing Ltd. using the strengthened vials and will become available through the course of 2024.</p>
<p><b>IMPACT ON PATIENT RESULTS / USER SAFETY:</b></p>	<p>No patient results are impacted by this issue as there is no change to the product formulation.</p> <p>There is the potential for minor injuries or biohazardous exposure due to the possibility of vials breaking during thawing or handling.</p>
<p><b>NECESSARY ACTION TO BE TAKEN BY CUSTOMER:</b></p>	<p>Continue to adhere to the following points from the Procedures for Use in the respective Instructions for Use:</p> <ul style="list-style-type: none"> <li>➤ Allow the control vial to stand at room temperature (18 to 25°C) for 30 minutes or until completely thawed. The control must be completely solubilized before continuing to the next step.</li> </ul>




**CLINICAL  
DIAGNOSTICS**

Technopath Manufacturing Ltd  
Fort Henry Business Park  
Ballina Co. Tipperary  
Ireland  
SRN: IE-MF-000001822  
Ref: CC0159 Rev 01

	<p>➤ Place the thawed control vial on a suitable mixing device, e.g. roller mixer for 15 minutes, to achieve a homogenous solution.</p> <p>If you identify that a vial has broken or cracks are visible, safely dispose of the vial in line with Safety Data Sheet (SDS).</p> <p>Continue to use caution to prevent operator exposure and wear gloves when handling the control, as it contains material of human and animal origin. In addition, refer to the respective Safety Data Sheet (SDS) for more information regarding exposure controls and personal protection.</p> <p>Please share this information with your laboratory staff who need to be aware of this safety notice and retain this notification as part of your laboratory Quality System documentation. If you have forwarded any of the impacted product(s) listed above to another laboratory, please provide them a copy of this letter.</p>
<b>CONTACT</b>	<p>If you or any of the health care providers you serve have questions regarding this notice, customers please contact <a href="mailto:fieldaction@lgcgroup.com">fieldaction@lgcgroup.com</a></p>

We sincerely apologize for any inconvenience this may cause your laboratory and hope for your understanding and cooperation.

Yours Sincerely



Technopath is part of LGC Clinical Diagnostics

Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, Ireland. t: +353 61 525 700 | e: [technopath.info@lgcgroup.com](mailto:technopath.info@lgcgroup.com) | w: [www.technopathed.com](http://www.technopathed.com)

## ACKNOWLEDGEMENT

### URGENT FIELD SAFETY NOTICE – FSN-RE23001 - Action Required

Field Safety Notice Ref #RE23001

Date: 31<sup>st</sup> Aug 2023

Type of Action: Field Safety Corrective Action (FSCA)

Kindly complete and return this form to [fieldaction@lgcgroup.com](mailto:fieldaction@lgcgroup.com) before **21<sup>st</sup> September 2023**.

Product Catalogue Number:

Customer Name and Dept:

Address:

Are the above contact details correct? Yes  No

(If no please insert correct details below)

Contact Name:

Department:

Address:

Telephone Number:

Email:

*Please acknowledge receipt of information and awareness of any required actions described within the accompanying Field Safety Notice. Please bring this notice to the attention of all personnel in your hospital or healthcare facility who need to be aware of this safety issue. If you have forwarded the affected product(s) listed above to another laboratory, please provide a copy of this notice to them.*

I acknowledge receipt of this Field Safety Notice and have read, understood, and implemented its contents.

Name:

Signed: \_\_\_\_\_

Date: \_\_\_\_\_