



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

Field Safety Notice (FSN)

FRED easyport plus

manufactured by

SCHILLER AG, Altgasse 68 CH-6341 Baar, Switzerland

www.schiller.ch

SRN: CH-MF-000012722 / CHRN: CHRN-MF-20000327

Date: 2022-03-02

Attention: Schiller authorized distributors and their customers.

Issues related to the functionality of FRED easyport plus has been identified:

➔ Misleading Ready-To-Use Status and Analysis interruption.

This notice is intended to inform you about:

- what the problem is and under what circumstances it can occur.
- the actions that you as a distributor/customer can take to minimize the effect of the problem.
- the actions planned by SCHILLER AG to correct the problem.

We kindly ask that you read this notice carefully and send us written acknowledgement **by 15th of April 2022** that you have read and understood the contents of this notice. Written acknowledgement can be sent to SCHILLER AG via the contact vigilance@schiller.ch.

If you need any further information or support concerning this issue, please do not hesitate to contact SCHILLER AG Customer Services: support@schiller.ch

SCHILLER AG apologizes for any inconveniences caused by this problem.

Sincerely,

Valentina Shcherba
Head of Regulatory Affairs
Altgasse 68, CH-6341 Baar, Switzerland
vigilance@schiller.ch
T: +41 41 766 42 42



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

1. INFORMATION ON AFFECTED DEVICES	
COMMERCIAL NAME(S):	FRED easyport plus
PRIMARY CLINICAL PURPOSE OF DEVICE(S)*	The FRED easyport plus® is a defibrillator with the possibility to deliver a shock in semi-automatic, fully automatic or manual mode. FRED easyport plus® is intended to be used to terminate cardiac arrhythmia such as Ventricular Fibrillation (VF) or Ventricular Tachycardia (VT) with a defibrillation shock.
MODEL/CATALOGUE/ REF NUMBER(S):	REF : 3.940060, 3.940066, 3.940063 Catalogue number: 0A.900000
SOFTWARE VERSION:	All software versions below 1.2.2
AFFECTED SERIAL OR LOT NUMBER RANGE :	See list of affected products in Annex III
UNIQUE DEVICE IDENTIFIER(S) (UDI-DI):	07613365001921
DEVICE TYPE:	Automated-External defibrillator (AED)

2. REASON FOR FIELD SAFETY CORRECTIVE ACTION (FSCA)	
PROBLEM DESCRIPTION	The device sets the RTU (Ready-To-Use) state to OK by default unless the self-test fails. This results in the RTU state being incorrectly set to OK if the device crashes during the self-test or shuts off (e.g., because of insufficient battery capacity).
	If the device is operated with the volume set to the maximum the audio amplifier consumes more energy and generates more noise on the internal supply voltage. The devices can (depending on absolute patient impedance and language setting) recognize that noise as an ongoing chest compression and therefore can interrupt the analysis and does not recommend delivering a shock.
Hazard giving rise to the FSCA	Possible loss or reduction of the device's functionality.
PROBABILITY OF PROBLEM ARISING	The probability of the occurrence is unlikely, but expected for the devices during lifetime, when the battery gets a very low battery capacity and will be not exchanged as required by maintenance instructions or in case of specific languages where a language audio file volume is set to maximum.
PREDICTED RISK TO PATIENT/USERS	Delay to treatment is possible in case the user is misled by the device readiness status. An interruption of the analysis and no shock delivery might affect probability of survival.
	However, this risk of no shock delivery can be mitigated by mechanical resuscitation (chest compressions) by trained medical personnel / users of the device who will always need take into account the current clinical condition of the patient as per BLS instructions and guidance and perform clinical measures as needed, including CPR.



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

3. TYPE OF ACTION TO MITIGATE THE RISK	
ACTION TO BE TAKEN BY THE USER or AUTHORIZED DISTRIBUTOR / CUSTOMER	<ol style="list-style-type: none">1) Send the Annex I - Distributor Reply Form / Annex II - Customer Reply Form back to SCHILLER AG by 15th of April 2022 that the content of this notice was read and understood.2) Update the affected devices according to the Service Instructions by 15th of June 2022. Please note: The new Software version 1.2.2 and its Service Release Notes on how to perform the update are available on SCHILLER Service website: Software (schiller.ch) : https://extra.schiller.ch/products/Rescue/Pages/software.aspx?productId=3673) Send the Annex I - Distributor Reply Form / Annex II - Customer Reply Form including the confirmation about the updated devices back to SCHILLER AG by 15th of June 2022.
DATE FOR COMPLETION:	15th of June 2022
ACTIONS BEING TAKEN BY THE MANUFACTURER	<input checked="" type="checkbox"/> Software upgrade A new version of the software will be provided for the update of the device.
FURTHER INFORMATION AND SUPPORT	If you need any further information or support concerning this issue, please contact SCHILLER AG Customer Services: support@schiller.ch

Transmission of this Field Safety Notice

This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)
Please transfer this notice to other organisations on which this action has an impact. (As appropriate)
Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.
Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback. *

The responsible National Authority has been informed about this communication of this field safety notice.

Contact person of manufacturer:

Valentina Shcherba
Head of Regulatory Affairs
Altgasse 68, CH-6341 Baar, Switzerland
vigilance@schiller.ch
T: +41 41 766 42 42



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

ANNEX I
Distributor/Importer Reply Form

FRED easyport plus
SCHILLER AG, Altgasse 68 CH-6341 Baar Switzerland

Please complete, sign, and return your acknowledgement by 15th of April 2022 to vigilance@schiller.ch

1. Distributor/Importer Details	
Company Name*	
Account Number	
Address*	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

For the list of affected products, see ANNEX III

2. Distributors/Importers (Tick all that apply)		
<input type="checkbox"/>	*I confirm the receipt, the reading and understanding of the Field Safety Notice.	Distributor/Importer to complete or enter N/A
<input type="checkbox"/>	I have checked my stock and quarantined inventory	Distributor/Importer to enter quantity and date
<input type="checkbox"/>	I have identified customers that received or may have received this device	
<input type="checkbox"/>	I have attached customer list	
<input type="checkbox"/>	I have informed the identified customers of this FSN	Date of communication:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers	
<input type="checkbox"/>	I have returned affected devices - enter number of devices returned and date complete.	Add quantity, Lot/Serial Number, Date Returned (same information as requested by the Customer Reply form)
<input type="checkbox"/>	I have destroyed affected devices – enter number destroyed and date complete.	Add quantity, Lot/Serial Number, Date Returned/Destroyed (same information as requested by the Customer Reply form)
<input type="checkbox"/>	Neither I nor any of my customers has any affected devices in inventory	



SCHILLER

The Art of Diagnostics

SCHILLER AG | Altgasse 68 | 6341 Baar, Switzerland

Tel: +41 41 766 42 42 | Fax: +41 41 761 08 80

info@schiller.ch | www.schiller.ch

FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241

FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

Print Name*	Distributor/Importer print name here
Signature*	Distributor/Importer sign Here
Date *	

Mandatory fields are marked with *

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.
Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

ANNEX II
Customer Reply Form

FRED easyport plus

SCHILLER AG, Altgasse 68 CH-6341 Baar Switzerland

Please complete, sign, and return your acknowledgement by 15th of June 2022 to vigilance@schiller.ch

1. Customer Details	
Account Number	
Healthcare Organisation Name*	
Organisation Address*	
Department/Unit	
Shipping address if different to above	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

For the list of affected products, see ANNEX III

2. Customer action undertaken on behalf of Healthcare Organisation				
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content.*	Customer to complete or enter N/A		
<input type="checkbox"/>	I performed all actions requested by the FSN.*	Customer to complete or enter N/A		
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.*	Customer to complete or enter N/A		
<input type="checkbox"/>	I have returned affected devices - enter number of devices returned and date complete.	Qty:	Lot/Serial Number:	Date Returned (DD/MM/YY):
		Qty:	Lot/Serial Number:	Date Returned (DD/MM/YY):
		Or provide a separate list		
N/A	Comments:			
<input type="checkbox"/>	I have destroyed affected devices – enter number destroyed and date complete.	Qty:	Lot/Serial Number:	
		Qty:	Lot/Serial Number:	
		Or provide a separate list		
N/A	Comments:			



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

<input type="checkbox"/>	No affected devices are available for return/ destruction	Customer to complete or enter N/A
<input type="checkbox"/>	Other Action (Define):	
<input type="checkbox"/>	I do not have any affected devices.	Customer to complete or enter N/A
<input type="checkbox"/>	I have a query please contact me (e.g. need for replacement of the product).	Customer to enter contact details if different from above and brief description of query
Print Name*		Customer print name here
Signature*		Customer sign here
Date*		

Mandatory fields are marked with *

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.
Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

Annex III List of affected devices

Serial Number	Date of the Update / Comments
9000.000650	
9000.000653	
9000.000654	
9000.000655	
9000.000656	
9000.000657	
9000.000658	
9000.000659	
9000.000660	
9000.000661	
9000.000662	
9000.000663	
9000.000664	
9000.000815	
9000.000816	
9000.000817	
9000.000818	
9000.000819	
9000.000820	
9000.000821	
9000.000822	
9000.000823	
9000.000824	
9000.000825	
9000.000826	
9000.000827	
9000.001514	
9000.001515	
9000.001516	
9000.001517	
9000.001518	
9000.001519	
9000.001520	
9000.001521	
9000.001522	
9000.001523	
9000.001524	
9000.001525	
9000.001526	
9000.001527	
9000.001528	
9000.000337	
9000.000430	
9000.000431	
9000.000432	
9000.000433	



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

9000.000434	
9000.000435	
9000.000436	
9000.000437	
9000.000438	
9000.000439	
9000.000440	
9000.000441	
9000.000442	
9000.000443	
9000.000444	
9000.000445	
9000.000446	
9000.000447	
9000.000448	
9000.000449	
9000.000450	
9000.000451	
9000.000452	
9000.000453	
9000.000454	
9000.000455	
9000.000456	
9000.000457	
9000.000458	
9000.000459	
9000.000460	
9000.000461	
9000.000462	
9000.000463	
9000.000464	
9000.000465	
9000.000466	
9000.000467	
9000.000468	
9000.000469	
9000.000470	
9000.000471	
9000.000472	
9000.000473	
9000.000474	
9000.000475	
9000.000476	
9000.000477	
9000.000478	
9000.000479	
9000.000511	
9000.000512	
9000.000513	
9000.000514	
9000.000515	



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

9000.000516	
9000.000517	
9000.000518	
9000.000519	
9000.000520	
9000.000521	
9000.000522	
9000.000523	
9000.000524	
9000.000525	
9000.000526	
9000.000527	
9000.000528	
9000.000529	
9000.000530	
9000.000531	
9000.000532	
9000.000533	
9000.000534	
9000.001789	
9000.001791	
9000.001792	
9000.001794	
9000.001795	
9000.001796	
9000.001797	
9000.001798	
9000.001799	
9000.001800	
9000.001801	
9000.001802	
9000.001803	
9000.001804	
9000.001805	
9000.001806	
9000.001807	
9000.001808	
9000.001809	
9000.001810	
9000.001811	
9000.001812	
9000.001813	
9000.001814	
9000.001815	
9000.001816	
9000.001817	
9000.001818	
9000.001819	
9000.001820	
9000.001821	
9000.001822	



FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241
FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

9000.001823	
9000.001926	
9000.001927	
9000.001928	
9000.001929	
9000.001930	
9000.001931	
9000.001932	
9000.001933	
9000.001934	
9000.001935	
9000.001936	
9000.001937	
9000.001938	
9000.001939	
9000.001940	
9000.001941	
9000.002129	
9000.002130	
9000.002131	
9000.002132	
9000.002133	
9000.002134	
9000.002135	
9000.002136	
9000.002137	
9000.002138	
9000.002139	
9000.002140	
9000.002141	
9000.002142	
9000.002143	
9000.002144	
9000.002145	
9000.002146	
9000.002147	
9000.002148	
9000.002149	
9000.002150	
9000.002151	
9000.002152	
9000.002153	
9000.002154	
9000.002155	
9000.002156	
9000.002157	
9000.002158	
9000.002159	
9000.002160	
9000.002161	
9000.002162	



SCHILLER

The Art of Diagnostics

SCHILLER AG | Altgasse 68 | 6341 Baar, Switzerland

Tel: +41 41 766 42 42 | Fax: +41 41 761 08 80

info@schiller.ch | www.schiller.ch

FSCA Ref: CAPA SAGQI-482, SAGQI-504 // # 1837771, 1839070, 1839241

FSN_with acknowledgment form_Fred easyport plus_SAGQI-482_SAGQI-504_EN.docx

9000.000634	
9000.000801	
9000.001283	
9000.001284	
9000.001285	
9000.001286	
9000.001287	
9000.001720	
9000.001721	
9000.001832	
9000.001833	
9000.001834	