Content Customer information cementless Bicontact SD hip stems

Weight restriction for Bicontact SD Implants sizes 9 and 10 and Bicontact SD-Revision size 11 and instructions for intraoperative technique of femoral canal preparation

Dear user of the Bicontact hip system,

within the Bicontact hip system the implants of the SD- and N-Series offer a stem range for the treatment of severe dysplastic hip anatomies as well as for narrow femoral canal conditions. Within more than 300.000 Bicontact implantations the stems of the SD- and N-series represent today a percentage of 5%.

In case of heavy patients small implant sizes are exposed to high material loads. Especially under conditions as distal pressfit and/or insufficient proximal bone support combined with higher offset loads the prosthesis stem can be overloaded in vivo.

New preclinical testing conditions with hip simulators have been set with the ISO 7106-04. This new test set-up leads to a higher load applied to the hip implant stems. Three implants of the Bicontact hip system (Bicontact SD Implants sizes 9 and 10 and Bicontact SD-Revision size 11) did not pass these new loading conditions. Therefore these implants can be only offered with weight restrictions for patient with hip arthroplasty treatment.

Clinical experiences have shown, that in particular cases the SD prosthesis stems failed due to fatigue overloading. The analysis of these cases had always shown a strong distal press-fit and load transfer situation which represented the laboratory test situation of the new ISO test-set-up.

As Bicontact implants are intended to be implanted with a proximal load transfer situation, distal reaming is recommended in case of narrow distal femoral canal conditions. This surgical procedure has been always described in the Bicontact surgery manual, because the use of the A-Osteoprofilers can not ensure a distal "cutting-free" of the femoral canal in all cases of small bone conditions. An alternative to the use of small SD sizes also N implants which can be used in such indications because these stems do not have any weight restrictions.

Due to CE-approval and CE regulations we are committed to add the weight limitation to the instructions for use. This letter is intended to inform you about these contraindications and the content of the Bicontact SD instructions for use. Bicontact D implants are exclusively delivered to Japan and therefore also part of this instruction for use.

If you have any additional question please do not hesitate to contact your B. Braun Aesculap branch or contact Aesculap Tuttlingen per e-mail to thomas.guettler@aesculap.de.

Enclosure: Supplement to instructions for use Bicontact SD and D implants

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TA-Nr. 011574 10/04

Contraindications – Weight restriction

Implant failure due to excessive body weight of the patient and mainly distal fixation of the small BiCONTACT SD implants!

- Observe maximum body weight of the patient for the implants listed below.

Art. no.	Designation	Max. body weight
NK709T	BiCONTACT SD, size 9	50 kg
NK710T	BiCONTACT SD, size 10	50 kg
NK210T	BiCONTACT SD, Revision, size 11	65 kg
NJ208T	BiCONTACT D, size 8	60 kg
NJ209T	BiCONTACT D, size 9	60 kg

To avoid distal implant fixation, carry out the following additional surgical measures:

- Bore up the distal marrow cavity, using flexible marrow cavity drills.
- Ensure sufficiently deep prosthesis seat through proximal-medial support.
- Carry out intraoperative check of the prosthesis seat, using an imaging device.
- Alternatively, use BiCONTACT N implants.

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