



Field Safety Notice

Device: Uniglide Unicondylar Knee, Femoral Component
Part Number 514.0400, LOT Number KNRU, Quantity 7 units
514.0300, LOT Number KNTK, Quantity 9 units

Type of action: Voluntary device withdrawal

Date: 14 April 2009

To:

Details on affected devices:

This Field Safety Notice concerns the above devices ONLY. The devices are size 3 & 4 femoral components used in a non-cemented total unicondylar knee arthroplasty. Our records indicate that you are in receipt of qty devices from the above affected lots.

No other devices within the range are affected. Only these two part number/lot number combinations are affected.

Description of the problem:

The nonconforming devices may have been comingled during the manufacturing process.

1. LOT number KNRU, there may be devices in unit boxes with labels stating the device to be size 4, the device itself marked as size 4 whereas the device is actually size 3.
2. LOT number KNTK, there may be devices in unit boxes with labels stating the device to be size 3, the device itself marked as size 3 whereas the device is actually size 4.

Clinical Implications

It is highly probable that a surgeon would identify this issue before implanting.

Outcome for a size 4 device located on a femur prepared for size 3:

The surgeon would find that despite the trial being a snug fit, the prosthesis would be loose with a significant gap between posterior cut and prosthesis. There would also be poor radial fit between the modelled surfaces of the femur and the prosthesis.

Outcome for a size 3 device located on a femur prepared for size 4:

The surgeon would find that, despite the trial being a good fit, the prosthesis would not seat on the femur.

Instruction

1. Discontinue use of the above devices, remove them from your inventory and place in a quarantine area pending return to Corin Ltd.
2. Corin Ltd customer representatives shall contact you during the week starting 14th April regarding the removal of the incorrect device and arranging for a replacement device. Alternatively please contact us on the Customer Hotline Tel: +44(0)1285 649 231.
3. Please forward this safety notice to any health professional in your organisation that needs to be aware.


Corin Limited
The Corinium Centre
 Cirencester
 Gloucestershire
 GL7 1YJ
t: +44 (0)1285 659866
f: +44 (0)1285 658960
e: info@coringroup.com
www.coringroup.com

Registered in England, No. 1610553

Your Competent Authority has been notified of this voluntary withdrawal.

If you have any questions or would like assistance with this issue please contact either your local Corin sales executive or the abovementioned customer services hotline.

Yours sincerely

A large black rectangular redaction box covers the signature area. A small handwritten mark is visible above the top right corner of the box.

Mr. [REDACTED]
Quality Assurance Manager
Corin Ltd.