



FSN Ref: FSN\_ Titan Pump\_20230215

FSCA Ref: FSCA\_ Titan Pump\_20230215

Date: 03.02.2023

**Urgent Field Safety Notice – Recall of specific lot numbers**  
**Titan Inflatable Penile Prosthesis**

Contact details of local representative (name, e-mail, telephone, address etc.)*
----------------------------------------------------------------------------------

Kyra Sievert service@coloplast.com
------------------------------------

**Urgent Field Safety Notice (FSN)**  
**Titan Inflatable Penile Prosthesis**  
**Risk addressed by FSN**

<b>1. Information on Affected Devices*</b>			
<b>1</b>	<b>1. Device Type(s)*</b>		
.	<p>This FSN concerns the Titan Inflatable Penile Prosthesis (IPP) pumps, which are part of a Titan IPP hydraulic system. The Titan implant consists of two inflatable Bioflex® penile cylinders that are implanted in the corpora cavernosa of the penis. The cylinders are attached to a pump, which is placed in the patient's scrotum, and the pump is connected to a fluid reservoir that is implanted underneath the abdominal muscles. The fluid reservoir contains a Lock-Out™ valve, which is intended to minimize the opportunity for auto-inflation. The fluid reservoir is filled with a sterile saline solution. Repetitive squeezing of the pump bulb transfers fluid from the reservoir to the cylinders in the penis. As the penile cylinders fill with fluid, the penis enlarges and becomes erect, thereby facilitating intercourse.</p> <p>The Titan IPP is a three-piece inflatable penile prosthesis (IPP) consisting of a pump, cylinder assemblies and reservoir. The pump and cylinders are pre-connected, and the reservoir is connected to the pump inlet tubing prior to implantation, using components of the assembly kit.</p>		
<b>1</b>	<b>2. Commercial name(s)</b>		
.	Titan Inflatable Penile Prosthesis		
<b>1</b>	<b>3. Primary clinical purpose of device(s)*</b>		
.	The Titan IPP is a self-contained hydraulic system designed to be surgically implanted for the management of erectile dysfunction. This implant provides the patient with voluntary control over the erect and flaccid states of the penis.		
<b>1</b>	<b>4. Device Model/Catalogue/part number(s)*</b>		
.	ES2918, ES2922, ES2920, EN2814, ES2916, EN2816, ES2914, EN2911, EN2918		
<b>1</b>	<b>5. Affected serial or lot number range</b>		
.	<b>Lot Serial Number</b>	<b>Item Number</b>	<b>Expiration Date</b>
	8812159	ES2918	Oct 5, 2027
	8840683	ES2922	Oct 12, 2027
	8812162	ES2920	Oct 6, 2027
	8849596	ES2922	Nov 3, 2027
	8852985	ES2918	Oct 24, 2027
	8852984	ES2918	Oct 24, 2027
	8887598	EN2814	Nov 16, 2027
	8849590	ES2918	Oct 27, 2027
	8840681	ES2916	Oct 12, 2027
	8812163	ES2920	Oct 6, 2027
	8812164	ES2922	Oct 6, 2027
	8812161	ES2920	Oct 6, 2027
	8812160	ES2920	Oct 6, 2027
	8849612	ES2920	Nov 6, 2027
	8812165	ES2922	Oct 6, 2027
	8849568	ES2916	Oct 20, 2027

8849594	ES2920	Nov 3, 2027
8849593	ES2920	Nov 3, 2027
8840682	ES2918	Oct 12, 2027
8849607	ES2916	Oct 31, 2027
8853029	ES2920	Oct 24, 2027
8895171	EN2816	Nov 17, 2027
8849597	ES2922	Nov 3, 2027
8849595	ES2920	Nov 3, 2027
8840680	ES2914	Oct 12, 2027
8849608	ES2918	Nov 6, 2027
8849583	EN2911	Oct 26, 2027
8849567	ES2916	Oct 20, 2027
8849574	EN2918	Oct 24, 2027
8849592	ES2918	Nov 3, 2027
8849610	ES2920	Nov 6, 2027
8849613	ES2922	Nov 6, 2027
8849591	ES2918	Nov 3, 2027
8812158	ES2918	Oct 5, 2027
8849611	ES2920	Nov 6, 2027

<b>2 Reason for Field Safety Corrective Action (FSCA)*</b>	
<b>2</b>	<p><b>1. Description of the product problem*</b></p> <p>Coloplast has identified that Titan Touch Pumps manufactured between September 17, 2022 through December 2, 2022 have a decreased wall thickness (compared to the current standard) and are therefore, subject to this voluntary recall. A decreased wall thickness may result in difficulty inflating and/or deflating the device, pump failure, or a fracture of the pump.</p>
<b>2</b>	<p><b>2. Hazard giving rise to the FSCA*</b></p> <p>Pump failure may present as a fracture of the pump wall and may compromise or prevent device function. In the event of a pump fracture, solution used to fill the IPP device may leak. This fluid is physiological saline (sterile, isotonic, pyrogen-free Sodium Chloride U.S.P. Solution for Injection) and is not a source of risk or harm in the case of leakage within the body. Temporary swelling may follow the fracture. Additional potential risks from pump fracture are those associated with replacement surgery if this is the course of action determined by you and your patient. Pump failure may present as the inability to inflate or deflate the device due to loss of the fill solution or weakened wall. If you suspect a patient has a device affected by this issue, it is recommended that you manage the patient as you would in the normal course of clinical practice. If a replacement is necessary, please report the issue immediately to <a href="mailto:FRcomplaints@coloplast.com">FRcomplaints@coloplast.com</a> and return the explanted device.</p>
<b>2</b>	<p><b>3. Probability of problem arising</b></p> <p>This voluntary recall is a proactive measure conducted by Coloplast. There have been no reports of pump failure directly related to this issue as of now. Pump longevity is impacted by the number of inflate/deflate cycles the device undergoes over its lifespan, which may vary by user. Testing confirms Titan Touch pumps manufactured prior to September 17, 2022, and after December 2, 2022, are not affected.</p>
	<p><b>4. Background on Issue</b></p>



<b>4. General Information*</b>	
4.	1. FSN Type* <span style="float: right;">New</span>
4.	2. Further advice or information already expected in follow-up FSN? * <span style="float: right;">No</span>
4.	3. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)
	a. Company Name <span style="float: right;"><b>Coloplast A/S</b></span>
	b. Address <span style="float: right;"><b>Holtedam 1 3050 Humlebæk Denmark</b></span>
4.	4. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *
4.	5. List of attachments/appendices: <span style="float: right;"><b>Customer Reply Form</b></span>
4.	6. Name/Signature

<b>Transmission of this Field Safety Notice</b>	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback..*</p>

Note: Fields indicated by \* are considered necessary for all FSNs. Others are optional.

## Field Safety Notice Customer Reply Form

It is important that the actions detailed in the FSN are taken and the completed form is returned to Coloplast.

This reply is the evidence we need to monitor the progress of the corrective actions taken.

**If you have further distributed this device, please forward this form to whom you have distributed the device too. They can utilize this form to communicate status of the product they have received.**

1. Field Safety Notice (FSN) information		
FSN Reference number*	FSN_Titan Pump_20230215	
FSN Date*	March 02, 2023	
Product/ Device name*	Titan Inflatable Penile Prosthesis	
Product Code(s)	ES2918, ES2922, ES2920, EN2814, ES2916, EN2816, ES2914, EN2911, EN2918	
Batch/Serial/Lot Number (s)	8812159, 8840683, 8812162, 8849596, 8852985 8852984, 8887598, 8849590, 8840681, 8812163 8812164, 8812161, 8812160, 8849612, 8812165 8849568, 8849594, 8849593, 8840682, 8849607 8853029, 8895171, 8849597, 8849595, 8840680 8849608, 8849583, 8849567, 8849574, 8849592 8849610, 8849613, 8849591, 8812158, 8849611	
2. Customer Details		
Account Number		
Healthcare Organisation or Company Name*		
Address*		
Department/Unit		
Shipping address if different to above		
Contact Name*		
Title or Function		
Telephone number*		
3. Return acknowledgement to sender		
E-mail	Kyra Sievert service@coloplast.com	
Deadline for returning the Reply form*	April 14 <sup>th</sup> , 2022	
4. Customer/Distributor/Importer's actions undertaken (Tick all that apply)		
<input type="checkbox"/>	I confirm the receipt, the reading and understanding of the Field Safety Notice.	Complete or enter N/A
<input type="checkbox"/>	I have performed all actions requested by the FSN.	Complete or enter N/A
<input type="checkbox"/>	I have checked my stock and quarantined inventory	Complete or enter N/A
<input type="checkbox"/>	I have identified customers that received or may have received this device	Relevant for Distributor and Importer
<input type="checkbox"/>	I have attached customer list	Relevant for Distributor and Importer
<input type="checkbox"/>	I have informed the identified customers of this FSN	Date of communication:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers	
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.	Relevant for Customer
<input type="checkbox"/>	I have returned affected devices	Qty:                      Lot no.:                      Date Returned (DD/MM/YY): Comments:
<input type="checkbox"/>	I have destroyed affected devices	Qty:                      Lot no.:                      Date Returned (DD/MM/YY): Comments:
<input type="checkbox"/>	No affected devices are available for return/ destruction	Complete or enter N/A
<input type="checkbox"/>	I have used the affected devices.	Qty:                      Lot no.:                      Date Returned (DD/MM/YY):

		Comments:
<input type="checkbox"/>	Other Action (define):	
<input type="checkbox"/>	I have a query please contact me (e.g. need for replacement of the product).	Enter contact details if different from above and brief description of query
Print Name*		
Signature*		
Date *		

\*Mandatory fields