

### Video Processor / Illuminator

**Recall Number: RA2021-2640074**

**March 19, 2021**

**Attn: Recall Coordinator**

**Account number:**

**Account name:**

**Account Address:**



### Response required

Catalog number	UDI	Product description	Serial number(s)
PC9001	00858701006049	Video Processor/Illuminator (VPI)	See Attachment A

#### Product description

A video processor/illuminator (VPI), provides VIS/NIR (Visible/Near Infrared) illumination to the surgical laparoscope via a flexible light guide cable, as well as the image processing necessary to generate simultaneous, real-time HD (High Definition) video color and NIR (Near Infrared) fluorescence images.

#### Intended Use

Since the VPI is a component that is used with both the SPY-PHI and PINPOINT systems, the intended use for both imaging systems are being documented here:

##### *PINPOINT SYSTEM:*

PINPOINT Endoscopic Fluorescence Imaging System is intended to provide real-time endoscopic visible and near infrared fluorescence imaging. This imaging system enables surgeons to perform routine visible light endoscopic procedures as well as further visually assess circulation including blood flow in vessels and microvessels, tissue and organ perfusion, and lymphatics and perfusion associated with tumors and tumor margins with near infrared fluorescence imaging during minimally invasive surgery.

##### *SPY-PHI SYSTEM:*

The SPY Portable Handheld Imaging System (SPY-PHI System) is an active device used to visualize circulation, including lymphatics and blood vessels, as well as related tissue perfusion with near infrared fluorescence imaging during a variety of surgical procedures.

#### Product issue

It was discovered that a portion of the VPIs shipped between 10-Jan-2020 and 02-Feb-2021 are exhibiting flickering, loss of image or unintended VPI reboots.

#### Potential risks and hazards

When used with SPY-PHI, if loss of image or unintended VPI reboots occur during surgery, the potential harm is a delay of surgery. When used with PINPOINT, if loss of image or unintended VPI reboots occurs during surgery, the potential harm is conversion to open. The potential of harm is remote.

## **Actions needed**

1. Immediately check inventory to locate the product listed on Attachment A and **remove them from their point of use**.
2. Return the enclosed business reply form by email to [xxxxxxx@stryker.com](mailto:xxxxxxx@stryker.com) to confirm receipt of this notification/document product segregation. **RESPONSE IS REQUIRED.**
  - a. If affected product is found, segregate the product and call Stryker customer service or email [xxxxxxx](mailto:xxxxxxx) to arrange for product return and receive a replacement (upon availability) and complete Business Reply form, Attachment B on page 4.
  - b. If no product is found, complete Business Reply form, Attachment B on page 4 of this recall letter.
3. Maintain awareness of this communication internally until all required actions have been completed.
4. Inform Stryker if any of the subject devices have been distributed to other organizations.
  - a. Please provide contact details so that Stryker can inform the recipients appropriately.
  - b. If you are a Distributor, note that you are responsible for notifying your affected customers.
5. Please inform Stryker of any adverse events concerning the use of the subject devices.
  - a. Please comply with any local laws or regulations concerning the notification of adverse events to your National Competent Authority.
6. Complete the attached customer response form. It may be that you no longer have any physical inventory on site. Completing this form will allow us to update our records and will also negate the need for us to send any further unnecessary communications on this matter. Therefore, please complete even if you no longer have any of the subject devices in your physical inventory.
7. Return the completed form to your nominated Stryker Representative (indicated below) for this PFA
  - a. On receipt of the form, a Stryker Representative will contact you to organize any applicable ongoing actions.

We request that you respond to this notice within XXX calendar days from the date of receipt.

The target date for completion of this action is XXX and your timely response will enable us to ensure that we meet this target/

This will enable us to update our records and negate the need to send unnecessary reminder letters.

Your designated contact person for this action is given below. Should you have any queries concerning this matter please do not hesitate to contact them directly.

**Name:** \_\_\_\_\_ **Position:** \_\_\_\_\_ **email:** \_\_\_\_\_

In line with the recommendations of the Meddev Vigilance Guidance document Ref 2.12-1, we can confirm that this FSCA has been notified appropriately to the National Competent Authority for your country.

On behalf of Stryker we thank you sincerely for your help and support in completing this action within the target date and regret any inconvenience that may be caused. We would like to reassure you that Stryker is committed to ensuring that only conforming devices, meeting our high internal quality standards, remain on the market.

Sincerely,

## Attachment A

### Recalled Serial Numbers

**Device Description and Part Number:** PC9001 – Video Processor/Illuminator (VPI)

**Impacted Serial Numbers:**

PI20362910	PI20082806	PI20112868	PI20322887	PI20372920	PI16180922	PI20052735	PI20362908
PI20362911	PI20092807	PI20112869	PI20322888	PI20372921	PI17381211	PI20052740	PI20362909
PI20362913	PI20102835	PI20112870	PI20322889	PI20372922	PI17381214	PI20062759	PI20463033
PI20362912	PI20102836	PI20112871	PI20322890	PI20372923	PI18041369	PI20062766	PI20463038
PI20022723	PI20102839	PI20112872	PI20322891	PI20052734	PI18041410	PI20062768	PI20463041
PI20022724	PI20102846	PI20112873	PI20362892	PI20052736	PI18181490	PI20062771	PI20463042
PI20022725	PI20102851	PI20112874	PI20362893	PI20062747	PI18241579	PI20062772	PI20463046
PI20022726	PI20112852	PI20112875	PI20362894	PI20062754	PI19182075	PI20072773	PI20463049
PI20062749	PI20112860	PI20112876	PI20362895	PI20082802	PI19182080	PI20072780	PI20463051
PI20082795	PI20112861	PI20282879	PI20362896	PI20092808	PI19242192	PI20072781	PI20463052
PI20082797	PI20112862	PI20282880	PI20372914	PI20022701	PI19342372	PI20072784	PI20463054
PI20082798	PI20112863	PI20312881	PI20372915	PI20082788	PI19372413	PI20072787	PI20463058
PI20082799	PI20112864	PI20312882	PI20372916	PI20463062	PI19432511	PI20362904	PI20463036
PI20082801	PI20112865	PI20312883	PI20372917	PI14350291	PI20052727	PI20362905	
PI20082804	PI20112866	PI20322885	PI20372918	PI15150533	PI20052728	PI20362906	
PI20082805	PI20112867	PI20322886	PI20372919	PI15240575	PI20052732	PI20362907	

# Business Reply Form

Account number:  
 Account name:  
 Account Address:

## Attachment B - Business Reply Form

### Video Processor / Illuminator

**Recall Number: RA2021-2640074**

**March 19, 2021**

Please complete and sign this form. Email the completed form to [endorecall@stryker.com](mailto:endorecall@stryker.com) by **01-June-2021**.

**Note:** Your signature indicates that you have received and understand the enclosed notification and that you have performed all actions requested.

Catalog number	Product	Serial number(s)	Quantity on hand
PC9001	VIDEO PROCESSOR/ILLUMINATOR (VPI)	List here	

**Form completed by:**

Printed Name		Title	
Signature		Phone	
Date		Email	

- If you no longer have affected product on hand, please check here.
  - Please state disposition of product no longer on hand: \_\_\_\_\_

If you have further distributed any affected product, please indicate to whom:

Product(s) Distributed		Quantity Distributed	
Facility Name		Contact Person	
Full Address			

Return completed Business Reply Form to [XXXXXXXXXXXXXXXXXXXX](#)