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Dear User,

We would like to inform you that our a.m. product, Artikon.genu-bi-lock (System-KneeJoint Set with locking mechanism) possibly carries a safety risk.

Therefore, please kindly note the following information and safety notes. The serial numbers concerned are indicated on page 5 of this information.

Please accept our sincere apologies for any inconvenience caused and do not hesitate to contact us anytime, if you have any questions.

Problem description

The Artikon.genu-bi-lock system-knee joint set with locking mechanism is a component used in knee-ankle-foot-orthoses (KAFO) and is designed to prevent unintended knee flexion. We have discovered that under certain circumstances the locking mechanism may fail, which in the worst case could cause the user to fall and consequently lead to injuries.

As a precautionary measure, we have decided to recall all

Artikon.genu-bi-lock system-knee joint sets with locking mechanism intended for system widths 16 and 20 mm immediately, and to stop supplying them until the problem has been cleared.

This recall concerns all products supplied during the period from 2018-05-14 until 2019-01-14.

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Instructions regarding actions to be taken by the customer

All users of knee-ankle-foot-orthoses (KAFO) fitted with an Artikon.genu-bi-lock system-knee joint set with locking mechanism intended for system widths 16 and 20 mm are asked to contact their orthopaedic workshop immediately.

The Artikon.genu-bi-lock system-knee joint must be removed from the orthosis; and the leg orthoses concerned must be equipped with any other commercial knee joint with locking mechanism. In order to eliminate any risk of falling, we urgently ask you not to use the leg orthosis any more until the joint has been exchanged.

Subsequently, we kindly ask all orthopaedic workshops and orthotists to return the affected Artikon.genu-bi-lock system-knee joint sets with locking mechanism to Streifeneder ortho.production GmbH, Moosfeldstrasse 10, 82275 Emmering, Germany.

Streifeneder ortho.production GmbH would like to support you during the modification of the knee-ankle-foot-orthoses (KAFO) with alternative knee joints with locking mechanism and we would like to give, as always, individual recommendations in order to find the best possible solution for each patient. Of course, we will provide support for the costs of the necessary products exchanged.

This precautionary measure and recall is valid from 2019-01-14 for an indefinite period.

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Customer reply

Product:

Artikon.genu-bi-lock System-Knee Joint Set with locking mechanism

Item number

- 2J30/16L
- 2J30/20L
- 2J30/16R
- 2J30/20R
- 2J31/16
- 2J31/20

Dear Sir or Madam,

According to the safety information, please complete this reply form and return it immediately to Streifeneder ortho.production GmbH, preferably by e-mail (d.rudolph@streifeneder.de).

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☐ Ja ☐ Nein

I have read and understood the recall instruction provided in the safety information.

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Record item number and serial number of the product returned.
Please use a separate sheet if necessary.

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Wird von Streifeneder ortho.production GmbH ausgefüllt / Will be filled out by Streifeneder ortho.production GmbH
Unterschrift Mitarbeiter Streifeneder ortho.production GmbH