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### Betreffend

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Dear User,

We would like to bring to your attention information  
concerning the above-mentioned product, GO.smart,  
Carbon-fiber foot, weight category 60 kg, which has  
presented safety concerns under certain circumstances.

**Please observe and adhere to the following infor-  
mation and instructions.**

We apologize for any inconvenience this may cause.  
If you have any further enquiries or need additional  
information please do not hesitate to contact us.

**Problem description**

The GO.smart Carbon-fiber Foot, weight category  
60 kg, is a foot component which is used in lower limb  
prosthesis.

This prosthetic foot is approved for users with activity  
levels 1 and 2 according to the Streifeneder Walky-system  
and for maximum body weight of up to 60 kg, including  
carried and supported loads.

Internal test results show that the elastomer bumper/  
wedge used is not sufficiently stable under certain  
conditions and can tear when bearing overload.

A defective foot could lead to a loss of functionality  
and stability, and should be replaced by an orthopedic  
workshop.

As a precautionary measure, all affected GO.smart  
Carbon-fiber Feet, weight category 60 kg are recalled  
and will be replaced by a new version of the GO.smart  
Carbon-fiber Foot, weight category 60 kg.

Affected are all GO.smart Carbon-fiber Feet, weight  
category 60 kg, which were delivered in the period from  
2018-06-01 to 2020-02-01.

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### **Instructions regarding actions to be taken by the customer**

All users of lower limb prostheses, which are fitted with GO.smart, Carbon-fiber Feet with weight category 60 kg, are asked to immediately contact their orthopedic workshop.

The affected lower limb prostheses have to be refitted with a replacement-GO.smart Carbon-fiber Foot with weight category 60 kg. We will provide the replacement foot free of charge if the customer reply is completely filled out attached below.

Subsequently we kindly ask all orthopedic workshops and prosthetists to send back the affected GO.smart Carbon-fiber Feet with weight category 60 kg to Streifeneder ortho.production GmbH, Moosfeldstraße 10, 82275 Emmering, Germany. A return label will be provided.

Streifeneder ortho.production GmbH would like to support you during the modification of the lower limb prosthesis with a replacement GO.smart, Carbon-fiber Foot.

This precautionary measures and recall is valid from 2020-08-18 to an indefinite time point.

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5A410/...

Dear Sir or Madam,

According to the safety information, please complete  
this reply form and return it immediately to  
Streifeneder ortho.production GmbH,  
preferably by e-mail (d.rudolph@streifeneder.de).

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I have read and understood the recall instruction provided  
in the safety information.

Yes      No

All affected users have received a  
replacement product.

Yes      No

Record item number and serial number of the product  
returned.  
Please use a separate sheet if necessary.

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Wird von Streifeneder ortho.production GmbH ausgefüllt / Will be filled out by Streifeneder ortho.production GmbH
Unterschrift Mitarbeiter Streifeneder ortho.production GmbH