



Meditrade GmbH · Medipark 1 · 83088 Kiefersfelden / Germany

Date
12.08.2021

Emergency - Safety Notice Nitril Sensory white and Nitril Sensory blue, LOT 3220003750820 - Holes in gloves

Dear customer,

Meditrade GmbH is publishing a field safety notice for the following products:
Examination gloves Nitril® Sensory® white, REF 2280, LOT 3220003750820
Examination gloves Nitril® Sensory® blue, REF 2283, LOT 3220003750820

The gloves in this batch have significantly more holes than would be permissible according to an AQL of 1.5. As a result, protection against contamination or chemicals is not adequately ensured.

The following actions should be taken immediately:

Stop using the gloves, Nitrile Sensory, of LOT 3220003750820.

Contact Meditrade or your distributor if you still have this product.

Further information:

If you require any further information or assistance in relation to this issue, please contact:

reklamation@meditrade.de

Or by telephone: 08033 9760 0

This field safety notice has also been reported to the relevant authorities.

Meditrade regrets any inconvenience caused by the issue.

Yours sincerely

