

.....
.....
- OP- Leitung -
.....
.....

Ihre Zeichen

Unsere Zeichen
QMV / IL

Tel. (0 74 61) 95-
2922


Fax (0 74 61) 95-
2177

Datum
2006-02-03

Recall FF013P Scalpfix – manufacturing lots 51291714 and 51280476

Dear Sirs,

through feedback from one of our customers, we received an information that the blister package of the above mentioned product could be damaged.



Aesculap.
All it takes to operate.

Vorsitzender des Beirates:
Dr. Fritz Lamb
Geschäftsleitung:
Prof. Dr. Dr. Dr. h.c. Michael Ungethüm
(Vorsitzender und
pers. haftender Gesellschafter)
Christof Hennigfeld
Dr. Harald Stallforth

Sitz der Gesellschaft: Tuttlingen
Reg. Gericht: Tuttlingen HRA 902
Steuernummer: 21060/00028

Bankverbindungen:
Deutsche Bank AG Tuttlingen
BLZ 653 700 75, Konto 21 22 000
IBAN DE44 6537 0075 0212 2000 00
Dresdner Bank AG Tuttlingen
BLZ 643 800 11, Konto 27 23 001
IBAN DE05 6438 0011 0272 3001 00

Anschrift:
Aesculap AG & Co. KG
Am Aesculap-Platz
D-78532 Tuttlingen

...

Our analysis revealed, that in specific cases (manufacturing lots 51291714 and 51280476) the sealed seam on the blister package is insufficient and there might be a possibility that the sterility of the product is not guaranteed.

To avoid the risk of non sterility, we decided to recall the product concerned.

Please send us the product back with enclosed recall form. You will receive replacement free of charge.

We apologize for any inconvenience it may have caused. For further questions please do not hesitate to contact us.

Best regards

Aesculap AG & Co. KG

i. V.

i. A.

Susanne Napravnik
(Product Management Neurosurgery)

Ines Lucchese
(Recall Coordinator Quality Management)