



SURGE

MEDICAL SOLUTIONS, LLC

April 17, 2007

[REDACTED]
[REDACTED]
[REDACTED]

GERMANY

RE: Medical Device Recall

We have initiated a recall of the following Surge Medical Solutions Self-Inflating Retrograde Cannulae, with lot numbers beginning with digits **0610, 0611, 0612, 0701, or 0702:**

RSG-M014S, RSG-M014S-L, RSG-T014S-L, RSH-L014S, RSH-M014S, RSH-M014S-L, RSH-S014S, RSH-T014S, RSH-T014S-L, RGR-L014S, RGR-M014S, RGR-M014S-L, RSR-T014S

The products are being recalled because of a potential leak in the pressure monitoring line, which could result in inaccurate pressure readings.

Enclosed is a Return Goods Authorization (RGA) form. It lists the products included in the recall that have been shipped to you, or in some instances, directly to your customer. We request your cooperation in determining the products listed that you or your customers have on hand and completing the RGA form. Please note that the instructions on the form apply both to you and any customers to whom these products were delivered.

If you have any questions, please contact me at [REDACTED] or [REDACTED], or [REDACTED]@surgemedical.com.

Thank you for your cooperation.

Sincerely,

SURGE MEDICAL SOLUTIONS, LLC

[REDACTED]
[REDACTED]
[REDACTED]

V.P. - Operations

enclosure: RGA form