

TO WHOM IT MAY CONCERN

Your reference

Our reference

Tel. +49/56 61/71-

Fax +49/56 61/71-

Date

October 17, 2013

Intrafix Primeline / Intrafix Air P, Recall of several batches as a preventive measure**Intrafix Primeline, REF: 4062957****Batches: 3F07018901, 3F08018901, 3F10018901, 3F11018901, 3F12018901, 3F13018901, 3F14018901****Intrafix Air P, REF: 4062955 (Concerned batches only distributed in Germany)****Batches: 3E31018901, 3F01018901, 3F03018901, 3F04018901, 3F05018901, 3F06018901**

Dear Sirs,

In the course of internal quality checks we have found that with a low probability of occurrence these IV-sets may leak due to an assembly issue at the tubing, near to the patient connector. Therefore we have decided to recall the concerned batches as a preventive measure.

Our records show that your hospital has received **Intrafix Primeline / Intrafix Air P** sets from the specific batches.

Please do not use devices of these batches and send them back to us. A replacement has already been initiated.

Please forward this information to the users within your hospital immediately.

Please accept our apologies for any inconveniences.

We kindly ask you to confirm the receipt of this information. Should you have any questions, please call your local B. Braun contact person under Tel xxx-xxxx.

B | BRAUN

Page 2 to the letter of August 28, 2013 to
TO WHOM IT MAY CONCERN

Yours sincerely,