

# URGENT PRODUCT RECALL

May X, 2014

Dear Peritoneal Dialysis Provider:

Baxter Healthcare Corporation is voluntarily recalling the below reported lot of MiniCap Extended Life PD Transfer Set with Twist Clamp.

<b>Affected Products</b>	MiniCap Extended Life PD Transfer Set with Twist Clamp Product Code: R5C4482      Lot: H12F18054
<b>Problem Description</b>	Baxter has received complaints reporting connection issues between the PD transfer set catheter connector and an adapter. Baxter is currently investigating the issue.
<b>Hazard Involved</b>	Connection issues between the PD transfer set catheter connector and an adapter may lead to leaks and contamination of the fluid pathway which may result in peritonitis.
<b>Action to be taken if product was purchased directly from Baxter (Section may be adapted locally)</b>	<p>Baxter is requesting that you take the following actions:</p> <ol style="list-style-type: none"><li>1. Locate and remove all affected product from your facility. The product code and lot number can be found on the individual product package or shipping carton.</li><li>2. If you can identify patients who are currently using the affected lot, Baxter recommends replacing the transfer sets for those patients.</li><li>3. Contact Baxter's Customer Service to arrange for return and credit. Customer Service can be reached at 1-888-229-0001. Please have your Baxter eight digit ship-to account number ready when calling. <b>To be adapted locally</b></li><li>4. If you received your product directly from Baxter, please complete the attached Customer Reply Form (Attachment 1), and return it to Baxter by either fax or scanned e-mail. <b>To be adapted locally</b></li><li>5. If you are a dealer, wholesaler, or distributor/reseller that distributed any product to other facilities, please conduct a recall with your end user customers in accordance with your customary procedures.</li></ol>

**Further Information and support** For general questions regarding this Urgent Product Recall, please contact The Center for One Baxter at \_\_\_\_\_.  
(To be adapted locally)

We apologize for any inconvenience this communication may cause you, your staff, and your peritoneal dialysis patients.

Local MOH has been notified of this action. (To be adapted locally)

Sincerely,

(To be adapted locally)

cc: Attachment 1: Customer Reply Form