

URGENT FIELD SAFETY NOTICE

Commercial Name: HandyVac™, Guedel Airways, Universal Connector, Exchange bags for UnoMeter™ 500 (see Attachment 2 for commercial name and full details)

Issue Date: 02nd October 2019
 Original Notice Revised Notice Revision No.: Rev. 1

Reference/REF No: See Attachment 2: Product List

FSCA Ref.: 2019-008-011MA & 013MA

Type of Action: Recall/Product Disposal

Please note that this action only applies to specific product codes of the above-mentioned products.

| | |
|-----------------------------------|---|
| DESCRIPTION OF THE PROBLEM | <p>ConvaTec has voluntarily initiated a recall of specific product codes of the above-mentioned products.</p> <p>Internal assessment of product packaging integrity has confirmed that these devices are not meeting our expectations or those of our customers. Transportation testing conducted on the product in its current market packaging failed, confirming the potential for a breach in the sterile barrier.</p> <p>ConvaTec has not received any reports of incidents related to the packaging seal issue.</p> <p>Using a non-sterile device on a patient may expose the patient to infectious agents increasing the patient risk of developing infection.</p> |
|-----------------------------------|---|

| | |
|--|---|
| INFORMATION ON AFFECTED DEVICES | <p>HandyVac™ HandyVac™ is a single-use sterile low vacuum wound drainage system for post-operative surgical wound drainage, to drain wound fluids and prevent formation of hematomas, which can delay healing, cause pain or support infection. HandyVac™ consists of a drainage tube connected to a low-vacuum bellows and a drainage bag to collect the drained fluids. It is supplied as a standalone system or together with wound drains or wound drains that are pre-connected to a trocar.</p> <p>Guedel Airways Single-use sterile oropharyngeal airways to maintain a patent airway in unconscious or anaesthetized patients by keeping the airway passage free.</p> <p>Universal Connector Sterile connector with CH 04-18 Christmas Tree connector on one side and Luer slip port on the other side.</p> <p>Exchange Bags for UnoMeter™ 500 Sterile exchange bags for use with the UnoMeter™ 500 hourly diuresis monitoring system. The system is used to monitor urine output in critically ill patients. The affected exchange bags are used to collect larger amounts of urine over a prolonged period. Once the bag is full, it is usually replaced with a new bag if the system is not replaced as a whole.</p> |
|--|---|

| | |
|---|--|
| PRODUCT IDENTIFICATION PROCEDURE | <p><u>Confirmation of Specific Product Codes and Lot Numbers:</u></p> <p>This issue is limited to specific product codes as identified in Attachment 2.</p> <p>Only the identified product part codes within this notice may have a potential breach in the sterile barrier packaging.</p> <p>For this reason and to address any potential risk of harm, all the affected products should not be used.</p> <p>The only way to identify affected product is by comparing Product code/REF and LOT/Batch number to the recalled product list (see Attachment 1). There is no other discernable difference between affected and unaffected product.</p> <p>See Attachment 1 for examples of package labelling that highlights the location of the product code and LOT/Batch number on the device labels. These are located on the primary packaging and the shipping carton. The product code (reference number) is preceded by the word 'REF' and the LOT/Batch number is preceded by the symbol, LOT.</p> |
|---|--|

DISTRIBUTOR ACTIONS

| | |
|---|---|
| 1 | Verify if you have affected product in your inventory and identify which of your customers have received potentially affected product using Attachment 2: Product List . |
| 2 | <p>If you have distributed this product:</p> <ul style="list-style-type: none"> • To other Wholesalers: <ul style="list-style-type: none"> ○ Please forward this Field Safety Notice, including attachments, to them and ask that they follow these Distributor Actions. • To End Users <ul style="list-style-type: none"> ○ Please forward this Field Safety Notice, including attachments, to your end-users and ask that they follow the End User Actions Section (below). |
| 3 | <p>If you have received returned product from your customers that they have identified or if you have any stock on hand which you have verified as included in the scope of this Field Safety Notice (as described in Attachment 2: Product List).</p> <ul style="list-style-type: none"> • Immediately stop distributing and quarantine all affected products. • Contact ConvaTec Customer Service for further instructions on how to return affected product and to arrange credit. <p>Products will not be automatically replaced. For replacement product(s), please place new orders with ConvaTec using the usual process. ConvaTec will supply new products as soon as possible.</p> |
| 4 | <p>Please complete the Field Safety Notice <u>Distributor Response Form</u> and return to your Customer Service as soon as possible. Include details of product to be returned, if applicable.</p> <p>Please send a complete list of all consignees to your Customer Service as soon as possible. Please indicate if you have sent this Field Safety Notice, including attachments, to your consignees. This information is not required but important to allow ConvaTec to perform effectiveness checks.</p> <p>Please contact your Customer Service if you have any additional questions.</p> |

PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

RETAILER ACTIONS

| | |
|---|--|
| 1 | <p>For Retailers with a physical location (e.g. a “brick and mortar” store)</p> <ul style="list-style-type: none"> • Immediately stop distributing and quarantine all affected products • Please post the Field Safety Notice including Attachment 2 in a conspicuous location in your store. <p>For Retailers, which have a digital location (e.g. an eRetailer)</p> <ul style="list-style-type: none"> • Please forward Attachment 2: Product List showing affected Product Codes/REF No. and LOT Numbers via physical mail or e-mail to End-users. |
| 2 | <p>If you have received returned product from your customers that they have identified or if you have any stock on hand which you have verified as included in the scope of this Field Safety Notice (as described in Attachment 2: Product List, contact Customer Services Centre for further instructions on how to return affected product and to arrange credit.</p> <p>Products will not be automatically replaced. For replacement product (s), please place new orders with ConvaTec using the usual process. ConvaTec will supply new products as soon as possible.</p> |
| 3 | <p>Please complete the Field Safety Notice Distributor Response Form and return to your Customer Services as soon as possible. Include details of product to be returned, if applicable.</p> <p>Please send a complete list of all consignees to your Customer Service as soon as possible. Please indicate if you have sent this Field Safety Notice, including attachments, to your consignees.</p> <p>This information is not required but important to allow ConvaTec to perform effectiveness checks. Please contact your Customer Service if you have any additional questions.</p> |

PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

END USER ACTIONS (CLINICIANS, HOSPITALS, PATIENTS, OTHERS)

| | |
|---|--|
| 1 | <p>Our records show that you have taken delivery of affected product. Please stop the use of all affected devices as defined in this document.</p> |
| 2 | <p>This issue is limited to specific product codes and lots. Check stock and ensure that all affected devices that you have in stock are quarantined.</p> <p>If you do not know where the product code and/or lot number are located, please refer to Attachment 1: Examples of Labelling.</p> <p>To verify whether your lot is associated with the product included in the scope of this Field Safety Notice use Attachment 2: Product List showing affected Product Codes/REF No. and LOT Numbers.</p> |
| 3 | <p>Please complete the enclosed, Field Safety Notice ‘End User Response Form’ and forward it to your distributor as soon as possible.</p> <p>Include details of product to be returned, if applicable.</p> <p>This information is required to allow ConvaTec to perform effectiveness checks.</p> |
| 4 | <p>Contact your distributor (or ConvaTec Customer Service if you have purchased the product directly from ConvaTec) to arrange return of affected products, if applicable, and to arrange credit. For replacement products, please place new orders with your distributor or ConvaTec using the usual route and procedure.</p> <p>Please contact your Customer Service if you have any additional questions.</p> |

PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

| | |
|---|---|
| TRANSMISSION OF THIS FIELD SAFETY NOTICE | <p>This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p> |
|---|---|

ConvaTec is committed to providing quality products and services to our customers and we sincerely apologise for any inconvenience this notice may cause.

The relevant National Authorities have been advised about this Field Safety Corrective Action.

AUTHORISATION

| | | | |
|--------------------------------|----------------------------|---|--|
| <u>Name</u> [REDACTED] | <u>Title</u> [REDACTED] | <u>Address</u> ConvaTec Limited, First Avenue, Deeside Industrial Park, Deeside, CH5 2NU, U.K. | |
| <u>Date</u> 02 October 2019 | | <u>Signature</u> [REDACTED] | |



Regional ConvaTec Customer Service Contact:

Austria, Belgium, Czech Republic, Estonia, France, Germany, Ireland, Israel, Netherlands, Poland, Spain, Switzerland & Turkey

Tel: +41 (0) 52 630 54 01

Email: ccc.customerservice@convatec.com

Denmark

Tel: +45 48 167 030

Email: customerservicenordic@convatec.com

Finland

Tel: +358 20 7659 630

Email: mail.fi@convatec.com

Italy

Tel: +39 06 89789928

Email: CustomerServiceItaly@Convatec.com

Norway

Tel: +47 22 686 095

Email: customerservicenordic@convatec.com

Sweden

Tel: +46 (0)42 332010

Email: customerservicenordic@convatec.com

United Kingdom

Tel: +44 (0) 1244 832206

Email: Unomedical-uk.customerservice@convatec.com



FIELD SAFETY NOTICE DISTRIBUTOR RESPONSE FORM
PLEASE COMPLETE AND RETURN by Fax/Email

Consignee of the device:

| | |
|------------------------------|--|
| Consignee Account No: | |
| Consignee Name: | |
| Consignee Address: | |

The following products, HandyVac™, Guedel Airways, Universal Connector and Exchange bags For UnoMeter™ 500 have been distributed to your facility:

| Invoice # | Sales Order # | Product Code / REF No. | SAP Code | LOT No. | Quantity Delivered (pieces) |
|-----------|---------------|------------------------|----------|---------|-----------------------------|
| | | | | | |

| Distributors (Tick all that apply and give details, where applicable) | | |
|---|--|------------------------|
| <input type="checkbox"/> | I confirm the receipt, the reading and understanding of the Field Safety Notice. | |
| <input type="checkbox"/> | I have checked my stock and quarantined inventory | Add details to Table 1 |
| <input type="checkbox"/> | I have identified customers that received or may have received this device | |
| <input type="checkbox"/> | I have attached customer list | Add details to Table 2 |
| <input type="checkbox"/> | I have informed the identified customers of this FSN | Date sent: |
| <input type="checkbox"/> | I have received confirmation of reply from all identified customers | Attach responses |
| <input type="checkbox"/> | Neither I nor any of my customers has any affected devices in inventory | |

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN. Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.



Table 1. Quarantined Inventory: Record quantity (pieces) for each LOT to be returned.

| LOT No. | Units on Hand | LOT No. | Units on Hand | LOT No. | Units on Hand | LOT No. | Units on Hand |
|---------|---------------|---------|---------------|---------|---------------|---------|---------------|
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |

Table 2. Customer List: Please provide details of affected HandyVac™, Guedel Airways, Universal Connector and Exchange bags For UnoMeter™ 500 products that were distributed to your customers.

| Customer Name | Product Code / REF No. | SAP Code | LOT No. | Quantity (pieces) |
|---------------|------------------------|----------|---------|-------------------|
| | | | | |
| | | | | |
| | | | | |
| | | | | |

| FORM Completed and Returned From: | |
|-----------------------------------|--|
| Name (CAPITAL LETTERS): | |
| Position: | |
| Company Name: | |
| Address: | |
| Phone No: | |
| Signature: | |
| Date (dd/mmm/yyyy): | |



FIELD SAFETY NOTICE CUSTOMER RESPONSE FORM
PLEASE COMPLETE AND RETURN by Fax/Email

Consignee of the device:

| | |
|------------------------------|--|
| Consignee Account No: | |
| Consignee Name: | |
| Consignee Address: | |

The following products, HandyVac™, Guedel Airways, Universal Connector and Exchange bags For UnoMeter™ 500 have been distributed to your facility:

| Invoice # | Sales Order # | Product Code / REF No. | SAP Code | LOT No. | Quantity Delivered (pieces) |
|-----------|---------------|------------------------|----------|---------|-----------------------------|
| | | | | | |

| Customer action undertaken on behalf of Healthcare Organisation (Tick all that apply) | | |
|---|---|------------------------|
| <input type="checkbox"/> | I confirm receipt of the Field Safety Notice and that I read and understand its content. | |
| <input type="checkbox"/> | I performed all actions requested by the FSN. | |
| <input type="checkbox"/> | The information and required actions have been brought to the attention of all relevant users and executed. | |
| <input type="checkbox"/> | I have checked my stock and quarantined inventory | Add details to Table 1 |
| <input type="checkbox"/> | No affected devices are available for return | |

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN. Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.

Table 1. Quarantined Inventory: Record quantity (pieces) for each LOT to be returned.

| LOT No. | Units on Hand | LOT No. | Units on Hand | LOT No. | Units on Hand | LOT No. | Units on Hand |
|---------|---------------|---------|---------------|---------|---------------|---------|---------------|
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |











| | |
|------------------------------------|--|
| FORM Completed and Returned From: | |
| Name (CAPITAL LETTERS): | |
| Position: | |
| Company Name: | |
| Address: | |
| Phone No: | |
| Signature: | |
| Date (dd/mmm/yyyy): | |

List of Attachments

Attachment 1: Examples of Labelling and Affected Product

Attachment 2: EU Product List showing affected Product Codes/REF No. and LOT Numbers

Attachment 1: Examples of Labelling and Affected Product

| For the AFFECTED product | |
|---|---|
| <p>...the shipper labelling will include the product code and lot number.</p> <p>See example labels to the right.</p> <p>Product Code (REF/SKU)</p> <p>Lot Number</p> | <div style="border: 2px solid yellow; padding: 5px; display: inline-block; margin-bottom: 10px;"> REF 84061181 </div> <p> GB Unometer bag 2000 S Unometerpåse DK Unometerpose 2000 FIN Unometerpussi NO Unometerpose 2000 Made in Belarus </p> <p style="text-align: right; font-size: 24px;">2 X 100</p> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="text-align: center;">  </div> <div style="text-align: center;">  </div> <div style="font-size: 8px;"> ConvaTec Limited First Avenue Deeside Industrial Park Deeside, Flintshire, CH5 2NU, UK </div> </div> <div style="display: flex; justify-content: space-between; align-items: center; margin-top: 10px;"> <div style="text-align: center;">  2012-04  2017-03 </div> <div style="border: 2px solid blue; padding: 5px; display: inline-block;"> LOT 632574 </div> </div> <p style="text-align: center; font-size: 8px;">EAN No. 05705243208804</p> <div style="text-align: center;">  <small>(01)05705243208804(17)170331(10)632574 v.04</small> </div> |
| <p>...the market unit labelling will include the product code and lot number.</p> | <div style="border: 2px solid yellow; padding: 5px; display: inline-block; margin-bottom: 10px;"> 100 REF 84061181 </div> <p> GB Unometer bag 2000 S Unometerpåse DK Unometerpose 2000 FIN Unometerpussi NO Unometerpose 2000 </p> <div style="display: flex; justify-content: space-between; align-items: center; margin-top: 10px;"> <div style="border: 1px solid black; padding: 2px;"> EC REP UNOMEDICAL A/S Aaholmvej 1-3, Østød 4320 LEJRE DENMARK </div> <div style="text-align: center;">  </div> <div style="font-size: 8px;"> Conva Tec Limited First Avenue Deeside Industrial Park Deeside, Flintshire, CH5 2NU, UK </div> </div> <div style="display: flex; justify-content: space-between; align-items: center; margin-top: 10px;"> <div style="text-align: center;"> <small>EAN No. 05705243208798</small>  <small>(01)05705243208798(17)240731(10)632574</small> </div> <div style="text-align: center;">  2019-08  2024-07 </div> <div style="text-align: center;">  </div> </div> <p style="text-align: right; font-size: 8px;">v.05</p> |

Attachment 2: Product List

PURPOSE

This attachment is a listing of EU products included in the scope of this Field Safety Notice, which has been organized by Product Code. Within each product code, the lot numbers are organized numerically, smallest to largest.

Please contact your ConvaTec Customer Service Team Member for assistance.

| Column Header | Details | What Level of Packaging includes this information? |
|---------------------|--|---|
| Product Code | The Product Code is also known as the Catalog Number, or the SKU. | This number is included on the Shipper Label and the Market Unit Label. |
| Product Description | <p>The Product Description is a description of the finished product from the internal ConvaTec system, SAP. It is a brief, coded description of the product.</p> <p>Utilizing the first line item example, Product Code 84062182, Batch: 156697, UNOMETER BAG 3L BO (50/100) INT:</p> <ul style="list-style-type: none"> • REF: Product code • Batch: Lot number • UNOMETER BAG: UnoMeter™ 500 • BO(50/100): with Bottom Outlet (pack of 50/100) | A description of the material is on the Shipper Label and on the Market Unit Labeling; however, the description on these labels may not exactly match (verbatim) the Material Description as listed in this attachment. |
| Lot Number | A traceable number which is specific to the product. | This number is included on the Shipper Label and the Market Unit Label. |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 167273 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 169311 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS | 170793 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS | 172068 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 172069 |
| 40097181 | HANDYVAC M/700 ML. POSE | 172070 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 172071 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 172072 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 172072 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 172073 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 172074 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 172075 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 172075 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 172076 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 172076 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 172077 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 172078 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 172489 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 173172 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 173173 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS | 174559 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 174560 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 174561 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS | 175474 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT | 175474 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 175475 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 176275 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 176276 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14 | 179215 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 179216 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 179217 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT | 179218 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)AUS | 179219 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 179219 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 179220 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 179221 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 179222 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 179222 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 179232 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 179233 |
| 22022181 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 179233 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 179577 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 179583 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 179591 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT | 179924 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 179925 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 179926 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 179927 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 179928 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 179928 |
| 40234182 | UNO HANDYVAC CH12 MTL TRCR WOBO(1/20)INT | 179929 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 179930 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 179931 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 179932 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 179933 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 179949 |
| 22034022 | UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA | 179950 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 179955 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 179960 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 181110 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 181129 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS | 181132 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 181140 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 181141 |
| 22019182 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 181150 |
| 22021182 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 181151 |
| 22023182 | UNO GUEDEL AIRWAY SIZE 0 (100/300) INT | 181152 |
| 22024182 | UNO GUEDEL AIRWAY SIZE 00 (100/300) INT | 181153 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 181170 |
| 84010181 | KATETERANSATS, STERIL | 181173 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 181176 |
| 40105182 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 181330 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 181895 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 181925 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 181926 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 182121 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 182770 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT | 182774 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)AUS | 182775 |
| 40097181 | HANDYVAC M/700 ML. POSE | 182782 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 183061 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 183089 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 183424 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT | 183425 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 183578 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 183578 |
| 40097181 | HANDYVAC M/700 ML. POSE | 183595 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 183596 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 183597 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 183598 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 183598 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 183606 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 183609 |
| 22054023 | UNO GUEDEL AIRWAY SIZE 5 (1X50PK) IT | 183617 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 183836 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 183837 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 183865 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 183866 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 183871 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14 | 184213 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 184214 |
| 40097181 | HANDYVAC M/700 ML. POSE | 184230 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 184231 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 184240 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 184241 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 184256 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 184257 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 184257 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 184258 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 184259 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 184259 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 184263 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 184273 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 184292 |
| 84010181 | KATETERANSATS, STERIL | 184300 |
| 22054181 | TUNGEHOLDER, STR. 5, LYSEBLAA, STERIL | 184301 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 22054181 | NIELUPUTKI 5 st (50 kpl/ltk) | 184301 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 184770 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 184771 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 185141 |
| 40234182 | UNO HANDYVAC CH12 MTL TRCR WOBO(1/20)INT | 185142 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 185870 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 186406 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 186412 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 186413 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 187037 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 187039 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 187521 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 187521 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 187522 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 187522 |
| 40223183 | HANDYVAC + trocarCH14 + BOLSA cerrada 20ud | 187531 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 187540 |
| 22019182 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 187567 |
| 22021182 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 187569 |
| 22023182 | UNO GUEDEL AIRWAY SIZE 0 (100/300) INT | 187570 |
| 22024182 | UNO GUEDEL AIRWAY SIZE 00 (100/300) INT | 187571 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 187572 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 187578 |
| 84010181 | KATETERANSATS, STERIL | 187580 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 187589 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 187980 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 188289 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 188302 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 188303 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 188593 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 188599 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 188917 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 188925 |
| 22020182 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 188929 |
| 22021182 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 188930 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 189330 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 189331 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 189332 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 189332 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40125181 | HandyVac m/draen og trocar CH 14 | 189879 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 189879 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 189883 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 190075 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 190417 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 190421 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 190815 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT | 190822 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 190830 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 190830 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 190831 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 190832 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 190832 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 190840 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 191331 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 191361 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 191363 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 191368 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 191760 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 192121 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 192122 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 192502 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 192529 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 192529 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 192529 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 192530 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 192565 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 192853 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 192856 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 192856 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 193310 |
| 40227182 | HANDYVAC U/BUNDAFTAPNING CH12 | 193311 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 193311 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 193337 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 193337 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 193370 |
| 84010181 | KATETERANSATS, STERIL | 193381 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 193833 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 193837 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 193838 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 194258 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 194265 |
| 40224182 | UNO HANDYVAC CH18 MTL TRCR WOBO(1/20)AUS | 194272 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 194275 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 194294 |
| 22034022 | UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA | 194295 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 194302 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 194307 |
| 40097181 | HANDYVAC M/700 ML. POSE | 195031 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 195311 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 195452 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 195659 |
| 22024181 | NIELUPUTKI 00 st (100 kpl/ltk) | 195667 |
| 22024181 | UNO GUEDEL AIRWAYS SZ 00 (100/300) INT | 195667 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 195671 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 196102 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 196103 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 196131 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 196135 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 196135 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 196136 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 196538 |
| 22054181 | NIELUPUTKI 5 st (50 kpl/ltk) | 196546 |
| 22054181 | UNO GUEDEL AIRWAY SIZE 5 (50/200) INT | 196546 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 196869 |
| 22019184 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 196871 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 196872 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 196873 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 197370 |
| 40097181 | HANDYVAC M/700 ML. POSE | 197385 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 197386 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 197391 |
| 40107183 | UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT | 197400 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 197428 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 197429 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 197431 |
| 84010181 | KATETERANSATS, STERIL | 197432 |
| 22054181 | NIELUPUTKI 5 st (50 kpl/ltk) | 197435 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|---|-------------------|
| Product Code | Description | Lot Number |
| 22054181 | UNO GUEDEL AIRWAY SIZE 5 (50/200) INT | 197435 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 197436 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 197437 |
| 22037182 | UNO GUEDEL AIRWAY SIZE 000 (100/300) INT | 197444 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 197686 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 198011 |
| 84010181 | KATETERANSATS, STERIL | 198044 |
| 40097181 | HANDYVAC M/700 ML. POSE | 198991 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 198992 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 198992 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 199003 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 199018 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 199019 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál | 199027 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 199365 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 199368 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 199369 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 199371 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 199572 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 199726 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 199824 |
| 22023184 | UNO GUEDEL AIRWAY SIZE 0 (100/300) INT | 200469 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 200470 |
| 40097181 | HANDYVAC M/700 ML. POSE | 201050 |
| 40103181 | HANDYVAC M/DRAEN, CH 12/700 ML POSE | 201052 |
| 40103181 | HANDYVAC CH 12 (20/ltk) | 201052 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 201053 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 201054 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 201054 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 201055 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 201055 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 201075 |
| 40160181 | HANDYVAC CH12, TROAKAARI (20 KPL) | 201075 |
| 40160181 | UNO HANDYVAC CH12 MTL TROC BO(1/20) INT | 201075 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 201093 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 201113 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 201113 |
| 40223183 | HANDYVAC + trocar CH14 + BOLSA cerrada 20ud | 201119 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 201126 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 201127 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 201129 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 201133 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 201134 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál | 201148 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 201149 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCER CH14 | 201955 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 201989 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 201990 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 202031 |
| 22019182 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 202050 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 202055 |
| 84010181 | KATETERANSATS, STERIL | 202061 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 202306 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 202466 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 202470 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 202791 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 202792 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 202793 |
| 40228183 | HANDYVAC + trocCH16 + BOLSA cerrada 20ud | 202798 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 202799 |
| 40227183 | HANDYVAC + trocCH12 + BOLSA cerrada 20ud | 202800 |
| 40227183 | UNO HANDYVAC CH12 MTL TROC (1/20) INT | 202800 |
| 40105182 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 203040 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 203056 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 203491 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 203492 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál | 203503 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 203504 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 203711 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 204020 |
| 22019184 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 204024 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 204025 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 204026 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 204029 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 204320 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 204335 |
| 40233182 | UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT | 204342 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 204349 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 204351 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 204353 |
| 22054023 | UNO GUEDEL AIRWAY SIZE 5 (50/200) IT | 204357 |
| 22054023 | UNO GUEDEL AIRWAY SIZE 5 (1X50PK) IT | 204357 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 204358 |
| 84010181 | KATETERANSATS, STERIL | 204359 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál | 204361 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 204362 |
| 22021181 | NIELUPUTKI 3 st (6x50 kpl/ltk) | 204367 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 205219 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 205242 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 205242 |
| 22019184 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 205691 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 206057 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 206072 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 206073 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 206074 |
| 22019184 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 206082 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 206083 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 206472 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 206476 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 206477 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 206478 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 206968 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 206973 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 207363 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 207375 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 207375 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 207375 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 207376 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 207376 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 207376 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 207377 |
| 22022182 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 207426 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 207429 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 207430 |
| 22023184 | UNO GUEDEL AIRWAY SIZE 0 (100/300) INT | 207431 |
| 22055181 | NIELUPUTKI 6 st (50 kpl/ltk) | 207432 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál | 207433 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 207434 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 207674 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 207685 |
| 22024184 | UNO GUEDEL AIRWAY SIZE 00 (100/300) INT | 207686 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 207691 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 208231 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 208237 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 208249 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 208604 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 208605 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 208842 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 208852 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 208864 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 208877 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 208878 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 208878 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 208878 |
| 22021181 | TUNGEHOLDER, STR. 3, GUL, STERIL | 209032 |
| 22021181 | NIELUPUTKI 3 st (6x50 kpl/ltk) | 209032 |
| 22021181 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 209032 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 209033 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 209316 |
| 22019181 | TUNGEHOLDER, STR. 1, HVID, STERIL | 210157 |
| 22019181 | NIELUPUTKI 1 st (100 kpl/ltk) | 210157 |
| 22019181 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 210157 |
| 22021181 | NIELUPUTKI 3 st (6x50 kpl/ltk) | 210158 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 210159 |
| 40097181 | HANDYVAC M/700 ML. POSE | 210488 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 210489 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 210489 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 210490 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 210508 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 210515 |
| 22019181 | NIELUPUTKI 1 st (100 kpl/ltk) | 210525 |
| 22021181 | NIELUPUTKI 3 st (6x50 kpl/ltk) | 210526 |
| 84010181 | KATETERANSATS, STERIL | 210527 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál | 210529 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 210529 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 210530 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 210950 |
| 84021182 | UNO URIMETER BAG 2L (75/300) IT | 210953 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 210957 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 211094 |
| 40107183 | UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT | 211095 |
| 22024184 | UNO GUEDEL AIRWAY SIZE 00 (100/300) INT | 211116 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 211197 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCER CH14 | 211357 |
| 40097181 | HANDYVAC M/700 ML. POSE | 211375 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 211376 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 211376 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 211386 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 211414 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 211414 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 211415 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 211415 |
| 22031022 | UNO GUEDEL AIRWAY SIZE 1 (100/300) ITA | 211434 |
| 22022181 | NIELUPUTKI 4 st (4x50 kpl/ltk) | 211439 |
| 22022181 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 211439 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 211653 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 211654 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 211655 |
| 22024184 | UNO GUEDEL AIRWAY SIZE 00 (100/300) INT | 211656 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 211659 |
| 22034022 | UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA | 211917 |
| 22023181 | NIELUPUTKI 0 st (100 kpl/ltk) | 212079 |
| 22023181 | UNO GUEDEL AIRWAY SIZE 0 (100/300) INT | 212079 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 212409 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 212417 |
| 40171183 | UNO HANDYVAC CH16 BAG W/TAP (1X1PK) IT | 212929 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 212939 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 212939 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 212939 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 212940 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 212940 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 212940 |
| 40223183 | HANDYVAC + trocarCH14 + BOLSA cerrada 20ud | 212944 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC (1/20) IT | 213036 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 213036 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40125181 | HandyVac m/draen og trocar CH 14 | 213050 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 213050 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 213305 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 213307 |
| 40097181 | HANDYVAC M/700 ML. POSE | 213804 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 213805 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 213809 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál | 213854 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 213855 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 213856 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 213857 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 214202 |
| 40097181 | HANDYVAC M/700 ML. POSE | 214593 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 214594 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 214595 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 214595 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 214600 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 214609 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 214619 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 215190 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 215190 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 215199 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 215203 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 215426 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 215435 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 215435 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 216421 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 216428 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 216428 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 216428 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 216429 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 216432 |
| 22034022 | UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA | 216452 |
| 22024184 | UNO GUEDEL AIRWAY SIZE 00 (100/300) INT | 216457 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 216864 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 217375 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 217386 |
| 40107183 | UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT | 217387 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 217661 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 217681 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT | 217723 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 218119 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 218128 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 218128 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 218129 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 218130 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 218131 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 218131 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 218132 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 218133 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 218429 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 218441 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 218482 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 218636 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 218641 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 218645 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 218645 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 218646 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 218653 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 219139 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 219167 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 219168 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 219169 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 219175 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 219651 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 219664 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 219665 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 219665 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 219690 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 219690 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 219691 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 219695 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 219696 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 219698 |
| 22020181 | NIELUPUTKI 2 st (6x50 kpl/ltk) | 219725 |
| 22021181 | NIELUPUTKI 3 st (6x50 kpl/ltk) | 219731 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 219733 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 219734 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 219982 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 220356 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 220357 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 220362 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 220374 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 220377 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 220798 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 221105 |
| 40097181 | HANDYVAC M/700 ML. POSE | 221123 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 221124 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 221124 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 221124 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 221125 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 221125 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 221126 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 221126 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 221127 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 221144 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 221502 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 221855 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT | 222210 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 222548 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 222549 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 222549 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 222549 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 222554 |
| 40097182 | UNO HANDYVAC BAG 700ML BO(1/20) INT | 222555 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 222569 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 223012 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 223013 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 223014 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 223018 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 223452 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 223452 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 223453 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 223453 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 223458 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 223459 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 223460 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 223470 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 223479 |
| 40234182 | UNO HANDYVAC CH12 MTL TRCR WOBO(1/20)INT | 223482 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 223482 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 223490 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 223493 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 223496 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT | 224343 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 224343 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 224364 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 224364 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 224365 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 224365 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 224374 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 224400 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 224400 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 224401 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 224401 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 224404 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 224408 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 224415 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 224468 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 225138 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 225144 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 225319 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 225323 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 225721 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 225722 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 226131 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 226132 |
| 40107183 | UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT | 226133 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 226444 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 226446 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 226447 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 226589 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT | 226589 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 226876 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 226878 |
| 40097181 | HANDYVAC M/700 ML. POSE | 226893 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40127181 | HandyVac m/draen og trocar CH 18 | 226895 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 226895 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 226922 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 226923 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 226923 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 226928 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 226929 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 226932 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 226938 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 226944 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 226947 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 227255 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 227257 |
| 40097181 | HANDYVAC M/700 ML. POSE | 227841 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 227842 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 227842 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 227843 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 227843 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 227844 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 227851 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 227883 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 228299 |
| 40105182 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 228300 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 228333 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 228335 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 228336 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 228337 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 228787 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 228788 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 228791 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 228822 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 228824 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 229734 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 229734 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 229735 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 229735 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 229739 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 229741 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 229751 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 229767 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 229769 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 230224 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 230760 |
| 40097181 | HANDYVAC M/700 ML. POSE | 230776 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 230777 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 230777 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 230778 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 230778 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 230779 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 230787 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 231155 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 231633 |
| 40097181 | HANDYVAC M/700 ML. POSE | 232455 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 232457 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 232467 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 232468 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 232468 |
| 22019184 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 232760 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 232761 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 232762 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 232763 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 232764 |
| 40097181 | HANDYVAC M/700 ML. POSE | 233067 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 233068 |
| 22035022 | UNO GUEDEL AIRWAY SIZE 0 (100/300) ITA | 233096 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 233098 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 233102 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 233717 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 233721 |
| 22019184 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 233722 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 233723 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 233724 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 233725 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 233730 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 233733 |
| 22021183 | UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT | 233740 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 233745 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 233766 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|---|-------------------|
| Product Code | Description | Lot Number |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 233766 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 234147 |
| 40107183 | UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT | 234155 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 234160 |
| 40097181 | HANDYVAC M/700 ML. POSE | 234909 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 234910 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 234911 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 234911 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 234912 |
| 84061181 | UNOMETER BAG 2L (100/200) INT | 234926 |
| 22020182 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 234929 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 235130 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 235130 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 235326 |
| 22019184 | UNO GUEDEL AIRWAY SIZE 1 (100/300) INT | 235328 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 235329 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 235330 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 235331 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 235332 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT | 235336 |
| 22021183 | UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT | 235337 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 235729 |
| 40097181 | HANDYVAC M/700 ML. POSE | 236643 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 236644 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 236644 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 236644 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 236645 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 236646 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 236646 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 236652 |
| 40223183 | HANDYVAC + trocar CH14 + BOLSA cerrada 20ud | 236654 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 236676 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 237467 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 237474 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 237476 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 237477 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 300309 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 300333 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 300339 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 300339 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 300353 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 300355 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 300412 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 300478 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 300479 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 300480 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 300542 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 300588 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 300590 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 300602 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCER CH14 | 300608 |
| 40097181 | HANDYVAC M/700 ML. POSE | 300625 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 300626 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 300627 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 300627 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 300634 |
| 84010181 | KATETERANSATS, STERIL | 300642 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 300665 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 300678 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 300678 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 300678 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 300682 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 300686 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 300692 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 300693 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 300705 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 300705 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 300706 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 300708 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 300708 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 300808 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 300837 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 300837 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 300839 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 300879 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 300892 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 300892 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 300893 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 300893 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 300900 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 300912 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 300921 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 300922 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 300924 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 300926 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 300931 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 300934 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 301369 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 301393 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 301393 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 301394 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 301394 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 301395 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 301395 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 301396 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 301400 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 301401 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 301402 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 301403 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 301405 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 301407 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 301460 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 301489 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 301490 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 301497 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 301498 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 301551 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 301571 |
| 22021183 | UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT | 301594 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 301765 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 301766 |
| 84010181 | KATETERANSATS, STERIL | 301784 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 301785 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 301785 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 301786 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 301786 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 301787 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 301787 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 301841 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 302103 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 302125 |
| 40097181 | HANDYVAC M/700 ML. POSE | 302139 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 302140 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 302140 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 302140 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 302144 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 302234 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 302277 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 302286 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 302301 |
| 40227183 | UNO HANDYVAC CH12 MTL TROC (1/20) INT | 302448 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 302527 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 302528 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 302528 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 302533 |
| 22020182 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 302541 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 302557 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 302560 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 302560 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 302587 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 302760 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 302760 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 302786 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 302787 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 302797 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 302798 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 302809 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 302810 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 302855 |
| 22022183 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 302864 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 302945 |
| 40097181 | HANDYVAC M/700 ML. POSE | 302971 |
| 40103181 | HANDYVAC M/DRAEN, CH 12/700 ML POSE | 302972 |
| 40103181 | HANDYVAC CH 12 (20/ltk) | 302972 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 302973 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 302973 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40125181 | HandyVac m/draen og trocar CH 14 | 302974 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 302974 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 302975 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 302975 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 302989 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 302990 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 303009 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 303009 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 303009 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 303011 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 303012 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 303014 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 303176 |
| 40105182 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 303177 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 303178 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 303297 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 303321 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 303343 |
| 84021182 | UNO URIMETER BAG 2L (75/300) IT | 303382 |
| 84010183 | UNO CONN CH04-CH18(100/800) IT | 303383 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 303403 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 303407 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 303410 |
| 22034022 | UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA | 303414 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 303427 |
| 40227183 | HANDYVAC + trocCH12 + BOLSA cerrada 20ud | 303457 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 303471 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 303501 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 303530 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 303585 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 303585 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 303586 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 303586 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 303586 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 303683 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 303686 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 303692 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 303761 |
| 84010181 | KATETERANSATS, STERIL | 303769 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 22020182 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 303771 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 303772 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 303781 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 303784 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 303956 |
| 84010181 | KATETERANSATS, STERIL | 304094 |
| 40097181 | HANDYVAC M/700 ML. POSE | 304098 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 304099 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 304100 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 304100 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 304123 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 304123 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 304123 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 304184 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 304271 |
| 22055181 | NIELUPUTKI 6 st (50 kpl/ltk) | 304372 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 304391 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 304401 |
| 84010181 | KATETERANSATS, STERIL | 304442 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 304448 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 304449 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 304449 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 304449 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 304450 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 304451 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 304451 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 304452 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 304452 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 304465 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 304466 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 304468 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 304470 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 304471 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 304623 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 304702 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 304734 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 304737 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 304738 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 304741 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 304763 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 304798 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 305018 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 305035 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 305035 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 305035 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 305036 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 305037 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 305053 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 305054 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 305409 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 305409 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 305409 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 305410 |
| 40103181 | HANDYVAC M/DRAEN, CH 12/700 ML POSE | 305411 |
| 40103181 | HANDYVAC CH 12 (20/ltk) | 305411 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 305412 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 305412 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 305412 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 305413 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 305414 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 305415 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 305416 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 305416 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 305416 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 305417 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 305418 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 305418 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 305419 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 305425 |
| 84010181 | KATETERANSATS, STERIL | 305438 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 305441 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 305444 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 305446 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 305447 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 305737 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 305738 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 305738 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 305739 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 305739 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 305739 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 305740 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 305740 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 305741 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 305745 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14 | 305747 |
| 84010181 | KATETERANSATS, STERIL | 305756 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 305760 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 305766 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 306028 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 306059 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 306060 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 306061 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 306194 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 306195 |
| 84010181 | KATETERANSATS, STERIL | 306197 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 306221 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 306346 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 306346 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 306347 |
| 84010181 | KATETERANSATS, STERIL | 306358 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 306452 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 306452 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 306635 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 306635 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 306644 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 306645 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 306857 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 306857 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 306858 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 306858 |
| 40123181 | HANDYVAC CH10, TROAKAARI (20 KPL) | 306858 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 306859 |
| 40097181 | HANDYVAC M/700 ML. POSE | 306874 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 306875 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 306876 |
| 40097181 | HANDYVAC M/700 ML. POSE | 307109 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 307110 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 307110 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 307111 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 307112 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 307113 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 307113 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 307114 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 307115 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 307117 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 307131 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 307133 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 307134 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 307353 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 307353 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 307354 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 307354 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 307354 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 307355 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 307355 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 307355 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 307364 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 307366 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 307367 |
| 40105182 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 307379 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 307411 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 307412 |
| 40097181 | HANDYVAC M/700 ML. POSE | 307561 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 307562 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 307562 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 307563 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 307564 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 307565 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 307568 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 307569 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 307581 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 307582 |
| 40097181 | HANDYVAC M/700 ML. POSE | 307757 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 307758 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 307758 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 307759 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 307759 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 307759 |
| 84010181 | KATETERANSATS, STERIL | 307774 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 307778 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 307791 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 307792 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 307794 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 308155 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 308157 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 308169 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 308173 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 308174 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 308181 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 308181 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 308187 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 308309 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 308385 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 308385 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 308392 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 308446 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 308449 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 308451 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 308452 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 308452 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 308454 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 308461 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 308461 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 308626 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 308627 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 308628 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 308664 |
| 40097181 | HANDYVAC M/700 ML. POSE | 308807 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 308808 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 308809 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 308809 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 308810 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 308811 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 308814 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 308816 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 308816 |
| 84010181 | KATETERANSATS, STERIL | 308822 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 308825 |
| 22034022 | UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA | 308837 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 308846 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 308972 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 308972 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 308976 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 308979 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 309003 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 309003 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 309129 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 309181 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 309181 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 309182 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 309183 |
| 40223183 | HANDYVAC + trocarCH14 + BOLSA cerrada 20ud | 309185 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 309186 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 309193 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 309193 |
| 40097181 | HANDYVAC M/700 ML. POSE | 309226 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 309227 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 309227 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 309227 |
| 40103181 | HANDYVAC M/DRAEN, CH 12/700 ML POSE | 309228 |
| 40103181 | HANDYVAC CH 12 (20/ltk) | 309228 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 309229 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 309229 |
| 84010181 | KATETERANSATS, STERIL | 309236 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 309422 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 309423 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 309424 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 309425 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 309454 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 309563 |
| 45 951 1S-50 | UNO PLEURA PUNCTURE (1/20) INT | 309567 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 309599 |
| 40097181 | HANDYVAC M/700 ML. POSE | 309665 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 309666 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 309667 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 309667 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 309667 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 309668 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 309668 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 309669 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 309670 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 309670 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 309671 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 309671 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 309671 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 309674 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 309676 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 309676 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 309677 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 309680 |
| 84010181 | KATETERANSATS, STERIL | 309691 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 309694 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 309696 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 309697 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 309697 |
| 40097181 | HANDYVAC M/700 ML. POSE | 309885 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 309886 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 309887 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 309887 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 309887 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 309916 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 309917 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 309918 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 309919 |
| 84010181 | KATETERANSATS, STERIL | 309923 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 309932 |
| 22034022 | UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA | 309933 |
| 40097181 | HANDYVAC M/700 ML. POSE | 310091 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 310092 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 310092 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 310093 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 310094 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 310095 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 310096 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 310099 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 310100 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 310112 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 310117 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 310119 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 310120 |
| 40184182 | UNO HANDYVAC CH16 MTL TROCAR BO(1/20) IT | 310144 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 310391 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 310392 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 310392 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 310392 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 310403 |
| 40097181 | HANDYVAC M/700 ML. POSE | 310405 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 310406 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 310423 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 310423 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 310578 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 310597 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 310598 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 310598 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 310603 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 310605 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 310606 |
| 40097181 | HANDYVAC M/700 ML. POSE | 310786 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 310787 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 310787 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 310788 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 310788 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 310789 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 310789 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 310789 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 310790 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 310790 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 310793 |
| 84010181 | KATETERANSATS, STERIL | 310805 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 310824 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 310825 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 310976 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 310977 |
| 40097181 | HANDYVAC M/700 ML. POSE | 311012 |
| 40103181 | HANDYVAC CH 12 (20/ltk) | 311013 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 311014 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 311019 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 311022 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 311023 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 311032 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 311186 |
| 84021182 | UNO URIMETER BAG 2L (75/300) IT | 311190 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 311201 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 311326 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 311327 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 311367 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 311367 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 311367 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 311368 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 311368 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 311368 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 311369 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 311437 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 311437 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 311437 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 311571 |
| 22020184 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 311680 |
| 22021184 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 311681 |
| 22022184 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 311682 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 311699 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 311700 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 311702 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 311717 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 311717 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 311717 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 311718 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 311735 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 311737 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 312059 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 312059 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 312060 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 312062 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 312063 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 312063 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 312067 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 312069 |
| 40103181 | HANDYVAC CH 12 (20/ltk) | 312108 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 312205 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 312251 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 312251 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 312251 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 312360 |
| 22021183 | UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT | 312370 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 312386 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 312390 |
| 40171183 | UNO HANDYVAC CH16 BAG W/TAP (1X1PK) IT | 312475 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 312498 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 312498 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 312501 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT | 312629 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 312715 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 312715 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 312715 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 312716 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 312716 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 312717 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 312733 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 312870 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 312884 |
| 84010181 | KATETERANSATS, STERIL | 312891 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 313032 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 313032 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 313032 |
| 84010181 | KATETERANSATS, STERIL | 313038 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 313041 |
| 22021183 | UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT | 313054 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 313311 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 313343 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 313470 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 313470 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84010181 | KATETERANSATS, STERIL | 313479 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 313481 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 313482 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 313482 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 313500 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 313500 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 313798 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 313799 |
| 40097181 | HANDYVAC M/700 ML. POSE | 313825 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 313828 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 313828 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 313845 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 313845 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 313845 |
| 40103181 | HANDYVAC M/DRAEN, CH 12/700 ML POSE | 313846 |
| 40103181 | HANDYVAC CH 12 (20/ltk) | 313846 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 313847 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 313848 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 313849 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 313849 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 313849 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 314075 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 314075 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 314075 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 314077 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 314099 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 314100 |
| 22020182 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 314111 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT | 314112 |
| 22054182 | UNO GUEDEL AIRWAY SIZE 5 (50/200) INT | 314114 |
| 40097181 | HANDYVAC M/700 ML. POSE | 314271 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 314274 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 314274 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 314290 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 314290 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 314459 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 314583 |
| 22032022 | UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA | 314764 |
| 40099181 | HANDYVAC CH 8 (20/ltk) | 314863 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84010183 | UNO CONN CH04-CH18(100/800) IT | 314883 |
| 40105182 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 315047 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 315170 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 315170 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 315188 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 315244 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 315397 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 315541 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 315541 |
| 84010181 | KATETERANSATS, STERIL | 315554 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 315597 |
| 84021183 | UNOMETER 500 BOLSAS 1'5L. 75uds (sin val | 315672 |
| 84021183 | UNO URIMETER BAG 2L (75/300) INT | 315672 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 315682 |
| 84077183 | UNOMETER BAG 2L BO (75/150) INT | 315682 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT | 315745 |
| 22021183 | UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT | 315746 |
| 84021182 | UNO URIMETER BAG 2L (75/300) IT | 218437 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 316293 |
| 84021182 | UNO URIMETER BAG 2L (75/300) IT | 304154 |
| 84021182 | UNO URIMETER BAG 2L (75/300) IT | 308337 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 316306 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 316306 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 316306 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 316308 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 316442 |
| 22033022 | UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA | 316550 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 316569 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 316569 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 316569 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 316577 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 316769 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 316769 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 316769 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 316772 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 316773 |
| 40234182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT | 316775 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 316817 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 316818 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 316820 |
| 40160181 | UNO HANDYVAC CH12 MTL TROC BO(1/20) INT | 316820 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 316825 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 316825 |
| 40097181 | HANDYVAC M/700 ML. POSE | 316986 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 317016 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 317117 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 317117 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 317117 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 317118 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 317118 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 317118 |
| 40097181 | HANDYVAC M/700 ML. POSE | 317308 |
| 40223183 | HANDYVAC + trocarCH14 + BOLSA cerrada 20ud | 317504 |
| 84010181 | KATETERANSATS, STERIL | 317508 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (50/300) IT | 317716 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT | 317716 |
| 45 951 1S-50 | UNO PLEURA PUNCTURE (1/20) INT | 317791 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 317991 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 317992 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 317992 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 318069 |
| 40099181 | HANDYVAC CH 8 (20/ltk) | 318168 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 318456 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 318525 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 318525 |
| 22022183 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 318604 |
| 22021182 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 318758 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 318781 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 319365 |
| 22022183 | UNO GUEDEL AIRWAY SIZE 4 (50/200) INT | 319385 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 319449 |
| 84010181 | KATETERANSATS, STERIL | 319796 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 319816 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 319817 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 319818 |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 319818 |
| 40123181 | HANDYVAC CH10, TROAKAARI (20 KPL) | 319818 |
| 40123181 | HandyVac m/draen og trocar CH 10 | 319819 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40123181 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 319819 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 319820 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 319821 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 319822 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14 | 319825 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 319825 |
| 40223183 | HANDYVAC + trocCH14 + BOLSA cerrada 20ud | 319826 |
| 40227182 | HANDYVAC U/BUNDAFTAPNING CH12 | 319828 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 319828 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 319829 |
| 40101181 | HANDYVAC M/DRAEN, CH 10/700 ML POSE | 319830 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 319830 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 319830 |
| 40097181 | HANDYVAC M/700 ML. POSE | 319962 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 320024 |
| 40222183 | HANDYVAC + trocCH10 + BOLSA cerrada 20ud | 320284 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 320287 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT | 320337 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 320347 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 320437 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 320558 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 320572 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 320582 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 320679 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 320679 |
| 40103181 | HANDYVAC M/DRAEN, CH 12/700 ML POSE | 320728 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 320835 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 320846 |
| 22037182 | UNO GUEDEL AIRWAY SIZE 000 (100/300) INT | 320897 |
| 40105182 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 320994 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 320995 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 321027 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 321027 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 321027 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 321028 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 321163 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 321211 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 321212 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 321213 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 321214 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 321214 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 321214 |
| 40127182 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 321215 |
| 22020182 | UNO GUEDEL AIRWAY SIZE 2 (50/300) INT | 321407 |
| 22021182 | UNO GUEDEL AIRWAY SIZE 3 (50/300) INT | 321408 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 321487 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 321487 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 321487 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 321488 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 321488 |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 321505 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 321528 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 321528 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 321528 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 321640 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 321641 |
| 40235182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT | 321644 |
| 40097181 | HANDYVAC M/700 ML. POSE | 321690 |
| 22020183 | UNO GUEDEL AIRWAY SIZE 2 (50/300) IT | 321694 |
| 84010182 | UNO CONN CH04-CH18(100/800) INT | 321887 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 321924 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 321924 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 322228 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 322411 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 322418 |
| 40127181 | HandyVac m/draen og trocar CH 18 | 322428 |
| 40127181 | UNO HANDYVAC CH18 MTL TROC BO(1/20) INT | 322428 |
| 40127181 | HANDYVAC CH18, TROAKAARI (20 KPL) | 322428 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 322886 |
| 40227182 | HANDYVAC U/BUNDAFTAPNING CH12 | 322891 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 322891 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 322898 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 322917 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14 | 322941 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 322941 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 322972 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 323025 |
| 40097181 | HANDYVAC M/700 ML. POSE | 323288 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 84061182 | UNOMETER BAG 2L (100/200) IT | 323383 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 323552 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 323599 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 323599 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 323702 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 323702 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 323702 |
| 84077181 | UNOMETER BAG 3L BO (75/150) INT | 323705 |
| 40227183 | HANDYVAC + trocCH12 + BOLSA cerrada 20ud | 323710 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14 | 323893 |
| 40223182 | UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS | 323893 |
| 84062182 | UNOMETER BAG 3L BO (50/100) INT | 323897 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 324189 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 324192 |
| 40107182 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 324205 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 324206 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 324363 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 324363 |
| 40105181 | HANDYVAC CH 14 (20/ltk) | 324363 |
| 40125182 | UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR | 324373 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 324407 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 324459 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 324619 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 324619 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 324619 |
| 84010181 | KATETERANSATS, STERIL | 324627 |
| 84010181 | UNO CONN CH04-CH18(100/800) INT | 324627 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 324706 |
| 40107181 | HANDYVAC CH 18 (20/ltk) | 324737 |
| 84077183 | UNOMETER 500 BOLSAS valv. vac 2L. 75uds | 324857 |
| 40101181 | UNO HANDYVAC CH10 BAG W/TAP (1/20) IT | 325020 |
| 40101181 | HANDYVAC CH 10 (20/ltk) | 325020 |
| 84077182 | UNOMETER BAG 2L BO (75/150) IT | 325035 |
| 84077182 | UNO UNOMETER BAG 2L BO (75/150) IT | 325035 |
| 40103182 | UNO HANDYVAC CH12 BAG W/TAP(1/20) INT | 325052 |
| 40101182 | UNO HANDYVAC CH10 BAG W/TAP(1/20) INT | 325408 |
| 40125183 | UNO HANDYVAC CH14 MTL TROC BO(1/20) IT | 325409 |
| 40223182 | HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14 | 325410 |
| 40105181 | HANDYVAC M/DRAEN, CH 14/700 ML POSE | 325429 |

| ATTACHMENT 2: EU Product List | | |
|--------------------------------------|--|-------------------|
| Product Code | Description | Lot Number |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 325429 |
| 40121182 | UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT | 325818 |
| 40125181 | HandyVac m/draen og trocar CH 14 | 326205 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 326205 |
| 40125181 | HANDYVAC CH14, TROAKAARI (20 KPL) | 326205 |
| 40105181 | UNO HANDYVAC CH14 BAG W/TAP (1/20) INT | 326619 |
| 40123182 | UNO HANDYVAC CH10 METAL TROC BO(1/20) IT | 326778 |
| 40160181 | HANDYVAC M/BUNDAFTAPNING CH12 | 326780 |
| 40107181 | HANDYVAC M/DRAEN, CH 18/700 ML POSE | 327459 |
| 40107181 | UNO HANDYVAC CH18 BAG W/TAP(1/20) INT | 327459 |
| 40125181 | UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU | 327942 |
| 40227182 | UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS | 328440 |