

April 18, 2007

## Information on Potassium Iquirapid Cat.No.: 10118 LOT H8146

Dear Distributor,

During internal Quality Control we recognized that the component TPB Reagent of Potassium Iquirapid Cat. No. 10118, Lot H8146 may degrade.

Such degradation causes a color shift of the TPB Reagent from colorless-yellowish to dark brown. The problem was caused by a defective raw material. Due to the slow degradation process, we were not able to notice this problem during manufacture and end control. Even if the product may work, we recommend to stop selling kits of this lot H8146. If you already have delivered kits of this lot to customers, please inform them accordingly and instruct them to return or destroy the kits following local regulations.

Human will certainly give full replacement for all the kits of lot H8146 you received. We will deliver the replacement kits with your next shipment. If you should urgently need replacements, please inform us immediately about the required amounts.

We apologize for the inconvenience you have had in this issue and remain

With best regards,

Human GmbH

  
Customer Service Manager

Deutsche Bank  
BLZ 510 70021  
Kto.-Nr. 4221 974  
SWIFT-BIC: DEUT DE FF 510  
IBAN: DE60 5107 0021 0422 1974 00

Norddeutsche Landesbank  
BLZ 250 500 00  
Kto.-Nr. 122 0377 24  
SWIFT-BIC: NOLA DE 2H  
IBAN: DE76 2505 0000 0122 0377 24

Dresdner Bank  
BLZ 510 800 60  
Kto.-Nr. 202 783 500  
SWIFT-BIC: DRES DE FF 510  
IBAN: DE10 5108 0060 0202 7835 00

Geschäftsführer:  
Wolfgang Flemming  
Vorsitzender des Beirats:  
Dr. Thomas Guth  
Telefon +49 6122 9988-0  
Telefax +49 6122 9988 100  
e-Mail human@human.de  
ILN 40 33145 000006  
RG: Wiesbaden HRB 10782