

DADE BEHRING

Dade Behring Marburg GmbH
P. O. Box 11 49
35001 Marburg, Germany
www.dadebehring.com

Phone: +49 (6421) 39-6402
Fax: +49 (6421) 39-4594

To whom it may concern

Our Ref.
BR-06407

Ext.
06421-39-5402

Fax
06421-39-4594

Date
2007-08-14

**Field Correction
on
Dade® Actin® Activated Cephaloplastin Reagent (B4218-2),
Lot Number 537148**

Dear Valued Dade Behring Customer,

Our records indicate that you have received Dade® Actin® Activated Cephaloplastin Reagent, number 537148.

The purpose of this letter is to inform you to discontinue the use of and discard the affected lot as well as any remaining inventory.

Dade Behring has confirmed that control and patient sample results obtained using the above lot of Dade® Actin® Activated Cephaloplastin Reagent may be slightly shortened and controls may recover outside the assigned range. This may result in invalid runs.

If you have any questions or need assistance with obtaining replacement product, please contact your local Dade Behring representative.

Please distribute this information to all appropriate personnel in your laboratory, retain a copy in your files, and forward this information to all parties that may use this product, including others to whom you may have transferred this lot.

We apologize for any inconvenience this action may have caused you and we thank you for your continued support of Dade Behring Hemostasis products.

Sincerely,


Director
Quality Systems & Compliance


Marketing Manager
GPM Hemostasis