

January 02, 2008

Field Correction

Syva® Emit® II Plus Amphetamines Assay
Catalogue Number 9C039UL, Lot Z10

Dear Customer:

Our records indicate that your facility may have received the Syva® Emit® II Plus Amphetamines Assay, Catalogue number 9C039UL, lot Z10.

Investigation by Siemens Healthcare Diagnostics has shown that the Z10 lot of 9C039UL is exhibiting unacceptable calibration characteristics. You will observe problems with calibration, e.g. not able to meet the lot card guidelines for minimum separations or expected calibrator rates.

We recommend that you discard the aforementioned lot and contact your local Siemens Representative with your replacement needs of 9C039UL

Please forward this notice to anyone to whom you may have distributed these products. If you have any questions about the information in the attached letter, please contact the Technical Assistance Center or call your local Siemens Healthcare Diagnostics Representative

Siemens Healthcare Diagnostics regrets any inconvenience this situation may cause. Thank you for your patience and continued support.

Sincerely,


Product Manager, Syva®

As of November 6, 2007, Dade Behring Inc. is a wholly owned subsidiary of Siemens Healthcare Diagnostics. In the short term, while we are working toward a legal integration of the two companies, please be aware that customer communications and legal documents will use either the Siemens Healthcare Diagnostics or the Dade Behring name, as appropriate.