

Urgent FIELD SAFETY NOTICE

Oxoid Ready Poured Plate *Helicobacter pylori* Selective Agar PB0398A Lot 1264570

16th August 2013

Field Safety Corrective Action (F.S.C.A.)

Contamination has been identified in Oxoid Ready Poured Plate *Helicobacter pylori* Selective Agar PB0398A Lot 1264570.

As the contamination may not be apparent prior to incubation, it is possible that growth of the contaminant could mask recovery of a microorganism from a sample. Use of this batch could lead to incorrect results reporting or a delay in results reporting due to the necessity for retesting. You should consider reviewing the outcomes of tests already performed using this batch of medium and seek appropriate expert advice where necessary.


In keeping with our Quality Policy, we request that you destroy all unused stock of this batch. Replacement packs, or a credit note, will be arranged via the Thermo Fisher Scientific, Microbiology UK, Customer Services Department.

This notice needs to be passed on to all those who need to be aware within your organisation.

The Medicines and Healthcare products Regulatory Agency (M.H.R.A.) has been informed of this F.S.C.A.

If you experience any problems complying with this notice, or you require further information, please contact the Technical Support Department on 01256 694281.

Please confirm receipt of this F.S.C.A., and advice of the removal and destruction of this batch, to our Customer Services Administrator.




Director of Regulatory Affairs
Oxoid & Remel
Microbiology UK

FSN/V3/10/2009



Oxoid Ltd.
Wade Road
Basingstoke

Hants
RG24 8PW
UK

+44 (0) 1256 841144 tel
+44 (0) 1256 814626 fax

www.thermofisher.com
www.oxoid.com