

Dear Zenzy customers,

This email is intended to inform you that a new version - Zenzy 2.46 Build 3 - is now available online on our website at www.heni-software.de for all our clients with a maintenance contract.

This new version now contains an additional validation check for the patient ID at the step at which you save the data after entry into the form on the screen.

During regular system maintenance, we observed that the old version allowed for generation of two patient records with the same patient ID if two or more users simultaneously entered and saved new data using the same patient ID, e.g. through entering manually the patient ID.

Although such an incident might be extremely rare, we recommend that you download the new software updates at your earliest convenience.

We would like to apologize in advance for any inconveniences this issue might have caused.

Sincerely,

Your Dr. Heni SOFTWARE GmbH & Co. KG Team

Dr. Heni SOFTWARE GmbH & Co. KG

Im Gässle 10

79199 Kirchzarten

Hotline +49 (0)7661 919 9999 · Fax +49 (0)7661 919 9990

E-Mail support@heni-software.de · Internetpräsenz www.heni-software.de