

30th April 2021

URGENT – FIELD SAFETY NOTICE

Type of Action	Recall
Teleflex Reference	EIF-000464
Commercial Name	RUSCHELIT® Super Safety Clear Tracheal Tube
Product Code/Lot Number	Refer to Appendix 2

Dear Customer,

Details of affected devices

Teleflex has initiated a voluntary Field Safety Corrective Action (FSCA) for RUSCHELIT® Super Safety Clear Tracheal Tube; refer to Appendix 2 for a list of impacted product codes, manufacturing and expiration dates.

Description of the problem & immediate actions required

Teleflex is initiating a Field Safety Corrective Action due to reports that the endotracheal tube cuff could be inflated while the pilot balloon remained flat. In the event of the pilot balloon remaining flat there is the potential risk of injury to the patient as the pilot balloon no longer reflects the condition of the cuff. The risk persists in situ, therefore clinicians should evaluate the risk benefit to individual patients of extubation and reintubation.

Our records indicate you have received products that are subject to this Field Safety Advisory Notification.

Depending on your device location please adhere to the following Action list:

Device location	Action List Number
Medical facilities (hospitals, medical staff, etc.)	1
Distributors	2

Action list number 1 – Medical facilities

1. We request that you check your inventory for product within the scope of this FSCA. Users should cease use and distribution of impacted product and quarantine immediately.
2. If you do not have stock in scope of this FSCA mark the according checkbox on the Acknowledgement Form (Appendix 1) and return the form to the fax number or e-Mail address mentioned below.
3. If you do have stock in scope of this FSCA, mark the according checkbox on the Acknowledgement Form (Appendix 1) and contact customer service by calling the phone number mentioned below. Customer service will issue you with a return number. Write the return number into the respective field in the Acknowledgement Form and return this form immediately to Customer Service.

4. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

Action list number 2 – Distributors

1. Provide this field safety notice to all customers who have received product in scope of this FSCA. Your customer is then required to complete the acknowledgement form and return to you.
2. We request that you check your inventory for product within the scope of this FSCA. Cease use and distribution of impacted product and quarantine immediately. You may then return all product in scope, refer to Appendix 2 for the list of impacted codes & lots, to Teleflex.
3. As a distributor, you are then required to confirm to Teleflex that you have completed the field activity outlined above. Upon completion of your actions, please forward the completed Acknowledgement Form to Customer Service.
4. Please be aware that all European Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities in which Teleflex distribute directly will be notified by Teleflex.
5. If you have further distributed product outside of your country, please notify Teleflex by return email to the e-Mail address below.
6. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/TR region, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

Teleflex

Teleflex informs all customers, employees of Teleflex and distributors of this Field Safety Corrective Action.

Transmission of this Field Safety Notice

This notice should be passed on to all persons who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. Please consider end users, clinicians, risk managers, supply chain/distribution centres etc. in the circulation of this notice. Maintain awareness of this notice until all required actions have been completed in your organisation.

Contact reference person

Should you require any further information or support concerning this issue, please contact:


Contact: Customer Service

FAX: 0711 / 49 05 08 71

Telephone: 0711 / 20 90 80 00

Email: recalls.de@teleflex.com

Please be advised that all European Economic Area/Switzerland (EEA/CH) and Turkey Member State Competent Authorities to which Teleflex distribute directly will be notified by Teleflex. Teleflex is committed to providing high quality, safe and effective products. We sincerely apologise for any inconvenience this action may cause your operations. If you have any other questions, feel free to contact your local sales representative or Customer Service.



FIELD SAFETY CORRECTIVE ACTION **ACKNOWLEDGEMENT FORM**

PRODUCT FIELD ACTION BY TELEFLEX – IMMEDIATE ATTENTION REQUIRED

Ref. EIF-000464

RETURN COMPLETED FORM IMMEDIATELY TO:

FAX: +353 (0) 1 4370773 **E-Mail:** Recalls.Intl@teleflex.com

<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm that our inventory does NOT include products affected by this Field Action.	<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm our inventory DOES include products affected by this Field Action. The use and further distribution of the affected products is stopped. All products are put on hold and the amount below will be returned. Return Authorisation No _____
---	--

PLEASE PRINT PRODUCT QUANTITY NUMBERS CLEARLY

PRODUCT NUMBER	LOT NUMBER	QUANTITY (Returning)
<ul style="list-style-type: none"> Include a copy of the completed Acknowledgement Form in the returns package with the returned units Ensure the RAN number is clearly visible on the returns package Please label returns as "Field Safety Returns" 		

Complete this Acknowledgement form and return immediately by using fax number or e-Mail address above.

INSTITUTION NAME (EG NAME OF HOSPITAL, HEALTH CARE ORGANISATION)	
INSITIUTION ADDRESS	Phone/FAX
FORM COMPLETED BY:	Stamp
PRINT NAME: _____	
SIGNATURE: _____	
DATE	

Appendix 2: Scope of product for EIF-000464

Product Code	Manufacturing Dates	Expiration Dates
112480-000050	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112480-000055	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112480-000060	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112480-000065	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112480-000070	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112480-000075	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112480-000080	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112480-000085	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112480-000090	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112480-000095	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112480-000100	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112482-000050	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112482-000055	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482-000060	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482-000065	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482-000070	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482-000075	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482-000080	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112482-000085	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482-000090	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482-000095	Apr 2020 through to Dec 2020	Mar 2025 through to Nov 2025
112482-000100	Apr 2020 through to Jan 2021	Mar 2025 through to Dec 2025
112482050LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482055LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482060LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482065LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482070LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482075LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482080LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482085LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482090LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025
112482095LA	Apr 2020 through to Nov 2020	Mar 2025 through to Oct 2025