

November 2022

**URGENT – FIELD SAFETY NOTICE**

Type of Action	Recall Notice
Teleflex Reference	EIF-000513-01
Product Codes	Commercial Name
19012	Iso-Gard® Filter S with Expandi-Flex/Elbow
19012T	Iso-Gard® Filter S with Expandi-Flex/Elbow, Tethered Cap
191667-000100	Filter + Catheter Mount
19211; 19212	Iso-Gard® Filter S
19261; 19272	Iso-Gard® Filter S with Expandi-Flex
19261T; 19262T; 19272T	Iso-Gard® Filter S with Expandi-Flex, Tethered Cap
19211T; 19212T	Iso-Gard® Filter S, Tethered Cap

Dear Customer,

**Details of impacted devices**

Teleflex Medical Europe Limited, has initiated a voluntary Field Safety Corrective Action (“FSCA”) for the above listed products; refer to Appendix 2 for Unique Device Identifier (“UDI”) and lot number information. This recall is an expansion to a recent recall (Teleflex reference EIF-000513) and includes additional product codes and lots based on the outcome of further investigation.

If you did not receive a copy of the Field Safety Notice for EIF-000513, Teleflex identified that you only received product in scope of the expansion, EIF-000513-01. Regardless, if you previously responded to the initial notification (EIF-000513), or are just receiving notification for the first time, Teleflex requests that you review the letter and complete the actions identified on page 2.

**Description of the issue & immediate actions required**

Teleflex is initiating this voluntary FSCA for the above-mentioned products due to reports of the device splitting or detaching during use. The devices are intended for connection to respiratory equipment in intensive care units and operating theatres to protect the equipment from potential airborne contaminants. If the filter splits or detaches during use, the functionality and efficiency of the device may be compromised, potentially resulting in a leakage and insufficient air supply to the patient with potential desaturation.

As of October 2022, a total of 38 complaints reporting the filter splitting or detaching have been received for the products in scope of this FSCA. Of these 38 complaints, 4 complaints reported injuries including patient desaturation. No deaths or long-term patient injuries have been reported at this time.

Our records indicate you have received products that are in scope of this FSCA.

Depending on your location please adhere to the following Action list:

Device location	Action List Number
Medical facilities (hospitals, medical staff, etc.)	1
Distributors	2

**Action list number 1 – Medical facilities**

1. We request that you immediately check your inventory for product within the scope of this FSCA. Users should cease use and distribution of affected product and immediately quarantine the affected product.
2. If you have impacted product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and contact Teleflex Customer Service by calling the phone number provided below. Teleflex Customer Service will issue a return goods authorisation (RGA) number to you. Write the (RGA) number into the respective field in the Acknowledgement Form and promptly return this form by email to Teleflex Customer Service.
3. If you do not have impacted product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and return the form to Teleflex at the contact details provided.
4. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

**Action list number 2 – Distributors**

1. Provide this field safety notice to all customers who have received impacted product. Each of your customers is then required to complete the Acknowledgement Form and return it to you.
2. We request that you immediately check your inventory for impacted product. Cease use and distribution of impacted product and immediately quarantine the affected product. You may then return all product in scope.
3. As a distributor, you are then required to confirm to Teleflex that you have completed the field activity outlined in actions 1 and 2 of this Action List Number 2. Upon completion of your actions, please forward the completed Acknowledgement Form to Teleflex Customer Service. **Important** - Please ensure you only list batch numbers in scope of this Field Safety Notice when completing this form.
4. Please be aware that all European Economic Area/Switzerland, United Kingdom (EEA/CH/UK) and Turkey (TR) Competent Authorities in which Teleflex distribute directly will be notified by Teleflex.
5. If you have further distributed product outside of your country, please notify Teleflex Customer Service by return email to the email address below.
6. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/UK/TR area, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

Adverse reactions or quality problems experienced with the use of this product should be reported to Teleflex Customer Service at the contact information below.

**Teleflex**

Teleflex informs impacted customers, employees of Teleflex and distributors of this FSCA.

**Transmission of this Field Safety Notice**

This notice should be passed on to all persons who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please consider end users, clinicians, risk managers, supply chain/distribution centres, etc., in the circulation of this notice. Please maintain awareness of this notice until all required actions have been completed in your organisation.

**Contact reference person**

Should you require further information or support concerning this issue, please contact Customer Service using the details provided below:

**Customer Service:****Contact:** Customer Service**FAX:** 0711 / 49 05 08 71**Telephone:** 0711 / 20 90 80 00**Email:** [recalls.de@teleflex.com](mailto:recalls.de@teleflex.com)

Teleflex is committed to providing high quality, safe and effective products. We regret any inconvenience this action may cause your operations. If you have any other questions, please contact your local Teleflex sales representative or Teleflex Customer Service.

*For and on behalf of Teleflex,*

*Padraig Hegarty*

*Padraig Hegarty VP, Global QA (Manufacturing)*

Appendix 1

Customer No  
\_\_\_\_\_

**FIELD SAFETY CORRECTIVE ACTION**  
**ACKNOWLEDGEMENT FORM**

**PRODUCT FIELD ACTION BY TELEFLEX – IMMEDIATE ATTENTION REQUIRED**

Ref. EIF-000513-01

RETURN COMPLETED FORM IMMEDIATELY TO:

E-Mail: [recalls.de@teleflex.com](mailto:recalls.de@teleflex.com)

<input type="checkbox"/> We confirm receipt of this FSN and completion of the required actions contained therein. We further confirm that our inventory does <b>NOT</b> include products impacted by this Field Action.	<input type="checkbox"/> We confirm receipt of this FSN and completion of the required actions contained therein. We further confirm our inventory <b>DOES</b> include products impacted by this Field Action. The use and further distribution of the affected products is stopped. All products are put on hold and the amount below will be returned. <b>Return Goods Authorisation No</b> _____
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Complete this Acknowledgement Form and return the completed form immediately using the contact information above.

Product code	Lot/batch number	Quantity returning
<b>Important</b> - Please ensure you only list batch numbers in scope of this Field Safety Notice when completing this form.		
<ul style="list-style-type: none"> <li>Include a copy of the <b>completed Acknowledgement Form</b> in the returns package with the returned units</li> <li>Ensure the <b>RGA number is clearly visible</b> on the returns package</li> <li>Please label returns as <b>"Field Safety Returns"</b></li> </ul>		
<b>Note:</b> Non-FSCA product returns should be processed per standard product return processes.		

<b>INSTITUTION NAME (E.G., NAME OF HOSPITAL, HEALTH CARE ORGANISATION)</b>	
<b>INSTITUTION ADDRESS</b>	<b>PHONE/FAX</b>
<b>FORM COMPLETED BY</b>	<b>STAMP</b>
PRINT NAME: _____	
SIGNATURE: _____	
DATE	

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19012	KMZ21A0099	(01)04026704388325(17)251228(10)KMZ21A0099
19211	19LT34	(01)04026704347995(17)250528(10)19LT34
19212	201951	(01)04026704348008(17)250528(10)201951
19212	KMZ20H0705	(01)04026704348008(17)250728(10)KMZ20H0705
19012T	202025	(01)04026704400423(17)250528(10)202025
19012T	202027	(01)04026704400423(17)250628(10)202027
19012T	202028	(01)04026704400423(17)250628(10)202028
19012T	KMZ20J0079	(01)04026704400423(17)250828(10)KMZ20J0079
19012T	KMZ20J0861	(01)04026704400423(17)250828(10)KMZ20J0861
19012T	KMZ21A0100	(01)04026704400423(17)251228(10)KMZ21A0100
19012T	KMZ21A0127	(01)04026704400423(17)251228(10)KMZ21A0127
19012T	KMZ21A0238	(01)04026704400423(17)251228(10)KMZ21A0238
19012T	KMZ21A0425	(01)04026704400423(17)251228(10)KMZ21A0425
19012T	KMZ21A0427	(01)04026704400423(17)251228(10)KMZ21A0427
19012T	KMZ21A0522	(01)04026704400423(17)251228(10)KMZ21A0522
19012T	KMZ21A0535	(01)04026704400423(17)251228(10)KMZ21A0535
19012T	KMZ21A0617	(01)04026704400423(17)251228(10)KMZ21A0617
19012T	KMZ21A0712	(01)04026704400423(17)251228(10)KMZ21A0712
19012T	KMZ21A0828	(01)04026704400423(17)251228(10)KMZ21A0828
19012T	KMZ21A0854	(01)04026704400423(17)251228(10)KMZ21A0854
19012T	KMZ21A0860	(01)04026704400423(17)251228(10)KMZ21A0860
19012T	KMZ21A0909	(01)04026704400423(17)251228(10)KMZ21A0909
19012T	KMZ21A0926	(01)04026704400423(17)251228(10)KMZ21A0926
19012T	KMZ21B0141	(01)04026704400423(17)260128(10)KMZ21B0141
19012T	KMZ21G0174	(01)04026704400423(17)260628(10)KMZ21G0174
19012T	KMZ21G0180	(01)04026704400423(17)260628(10)KMZ21G0180
19012T	KMZ21G0181	(01)04026704400423(17)260628(10)KMZ21G0181
19012T	KMZ21G1067	(01)04026704400423(17)260628(10)KMZ21G1067
19012T	KMZ21H0277	(01)04026704400423(17)260728(10)KMZ21H0277
19012T	KMZ21J0069	(01)04026704400423(17)260828(10)KMZ21J0069
19012T	KMZ21J0386	(01)04026704400423(17)260828(10)KMZ21J0386
19012T	KMZ21K0208	(01)04026704400423(17)260928(10)KMZ21K0208
19012T	KMZ21K0648	(01)04026704400423(17)260928(10)KMZ21K0648
19012T	KMZ21K0707	(01)04026704400423(17)260928(10)KMZ21K0707
19012T	KMZ21K0856	(01)04026704400423(17)260928(10)KMZ21K0856
19012T	KMZ21L0309	(01)04026704400423(17)261028(10)KMZ21L0309
19012T	KMZ21M0061	(01)04026704400423(17)261128(10)KMZ21M0061
19012T	KMZ21M0492	(01)04026704400423(17)261128(10)KMZ21M0492
19012T	KMZ21M0493	(01)04026704400423(17)261128(10)KMZ21M0493
19012T	KMZ21M0813	(01)04026704400423(17)261128(10)KMZ21M0813
19012T	KMZ22B0105	(01)04026704400423(17)270128(10)KMZ22B0105
19012T	KMZ22B0117	(01)04026704400423(17)270128(10)KMZ22B0117
19012T	KMZ22B0118	(01)04026704400423(17)270128(10)KMZ22B0118
19012T	KMZ22B0508	(01)04026704400423(17)270128(10)KMZ22B0508
19012T	KMZ22C0622	(01)04026704400423(17)270228(10)KMZ22C0622
191667-000100	KMH20H0092	(01)04026704561544(17)250728(10)KMH20H0092
191667-000100	KMH20H0109	(01)04026704561544(17)250728(10)KMH20H0109
191667-000100	KMH20L0312	(01)04026704561544(17)251028(10)KMH20L0312
191667-000100	KMH20L0344	(01)04026704561544(17)251028(10)KMH20L0344
191667-000100	KMH20M0237	(01)04026704561544(17)251128(10)KMH20M0237
191667-000100	KMH21A0033	(01)04026704561544(17)251228(10)KMH21A0033
191667-000100	KMH21B0197	(01)04026704561544(17)260128(10)KMH21B0197

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
191667-000100	KMH21D0181	(01)04026704561544(17)260328(10)KMH21D0181
191667-000100	KMH21J0103	(01)04026704561544(17)260828(10)KMH21J0103
191667-000100	KMH21J0150	(01)04026704561544(17)260828(10)KMH21J0150
191667-000100	KMH22A0181	(01)04026704561544(17)261228(10)KMH22A0181
191667-000100	KMH22A0185	(01)04026704561544(17)261228(10)KMH22A0185
19211	20FT02	(01)04026704347995(17)250528(10)20FT02
19211	20FT03	(01)04026704347995(17)250528(10)20FT03
19211	20FT06	(01)04026704347995(17)250528(10)20FT06
19211	20FT07	(01)04026704347995(17)250528(10)20FT07
19211	20FT13	(01)04026704347995(17)250528(10)20FT13
19211	20FT14	(01)04026704347995(17)250528(10)20FT14
19211	20FT16	(01)04026704347995(17)250528(10)20FT16
19211	20FT20	(01)04026704347995(17)250528(10)20FT20
19211	20FT22	(01)04026704347995(17)250528(10)20FT22
19211	20FT38	(01)04026704347995(17)250528(10)20FT38
19211	20FT41	(01)04026704347995(17)250528(10)20FT41
19211	20FT42	(01)04026704347995(17)250528(10)20FT42
19211	20FT50	(01)04026704347995(17)250528(10)20FT50
19211	20FT51	(01)04026704347995(17)250528(10)20FT51
19211	20FT56	(01)04026704347995(17)250528(10)20FT56
19211	20FT57	(01)04026704347995(17)250528(10)20FT57
19211	20FT62	(01)04026704347995(17)250528(10)20FT62
19211	20FT63	(01)04026704347995(17)250528(10)20FT63
19211	20FT67	(01)04026704347995(17)250528(10)20FT67
19211	20GT02	(01)04026704347995(17)250628(10)20GT02
19211	20GT03	(01)04026704347995(17)250628(10)20GT03
19211	20GT15	(01)04026704347995(17)250628(10)20GT15
19211	20GT17	(01)04026704347995(17)250628(10)20GT17
19211	20GT20	(01)04026704347995(17)250628(10)20GT20
19211	20GT21	(01)04026704347995(17)250628(10)20GT21
19211	20GT27	(01)04026704347995(17)250628(10)20GT27
19211	20GT28	(01)04026704347995(17)250628(10)20GT28
19211	20GT39	(01)04026704347995(17)250628(10)20GT39
19211	20GT46	(01)04026704347995(17)250628(10)20GT46
19211	20GT48	(01)04026704347995(17)250628(10)20GT48
19211	20GT53	(01)04026704347995(17)250628(10)20GT53
19211	20GT55	(01)04026704347995(17)250628(10)20GT55
19211	KMH20G0002	(01)04026704347995(17)250628(10)KMH20G0002
19211	KMH20G0003	(01)04026704347995(17)250628(10)KMH20G0003
19211	KMH20G0004	(01)04026704347995(17)250628(10)KMH20G0004
19211	KMH20G0010	(01)04026704347995(17)250628(10)KMH20G0010
19211	KMH20G0011	(01)04026704347995(17)250628(10)KMH20G0011
19211	KMH20G0014	(01)04026704347995(17)250628(10)KMH20G0014
19211	KMH20G0017	(01)04026704347995(17)250628(10)KMH20G0017
19211	KMH20G0019	(01)04026704347995(17)250628(10)KMH20G0019
19211	KMH20G0026	(01)04026704347995(17)250628(10)KMH20G0026
19211	KMH20G0032	(01)04026704347995(17)250628(10)KMH20G0032
19211	KMH20H0005	(01)04026704347995(17)250728(10)KMH20H0005
19211	KMH20H0011	(01)04026704347995(17)250728(10)KMH20H0011
19211	KMH20H0015	(01)04026704347995(17)250728(10)KMH20H0015
19211	KMH20H0017	(01)04026704347995(17)250728(10)KMH20H0017
19211	KMH20H0019	(01)04026704347995(17)250728(10)KMH20H0019

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH20H0028	(01)04026704347995(17)250728(10)KMH20H0028
19211	KMH20H0035	(01)04026704347995(17)250728(10)KMH20H0035
19211	KMH20H0038	(01)04026704347995(17)250728(10)KMH20H0038
19211	KMH20H0053	(01)04026704347995(17)250728(10)KMH20H0053
19211	KMH20H0086	(01)04026704347995(17)250728(10)KMH20H0086
19211	KMH20H0099	(01)04026704347995(17)250728(10)KMH20H0099
19211	KMH20H0105	(01)04026704347995(17)250728(10)KMH20H0105
19211	KMH20H0112	(01)04026704347995(17)250728(10)KMH20H0112
19211	KMH20H0132	(01)04026704347995(17)250728(10)KMH20H0132
19211	KMH20H0159	(01)04026704347995(17)250728(10)KMH20H0159
19211	KMH20H0165	(01)04026704347995(17)250728(10)KMH20H0165
19211	KMH20H0167	(01)04026704347995(17)250728(10)KMH20H0167
19211	KMH20H0174	(01)04026704347995(17)250728(10)KMH20H0174
19211	KMH20H0179	(01)04026704347995(17)250728(10)KMH20H0179
19211	KMH20H0180	(01)04026704347995(17)250728(10)KMH20H0180
19211	KMH20H0182	(01)04026704347995(17)250728(10)KMH20H0182
19211	KMH20H0184	(01)04026704347995(17)250728(10)KMH20H0184
19211	KMH20H0189	(01)04026704347995(17)250728(10)KMH20H0189
19211	KMH20H0194	(01)04026704347995(17)250728(10)KMH20H0194
19211	KMH20H0229	(01)04026704347995(17)250728(10)KMH20H0229
19211	KMH20H0231	(01)04026704347995(17)250728(10)KMH20H0231
19211	KMH20H0232	(01)04026704347995(17)250728(10)KMH20H0232
19211	KMH20H0233	(01)04026704347995(17)250728(10)KMH20H0233
19211	KMH20H0235	(01)04026704347995(17)250728(10)KMH20H0235
19211	KMH20H0252	(01)04026704347995(17)250728(10)KMH20H0252
19211	KMH20H0254	(01)04026704347995(17)250728(10)KMH20H0254
19211	KMH20H0257	(01)04026704347995(17)250728(10)KMH20H0257
19211	KMH20H0266	(01)04026704347995(17)250728(10)KMH20H0266
19211	KMH20H0268	(01)04026704347995(17)250728(10)KMH20H0268
19211	KMH20H0270	(01)04026704347995(17)250728(10)KMH20H0270
19211	KMH20H0290	(01)04026704347995(17)250728(10)KMH20H0290
19211	KMH20H0295	(01)04026704347995(17)250728(10)KMH20H0295
19211	KMH20H0298	(01)04026704347995(17)250728(10)KMH20H0298
19211	KMH20H0299	(01)04026704347995(17)250728(10)KMH20H0299
19211	KMH20H0301	(01)04026704347995(17)250728(10)KMH20H0301
19211	KMH20H0302	(01)04026704347995(17)250728(10)KMH20H0302
19211	KMH20H0307	(01)04026704347995(17)250728(10)KMH20H0307
19211	KMH20H0312	(01)04026704347995(17)250728(10)KMH20H0312
19211	KMH20H0313	(01)04026704347995(17)250728(10)KMH20H0313
19211	KMH20H0316	(01)04026704347995(17)250728(10)KMH20H0316
19211	KMH20H0318	(01)04026704347995(17)250728(10)KMH20H0318
19211	KMH20H0319	(01)04026704347995(17)250728(10)KMH20H0319
19211	KMH20J0001	(01)04026704347995(17)250828(10)KMH20J0001
19211	KMH20J0018	(01)04026704347995(17)250828(10)KMH20J0018
19211	KMH20J0023	(01)04026704347995(17)250828(10)KMH20J0023
19211	KMH20J0030	(01)04026704347995(17)250828(10)KMH20J0030
19211	KMH20J0037	(01)04026704347995(17)250828(10)KMH20J0037
19211	KMH20J0038	(01)04026704347995(17)250828(10)KMH20J0038
19211	KMH20J0042	(01)04026704347995(17)250828(10)KMH20J0042
19211	KMH20J0047	(01)04026704347995(17)250828(10)KMH20J0047
19211	KMH20J0052	(01)04026704347995(17)250828(10)KMH20J0052
19211	KMH20J0054	(01)04026704347995(17)250828(10)KMH20J0054

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH20J0067	(01)04026704347995(17)250828(10)KMH20J0067
19211	KMH20J0071	(01)04026704347995(17)250828(10)KMH20J0071
19211	KMH20J0072	(01)04026704347995(17)250828(10)KMH20J0072
19211	KMH20J0075	(01)04026704347995(17)250828(10)KMH20J0075
19211	KMH20J0076	(01)04026704347995(17)250828(10)KMH20J0076
19211	KMH20J0079	(01)04026704347995(17)250828(10)KMH20J0079
19211	KMH20J0082	(01)04026704347995(17)250828(10)KMH20J0082
19211	KMH20J0086	(01)04026704347995(17)250828(10)KMH20J0086
19211	KMH20J0091	(01)04026704347995(17)250828(10)KMH20J0091
19211	KMH20J0094	(01)04026704347995(17)250828(10)KMH20J0094
19211	KMH20J0097	(01)04026704347995(17)250828(10)KMH20J0097
19211	KMH20J0098	(01)04026704347995(17)250828(10)KMH20J0098
19211	KMH20J0100	(01)04026704347995(17)250828(10)KMH20J0100
19211	KMH20J0111	(01)04026704347995(17)250828(10)KMH20J0111
19211	KMH20J0113	(01)04026704347995(17)250828(10)KMH20J0113
19211	KMH20J0119	(01)04026704347995(17)250828(10)KMH20J0119
19211	KMH20J0130	(01)04026704347995(17)250828(10)KMH20J0130
19211	KMH20J0131	(01)04026704347995(17)250828(10)KMH20J0131
19211	KMH20J0132	(01)04026704347995(17)250828(10)KMH20J0132
19211	KMH20J0136	(01)04026704347995(17)250828(10)KMH20J0136
19211	KMH20J0150	(01)04026704347995(17)250828(10)KMH20J0150
19211	KMH20J0154	(01)04026704347995(17)250828(10)KMH20J0154
19211	KMH20J0159	(01)04026704347995(17)250828(10)KMH20J0159
19211	KMH20J0169	(01)04026704347995(17)250828(10)KMH20J0169
19211	KMH20J0171	(01)04026704347995(17)250828(10)KMH20J0171
19211	KMH20J0195	(01)04026704347995(17)250828(10)KMH20J0195
19211	KMH20J0198	(01)04026704347995(17)250828(10)KMH20J0198
19211	KMH20J0202	(01)04026704347995(17)250828(10)KMH20J0202
19211	KMH20J0204	(01)04026704347995(17)250828(10)KMH20J0204
19211	KMH20J0214	(01)04026704347995(17)250828(10)KMH20J0214
19211	KMH20J0232	(01)04026704347995(17)250828(10)KMH20J0232
19211	KMH20J0252	(01)04026704347995(17)250828(10)KMH20J0252
19211	KMH20J0256	(01)04026704347995(17)250828(10)KMH20J0256
19211	KMH20J0259	(01)04026704347995(17)250828(10)KMH20J0259
19211	KMH20J0262	(01)04026704347995(17)250828(10)KMH20J0262
19211	KMH20J0271	(01)04026704347995(17)250828(10)KMH20J0271
19211	KMH20J0274	(01)04026704347995(17)250828(10)KMH20J0274
19211	KMH20J0284	(01)04026704347995(17)250828(10)KMH20J0284
19211	KMH20J0286	(01)04026704347995(17)250828(10)KMH20J0286
19211	KMH20J0288	(01)04026704347995(17)250828(10)KMH20J0288
19211	KMH20J0289	(01)04026704347995(17)250828(10)KMH20J0289
19211	KMH20J0292	(01)04026704347995(17)250828(10)KMH20J0292
19211	KMH20J0297	(01)04026704347995(17)250828(10)KMH20J0297
19211	KMH20J0299	(01)04026704347995(17)250828(10)KMH20J0299
19211	KMH20J0305	(01)04026704347995(17)250828(10)KMH20J0305
19211	KMH20J0314	(01)04026704347995(17)250828(10)KMH20J0314
19211	KMH20J0316	(01)04026704347995(17)250828(10)KMH20J0316
19211	KMH20J0318	(01)04026704347995(17)250828(10)KMH20J0318
19211	KMH20J0320	(01)04026704347995(17)250828(10)KMH20J0320
19211	KMH20J0321	(01)04026704347995(17)250828(10)KMH20J0321
19211	KMH20J0325	(01)04026704347995(17)250828(10)KMH20J0325
19211	KMH20J0336	(01)04026704347995(17)250828(10)KMH20J0336



**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH20J0339	(01)04026704347995(17)250828(10)KMH20J0339
19211	KMH20J0347	(01)04026704347995(17)250828(10)KMH20J0347
19211	KMH20K0004	(01)04026704347995(17)250928(10)KMH20K0004
19211	KMH20K0013	(01)04026704347995(17)250928(10)KMH20K0013
19211	KMH20K0015	(01)04026704347995(17)250928(10)KMH20K0015
19211	KMH20K0023	(01)04026704347995(17)250928(10)KMH20K0023
19211	KMH20K0031	(01)04026704347995(17)250928(10)KMH20K0031
19211	KMH20K0035	(01)04026704347995(17)250928(10)KMH20K0035
19211	KMH20K0038	(01)04026704347995(17)250928(10)KMH20K0038
19211	KMH20K0040	(01)04026704347995(17)250928(10)KMH20K0040
19211	KMH20K0043	(01)04026704347995(17)250928(10)KMH20K0043
19211	KMH20K0045	(01)04026704347995(17)250928(10)KMH20K0045
19211	KMH20K0049	(01)04026704347995(17)250928(10)KMH20K0049
19211	KMH20K0052	(01)04026704347995(17)250928(10)KMH20K0052
19211	KMH20K0054	(01)04026704347995(17)250928(10)KMH20K0054
19211	KMH20K0057	(01)04026704347995(17)250928(10)KMH20K0057
19211	KMH20K0059	(01)04026704347995(17)250928(10)KMH20K0059
19211	KMH20K0060	(01)04026704347995(17)250928(10)KMH20K0060
19211	KMH20K0062	(01)04026704347995(17)250928(10)KMH20K0062
19211	KMH20K0063	(01)04026704347995(17)250928(10)KMH20K0063
19211	KMH20K0065	(01)04026704347995(17)250928(10)KMH20K0065
19211	KMH20K0066	(01)04026704347995(17)250928(10)KMH20K0066
19211	KMH20K0071	(01)04026704347995(17)250928(10)KMH20K0071
19211	KMH20K0076	(01)04026704347995(17)250928(10)KMH20K0076
19211	KMH20K0079	(01)04026704347995(17)250928(10)KMH20K0079
19211	KMH20K0082	(01)04026704347995(17)250928(10)KMH20K0082
19211	KMH20K0083	(01)04026704347995(17)250928(10)KMH20K0083
19211	KMH20K0085	(01)04026704347995(17)250928(10)KMH20K0085
19211	KMH20K0091	(01)04026704347995(17)250928(10)KMH20K0091
19211	KMH20K0095	(01)04026704347995(17)250928(10)KMH20K0095
19211	KMH20K0098	(01)04026704347995(17)250928(10)KMH20K0098
19211	KMH20K0104	(01)04026704347995(17)250928(10)KMH20K0104
19211	KMH20K0110	(01)04026704347995(17)250928(10)KMH20K0110
19211	KMH20K0112	(01)04026704347995(17)250928(10)KMH20K0112
19211	KMH20K0116	(01)04026704347995(17)250928(10)KMH20K0116
19211	KMH20K0145	(01)04026704347995(17)250928(10)KMH20K0145
19211	KMH20K0147	(01)04026704347995(17)250928(10)KMH20K0147
19211	KMH20K0155	(01)04026704347995(17)250928(10)KMH20K0155
19211	KMH20K0159	(01)04026704347995(17)250928(10)KMH20K0159
19211	KMH20K0160	(01)04026704347995(17)250928(10)KMH20K0160
19211	KMH20K0167	(01)04026704347995(17)250928(10)KMH20K0167
19211	KMH20K0191	(01)04026704347995(17)250928(10)KMH20K0191
19211	KMH20K0194	(01)04026704347995(17)250928(10)KMH20K0194
19211	KMH20K0207	(01)04026704347995(17)250928(10)KMH20K0207
19211	KMH20K0238	(01)04026704347995(17)250928(10)KMH20K0238
19211	KMH20K0253	(01)04026704347995(17)250928(10)KMH20K0253
19211	KMH20K0259	(01)04026704347995(17)250928(10)KMH20K0259
19211	KMH20K0273	(01)04026704347995(17)250928(10)KMH20K0273
19211	KMH20K0281	(01)04026704347995(17)250928(10)KMH20K0281
19211	KMH20K0283	(01)04026704347995(17)250928(10)KMH20K0283
19211	KMH20K0288	(01)04026704347995(17)250928(10)KMH20K0288
19211	KMH20K0294	(01)04026704347995(17)250928(10)KMH20K0294

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH20K0296	(01)04026704347995(17)250928(10)KMH20K0296
19211	KMH20K0300	(01)04026704347995(17)250928(10)KMH20K0300
19211	KMH20K0305	(01)04026704347995(17)250928(10)KMH20K0305
19211	KMH20K0323	(01)04026704347995(17)250928(10)KMH20K0323
19211	KMH20K0337	(01)04026704347995(17)250928(10)KMH20K0337
19211	KMH20K0342	(01)04026704347995(17)250928(10)KMH20K0342
19211	KMH20K0360	(01)04026704347995(17)250928(10)KMH20K0360
19211	KMH20K0409	(01)04026704347995(17)250928(10)KMH20K0409
19211	KMH20K0414	(01)04026704347995(17)250928(10)KMH20K0414
19211	KMH20K0415	(01)04026704347995(17)250928(10)KMH20K0415
19211	KMH20K0423	(01)04026704347995(17)250928(10)KMH20K0423
19211	KMH20K0430	(01)04026704347995(17)250928(10)KMH20K0430
19211	KMH20K0431	(01)04026704347995(17)250928(10)KMH20K0431
19211	KMH20K0434	(01)04026704347995(17)250928(10)KMH20K0434
19211	KMH20K0438	(01)04026704347995(17)250928(10)KMH20K0438
19211	KMH20K0461	(01)04026704347995(17)250928(10)KMH20K0461
19211	KMH20L0004	(01)04026704347995(17)251028(10)KMH20L0004
19211	KMH20L0033	(01)04026704347995(17)251028(10)KMH20L0033
19211	KMH20L0057	(01)04026704347995(17)251028(10)KMH20L0057
19211	KMH20L0081	(01)04026704347995(17)251028(10)KMH20L0081
19211	KMH20L0082	(01)04026704347995(17)251028(10)KMH20L0082
19211	KMH20L0089	(01)04026704347995(17)251028(10)KMH20L0089
19211	KMH20L0099	(01)04026704347995(17)251028(10)KMH20L0099
19211	KMH20L0106	(01)04026704347995(17)251028(10)KMH20L0106
19211	KMH20L0109	(01)04026704347995(17)251028(10)KMH20L0109
19211	KMH20L0148	(01)04026704347995(17)251028(10)KMH20L0148
19211	KMH20L0200	(01)04026704347995(17)251028(10)KMH20L0200
19211	KMH20L0204	(01)04026704347995(17)251028(10)KMH20L0204
19211	KMH20L0238	(01)04026704347995(17)251028(10)KMH20L0238
19211	KMH20L0239	(01)04026704347995(17)251028(10)KMH20L0239
19211	KMH20L0240	(01)04026704347995(17)251028(10)KMH20L0240
19211	KMH20L0241	(01)04026704347995(17)251028(10)KMH20L0241
19211	KMH20L0242	(01)04026704347995(17)251028(10)KMH20L0242
19211	KMH20L0243	(01)04026704347995(17)251028(10)KMH20L0243
19211	KMH20L0252	(01)04026704347995(17)251028(10)KMH20L0252
19211	KMH20L0278	(01)04026704347995(17)251028(10)KMH20L0278
19211	KMH20L0280	(01)04026704347995(17)251028(10)KMH20L0280
19211	KMH20L0290	(01)04026704347995(17)251028(10)KMH20L0290
19211	KMH20L0296	(01)04026704347995(17)251028(10)KMH20L0296
19211	KMH20L0302	(01)04026704347995(17)251028(10)KMH20L0302
19211	KMH20L0303	(01)04026704347995(17)251028(10)KMH20L0303
19211	KMH20L0305	(01)04026704347995(17)251028(10)KMH20L0305
19211	KMH20L0308	(01)04026704347995(17)251028(10)KMH20L0308
19211	KMH20L0310	(01)04026704347995(17)251028(10)KMH20L0310
19211	KMH20L0317	(01)04026704347995(17)251028(10)KMH20L0317
19211	KMH20L0322	(01)04026704347995(17)251028(10)KMH20L0322
19211	KMH20L0324	(01)04026704347995(17)251028(10)KMH20L0324
19211	KMH20L0329	(01)04026704347995(17)251028(10)KMH20L0329
19211	KMH20L0337	(01)04026704347995(17)251028(10)KMH20L0337
19211	KMH20L0339	(01)04026704347995(17)251028(10)KMH20L0339
19211	KMH20L0348	(01)04026704347995(17)251028(10)KMH20L0348
19211	KMH20L0376	(01)04026704347995(17)251028(10)KMH20L0376

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH20L0421	(01)04026704347995(17)251028(10)KMH20L0421
19211	KMH20M0093	(01)04026704347995(17)251128(10)KMH20M0093
19211	KMH20M0098	(01)04026704347995(17)251128(10)KMH20M0098
19211	KMH20M0103	(01)04026704347995(17)251128(10)KMH20M0103
19211	KMH20M0114	(01)04026704347995(17)251128(10)KMH20M0114
19211	KMH20M0121	(01)04026704347995(17)251128(10)KMH20M0121
19211	KMH20M0124	(01)04026704347995(17)251128(10)KMH20M0124
19211	KMH20M0131	(01)04026704347995(17)251128(10)KMH20M0131
19211	KMH20M0146	(01)04026704347995(17)251128(10)KMH20M0146
19211	KMH20M0196	(01)04026704347995(17)251128(10)KMH20M0196
19211	KMH20M0203	(01)04026704347995(17)251128(10)KMH20M0203
19211	KMH20M0211	(01)04026704347995(17)251128(10)KMH20M0211
19211	KMH20M0221	(01)04026704347995(17)251128(10)KMH20M0221
19211	KMH20M0227	(01)04026704347995(17)251128(10)KMH20M0227
19211	KMH20M0231	(01)04026704347995(17)251128(10)KMH20M0231
19211	KMH20M0235	(01)04026704347995(17)251128(10)KMH20M0235
19211	KMH20M0239	(01)04026704347995(17)251128(10)KMH20M0239
19211	KMH20M0241	(01)04026704347995(17)251128(10)KMH20M0241
19211	KMH20M0243	(01)04026704347995(17)251128(10)KMH20M0243
19211	KMH20M0244	(01)04026704347995(17)251128(10)KMH20M0244
19211	KMH20M0247	(01)04026704347995(17)251128(10)KMH20M0247
19211	KMH20M0251	(01)04026704347995(17)251128(10)KMH20M0251
19211	KMH20M0252	(01)04026704347995(17)251128(10)KMH20M0252
19211	KMH20M0260	(01)04026704347995(17)251128(10)KMH20M0260
19211	KMH20M0263	(01)04026704347995(17)251128(10)KMH20M0263
19211	KMH20M0266	(01)04026704347995(17)251128(10)KMH20M0266
19211	KMH20M0267	(01)04026704347995(17)251128(10)KMH20M0267
19211	KMH20M0269	(01)04026704347995(17)251128(10)KMH20M0269
19211	KMH20M0278	(01)04026704347995(17)251128(10)KMH20M0278
19211	KMH20M0292	(01)04026704347995(17)251128(10)KMH20M0292
19211	KMH20M0298	(01)04026704347995(17)251128(10)KMH20M0298
19211	KMH20M0301	(01)04026704347995(17)251128(10)KMH20M0301
19211	KMH20M0316	(01)04026704347995(17)251128(10)KMH20M0316
19211	KMH20M0329	(01)04026704347995(17)251128(10)KMH20M0329
19211	KMH20M0426	(01)04026704347995(17)251128(10)KMH20M0426
19211	KMH20M0427	(01)04026704347995(17)251128(10)KMH20M0427
19211	KMH20M0461	(01)04026704347995(17)251128(10)KMH20M0461
19211	KMH20M0468	(01)04026704347995(17)251128(10)KMH20M0468
19211	KMH20M0481	(01)04026704347995(17)251128(10)KMH20M0481
19211	KMH20M0526	(01)04026704347995(17)251128(10)KMH20M0526
19211	KMH21A0006	(01)04026704347995(17)251228(10)KMH21A0006
19211	KMH21A0008	(01)04026704347995(17)251228(10)KMH21A0008
19211	KMH21A0021	(01)04026704347995(17)251228(10)KMH21A0021
19211	KMH21A0024	(01)04026704347995(17)251228(10)KMH21A0024
19211	KMH21A0027	(01)04026704347995(17)251228(10)KMH21A0027
19211	KMH21A0028	(01)04026704347995(17)251228(10)KMH21A0028
19211	KMH21A0034	(01)04026704347995(17)251228(10)KMH21A0034
19211	KMH21A0064	(01)04026704347995(17)251228(10)KMH21A0064
19211	KMH21A0087	(01)04026704347995(17)251228(10)KMH21A0087
19211	KMH21A0096	(01)04026704347995(17)251228(10)KMH21A0096
19211	KMH21A0102	(01)04026704347995(17)251228(10)KMH21A0102
19211	KMH21A0120	(01)04026704347995(17)251228(10)KMH21A0120

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21A0124	(01)04026704347995(17)251228(10)KMH21A0124
19211	KMH21A0128	(01)04026704347995(17)251228(10)KMH21A0128
19211	KMH21A0135	(01)04026704347995(17)251228(10)KMH21A0135
19211	KMH21A0147	(01)04026704347995(17)251228(10)KMH21A0147
19211	KMH21A0148	(01)04026704347995(17)251228(10)KMH21A0148
19211	KMH21A0151	(01)04026704347995(17)251228(10)KMH21A0151
19211	KMH21A0154	(01)04026704347995(17)251228(10)KMH21A0154
19211	KMH21A0159	(01)04026704347995(17)251228(10)KMH21A0159
19211	KMH21A0163	(01)04026704347995(17)251228(10)KMH21A0163
19211	KMH21A0168	(01)04026704347995(17)251228(10)KMH21A0168
19211	KMH21A0170	(01)04026704347995(17)251228(10)KMH21A0170
19211	KMH21A0189	(01)04026704347995(17)251228(10)KMH21A0189
19211	KMH21A0198	(01)04026704347995(17)251228(10)KMH21A0198
19211	KMH21A0202	(01)04026704347995(17)251228(10)KMH21A0202
19211	KMH21A0204	(01)04026704347995(17)251228(10)KMH21A0204
19211	KMH21A0216	(01)04026704347995(17)251228(10)KMH21A0216
19211	KMH21A0225	(01)04026704347995(17)251228(10)KMH21A0225
19211	KMH21A0321	(01)04026704347995(17)251228(10)KMH21A0321
19211	KMH21A0330	(01)04026704347995(17)251228(10)KMH21A0330
19211	KMH21A0333	(01)04026704347995(17)251228(10)KMH21A0333
19211	KMH21A0337	(01)04026704347995(17)251228(10)KMH21A0337
19211	KMH21A0344	(01)04026704347995(17)251228(10)KMH21A0344
19211	KMH21A0345	(01)04026704347995(17)251228(10)KMH21A0345
19211	KMH21A0349	(01)04026704347995(17)251228(10)KMH21A0349
19211	KMH21A0365	(01)04026704347995(17)251228(10)KMH21A0365
19211	KMH21A0379	(01)04026704347995(17)251228(10)KMH21A0379
19211	KMH21A0384	(01)04026704347995(17)251228(10)KMH21A0384
19211	KMH21A0400	(01)04026704347995(17)251228(10)KMH21A0400
19211	KMH21A0407	(01)04026704347995(17)251228(10)KMH21A0407
19211	KMH21A0410	(01)04026704347995(17)251228(10)KMH21A0410
19211	KMH21A0412	(01)04026704347995(17)251228(10)KMH21A0412
19211	KMH21A0418	(01)04026704347995(17)251228(10)KMH21A0418
19211	KMH21A0427	(01)04026704347995(17)251228(10)KMH21A0427
19211	KMH21A0428	(01)04026704347995(17)251228(10)KMH21A0428
19211	KMH21A0430	(01)04026704347995(17)251228(10)KMH21A0430
19211	KMH21A0434	(01)04026704347995(17)251228(10)KMH21A0434
19211	KMH21A0458	(01)04026704347995(17)251228(10)KMH21A0458
19211	KMH21A0459	(01)04026704347995(17)251228(10)KMH21A0459
19211	KMH21B0005	(01)04026704347995(17)260128(10)KMH21B0005
19211	KMH21B0007	(01)04026704347995(17)260128(10)KMH21B0007
19211	KMH21B0009	(01)04026704347995(17)260128(10)KMH21B0009
19211	KMH21B0048	(01)04026704347995(17)260128(10)KMH21B0048
19211	KMH21B0053	(01)04026704347995(17)260128(10)KMH21B0053
19211	KMH21B0069	(01)04026704347995(17)260128(10)KMH21B0069
19211	KMH21B0078	(01)04026704347995(17)260128(10)KMH21B0078
19211	KMH21B0084	(01)04026704347995(17)260128(10)KMH21B0084
19211	KMH21B0100	(01)04026704347995(17)260128(10)KMH21B0100
19211	KMH21B0101	(01)04026704347995(17)260128(10)KMH21B0101
19211	KMH21B0107	(01)04026704347995(17)260128(10)KMH21B0107
19211	KMH21B0109	(01)04026704347995(17)260128(10)KMH21B0109
19211	KMH21B0111	(01)04026704347995(17)260128(10)KMH21B0111
19211	KMH21B0118	(01)04026704347995(17)260128(10)KMH21B0118

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21B0126	(01)04026704347995(17)260128(10)KMH21B0126
19211	KMH21B0139	(01)04026704347995(17)260128(10)KMH21B0139
19211	KMH21B0143	(01)04026704347995(17)260128(10)KMH21B0143
19211	KMH21B0162	(01)04026704347995(17)260128(10)KMH21B0162
19211	KMH21C0090	(01)04026704347995(17)260228(10)KMH21C0090
19211	KMH21C0091	(01)04026704347995(17)260228(10)KMH21C0091
19211	KMH21C0100	(01)04026704347995(17)260228(10)KMH21C0100
19211	KMH21C0101	(01)04026704347995(17)260228(10)KMH21C0101
19211	KMH21C0105	(01)04026704347995(17)260228(10)KMH21C0105
19211	KMH21C0108	(01)04026704347995(17)260228(10)KMH21C0108
19211	KMH21C0109	(01)04026704347995(17)260228(10)KMH21C0109
19211	KMH21C0112	(01)04026704347995(17)260228(10)KMH21C0112
19211	KMH21C0114	(01)04026704347995(17)260228(10)KMH21C0114
19211	KMH21C0118	(01)04026704347995(17)260228(10)KMH21C0118
19211	KMH21C0122	(01)04026704347995(17)260228(10)KMH21C0122
19211	KMH21C0124	(01)04026704347995(17)260228(10)KMH21C0124
19211	KMH21C0125	(01)04026704347995(17)260228(10)KMH21C0125
19211	KMH21C0282	(01)04026704347995(17)260228(10)KMH21C0282
19211	KMH21C0286	(01)04026704347995(17)260228(10)KMH21C0286
19211	KMH21C0292	(01)04026704347995(17)260228(10)KMH21C0292
19211	KMH21C0296	(01)04026704347995(17)260228(10)KMH21C0296
19211	KMH21C0303	(01)04026704347995(17)260228(10)KMH21C0303
19211	KMH21C0306	(01)04026704347995(17)260228(10)KMH21C0306
19211	KMH21C0315	(01)04026704347995(17)260228(10)KMH21C0315
19211	KMH21C0318	(01)04026704347995(17)260228(10)KMH21C0318
19211	KMH21C0320	(01)04026704347995(17)260228(10)KMH21C0320
19211	KMH21C0337	(01)04026704347995(17)260228(10)KMH21C0337
19211	KMH21C0341	(01)04026704347995(17)260228(10)KMH21C0341
19211	KMH21C0348	(01)04026704347995(17)260228(10)KMH21C0348
19211	KMH21C0349	(01)04026704347995(17)260228(10)KMH21C0349
19211	KMH21C0357	(01)04026704347995(17)260228(10)KMH21C0357
19211	KMH21C0372	(01)04026704347995(17)260228(10)KMH21C0372
19211	KMH21C0373	(01)04026704347995(17)260228(10)KMH21C0373
19211	KMH21C0375	(01)04026704347995(17)260228(10)KMH21C0375
19211	KMH21C0380	(01)04026704347995(17)260228(10)KMH21C0380
19211	KMH21C0383	(01)04026704347995(17)260228(10)KMH21C0383
19211	KMH21C0384	(01)04026704347995(17)260228(10)KMH21C0384
19211	KMH21C0395	(01)04026704347995(17)260228(10)KMH21C0395
19211	KMH21C0396	(01)04026704347995(17)260228(10)KMH21C0396
19211	KMH21C0405	(01)04026704347995(17)260228(10)KMH21C0405
19211	KMH21C0411	(01)04026704347995(17)260228(10)KMH21C0411
19211	KMH21C0420	(01)04026704347995(17)260228(10)KMH21C0420
19211	KMH21C0423	(01)04026704347995(17)260228(10)KMH21C0423
19211	KMH21C0427	(01)04026704347995(17)260228(10)KMH21C0427
19211	KMH21C0437	(01)04026704347995(17)260228(10)KMH21C0437
19211	KMH21C0441	(01)04026704347995(17)260228(10)KMH21C0441
19211	KMH21C0445	(01)04026704347995(17)260228(10)KMH21C0445
19211	KMH21C0453	(01)04026704347995(17)260228(10)KMH21C0453
19211	KMH21D0045	(01)04026704347995(17)260328(10)KMH21D0045
19211	KMH21D0057	(01)04026704347995(17)260328(10)KMH21D0057
19211	KMH21D0060	(01)04026704347995(17)260328(10)KMH21D0060
19211	KMH21D0070	(01)04026704347995(17)260328(10)KMH21D0070

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21D0074	(01)04026704347995(17)260328(10)KMH21D0074
19211	KMH21D0076	(01)04026704347995(17)260328(10)KMH21D0076
19211	KMH21D0081	(01)04026704347995(17)260328(10)KMH21D0081
19211	KMH21D0082	(01)04026704347995(17)260328(10)KMH21D0082
19211	KMH21D0088	(01)04026704347995(17)260328(10)KMH21D0088
19211	KMH21D0097	(01)04026704347995(17)260328(10)KMH21D0097
19211	KMH21D0103	(01)04026704347995(17)260328(10)KMH21D0103
19211	KMH21D0108	(01)04026704347995(17)260328(10)KMH21D0108
19211	KMH21D0110	(01)04026704347995(17)260328(10)KMH21D0110
19211	KMH21D0114	(01)04026704347995(17)260328(10)KMH21D0114
19211	KMH21D0120	(01)04026704347995(17)260328(10)KMH21D0120
19211	KMH21D0126	(01)04026704347995(17)260328(10)KMH21D0126
19211	KMH21D0151	(01)04026704347995(17)260328(10)KMH21D0151
19211	KMH21D0157	(01)04026704347995(17)260328(10)KMH21D0157
19211	KMH21D0165	(01)04026704347995(17)260328(10)KMH21D0165
19211	KMH21D0170	(01)04026704347995(17)260328(10)KMH21D0170
19211	KMH21D0173	(01)04026704347995(17)260328(10)KMH21D0173
19211	KMH21D0179	(01)04026704347995(17)260328(10)KMH21D0179
19211	KMH21D0193	(01)04026704347995(17)260328(10)KMH21D0193
19211	KMH21D0198	(01)04026704347995(17)260328(10)KMH21D0198
19211	KMH21D0203	(01)04026704347995(17)260328(10)KMH21D0203
19211	KMH21D0207	(01)04026704347995(17)260328(10)KMH21D0207
19211	KMH21D0212	(01)04026704347995(17)260328(10)KMH21D0212
19211	KMH21D0219	(01)04026704347995(17)260328(10)KMH21D0219
19211	KMH21D0224	(01)04026704347995(17)260328(10)KMH21D0224
19211	KMH21D0226	(01)04026704347995(17)260328(10)KMH21D0226
19211	KMH21D0231	(01)04026704347995(17)260328(10)KMH21D0231
19211	KMH21D0233	(01)04026704347995(17)260328(10)KMH21D0233
19211	KMH21D0236	(01)04026704347995(17)260328(10)KMH21D0236
19211	KMH21D0240	(01)04026704347995(17)260328(10)KMH21D0240
19211	KMH21D0244	(01)04026704347995(17)260328(10)KMH21D0244
19211	KMH21D0247	(01)04026704347995(17)260328(10)KMH21D0247
19211	KMH21D0299	(01)04026704347995(17)260328(10)KMH21D0299
19211	KMH21D0300	(01)04026704347995(17)260328(10)KMH21D0300
19211	KMH21D0305	(01)04026704347995(17)260328(10)KMH21D0305
19211	KMH21D0307	(01)04026704347995(17)260328(10)KMH21D0307
19211	KMH21D0309	(01)04026704347995(17)260328(10)KMH21D0309
19211	KMH21D0323	(01)04026704347995(17)260328(10)KMH21D0323
19211	KMH21D0325	(01)04026704347995(17)260328(10)KMH21D0325
19211	KMH21D0326	(01)04026704347995(17)260328(10)KMH21D0326
19211	KMH21D0333	(01)04026704347995(17)260328(10)KMH21D0333
19211	KMH21D0354	(01)04026704347995(17)260328(10)KMH21D0354
19211	KMH21D0357	(01)04026704347995(17)260328(10)KMH21D0357
19211	KMH21D0361	(01)04026704347995(17)260328(10)KMH21D0361
19211	KMH21D0363	(01)04026704347995(17)260328(10)KMH21D0363
19211	KMH21E0001	(01)04026704347995(17)260428(10)KMH21E0001
19211	KMH21E0006	(01)04026704347995(17)260428(10)KMH21E0006
19211	KMH21E0118	(01)04026704347995(17)260428(10)KMH21E0118
19211	KMH21E0162	(01)04026704347995(17)260428(10)KMH21E0162
19211	KMH21E0193	(01)04026704347995(17)260428(10)KMH21E0193
19211	KMH21E0268	(01)04026704347995(17)260428(10)KMH21E0268
19211	KMH21F0007	(01)04026704347995(17)260528(10)KMH21F0007

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21F0017	(01)04026704347995(17)260528(10)KMH21F0017
19211	KMH21F0023	(01)04026704347995(17)260528(10)KMH21F0023
19211	KMH21F0176	(01)04026704347995(17)260528(10)KMH21F0176
19211	KMH21F0178	(01)04026704347995(17)260528(10)KMH21F0178
19211	KMH21F0182	(01)04026704347995(17)260528(10)KMH21F0182
19211	KMH21F0204	(01)04026704347995(17)260528(10)KMH21F0204
19211	KMH21F0212	(01)04026704347995(17)260528(10)KMH21F0212
19211	KMH21F0213	(01)04026704347995(17)260528(10)KMH21F0213
19211	KMH21F0220	(01)04026704347995(17)260528(10)KMH21F0220
19211	KMH21F0266	(01)04026704347995(17)260528(10)KMH21F0266
19211	KMH21F0288	(01)04026704347995(17)260528(10)KMH21F0288
19211	KMH21F0289	(01)04026704347995(17)260528(10)KMH21F0289
19211	KMH21F0290	(01)04026704347995(17)260528(10)KMH21F0290
19211	KMH21F0298	(01)04026704347995(17)260528(10)KMH21F0298
19211	KMH21F0300	(01)04026704347995(17)260528(10)KMH21F0300
19211	KMH21F0319	(01)04026704347995(17)260528(10)KMH21F0319
19211	KMH21F0352	(01)04026704347995(17)260528(10)KMH21F0352
19211	KMH21F0369	(01)04026704347995(17)260528(10)KMH21F0369
19211	KMH21F0371	(01)04026704347995(17)260528(10)KMH21F0371
19211	KMH21F0376	(01)04026704347995(17)260528(10)KMH21F0376
19211	KMH21F0378	(01)04026704347995(17)260528(10)KMH21F0378
19211	KMH21F0381	(01)04026704347995(17)260528(10)KMH21F0381
19211	KMH21F0387	(01)04026704347995(17)260528(10)KMH21F0387
19211	KMH21F0389	(01)04026704347995(17)260528(10)KMH21F0389
19211	KMH21F0396	(01)04026704347995(17)260528(10)KMH21F0396
19211	KMH21F0401	(01)04026704347995(17)260528(10)KMH21F0401
19211	KMH21F0404	(01)04026704347995(17)260528(10)KMH21F0404
19211	KMH21F0407	(01)04026704347995(17)260528(10)KMH21F0407
19211	KMH21F0410	(01)04026704347995(17)260528(10)KMH21F0410
19211	KMH21G0004	(01)04026704347995(17)260628(10)KMH21G0004
19211	KMH21G0009	(01)04026704347995(17)260628(10)KMH21G0009
19211	KMH21G0012	(01)04026704347995(17)260628(10)KMH21G0012
19211	KMH21G0016	(01)04026704347995(17)260628(10)KMH21G0016
19211	KMH21G0021	(01)04026704347995(17)260628(10)KMH21G0021
19211	KMH21G0025	(01)04026704347995(17)260628(10)KMH21G0025
19211	KMH21G0026	(01)04026704347995(17)260628(10)KMH21G0026
19211	KMH21G0030	(01)04026704347995(17)260628(10)KMH21G0030
19211	KMH21G0035	(01)04026704347995(17)260628(10)KMH21G0035
19211	KMH21G0036	(01)04026704347995(17)260628(10)KMH21G0036
19211	KMH21G0042	(01)04026704347995(17)260628(10)KMH21G0042
19211	KMH21G0046	(01)04026704347995(17)260628(10)KMH21G0046
19211	KMH21G0050	(01)04026704347995(17)260628(10)KMH21G0050
19211	KMH21G0051	(01)04026704347995(17)260628(10)KMH21G0051
19211	KMH21G0052	(01)04026704347995(17)260628(10)KMH21G0052
19211	KMH21G0054	(01)04026704347995(17)260628(10)KMH21G0054
19211	KMH21G0055	(01)04026704347995(17)260628(10)KMH21G0055
19211	KMH21G0056	(01)04026704347995(17)260628(10)KMH21G0056
19211	KMH21G0060	(01)04026704347995(17)260628(10)KMH21G0060
19211	KMH21G0064	(01)04026704347995(17)260628(10)KMH21G0064
19211	KMH21G0070	(01)04026704347995(17)260628(10)KMH21G0070
19211	KMH21G0074	(01)04026704347995(17)260628(10)KMH21G0074
19211	KMH21G0082	(01)04026704347995(17)260628(10)KMH21G0082

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21G0086	(01)04026704347995(17)260628(10)KMH21G0086
19211	KMH21G0088	(01)04026704347995(17)260628(10)KMH21G0088
19211	KMH21G0091	(01)04026704347995(17)260628(10)KMH21G0091
19211	KMH21G0093	(01)04026704347995(17)260628(10)KMH21G0093
19211	KMH21G0098	(01)04026704347995(17)260628(10)KMH21G0098
19211	KMH21G0104	(01)04026704347995(17)260628(10)KMH21G0104
19211	KMH21G0114	(01)04026704347995(17)260628(10)KMH21G0114
19211	KMH21G0119	(01)04026704347995(17)260628(10)KMH21G0119
19211	KMH21G0126	(01)04026704347995(17)260628(10)KMH21G0126
19211	KMH21G0128	(01)04026704347995(17)260628(10)KMH21G0128
19211	KMH21G0135	(01)04026704347995(17)260628(10)KMH21G0135
19211	KMH21G0139	(01)04026704347995(17)260628(10)KMH21G0139
19211	KMH21G0142	(01)04026704347995(17)260628(10)KMH21G0142
19211	KMH21G0146	(01)04026704347995(17)260628(10)KMH21G0146
19211	KMH21G0151	(01)04026704347995(17)260628(10)KMH21G0151
19211	KMH21G0157	(01)04026704347995(17)260628(10)KMH21G0157
19211	KMH21G0158	(01)04026704347995(17)260628(10)KMH21G0158
19211	KMH21G0159	(01)04026704347995(17)260628(10)KMH21G0159
19211	KMH21G0172	(01)04026704347995(17)260628(10)KMH21G0172
19211	KMH21G0184	(01)04026704347995(17)260628(10)KMH21G0184
19211	KMH21G0186	(01)04026704347995(17)260628(10)KMH21G0186
19211	KMH21G0247	(01)04026704347995(17)260628(10)KMH21G0247
19211	KMH21G0262	(01)04026704347995(17)260628(10)KMH21G0262
19211	KMH21G0267	(01)04026704347995(17)260628(10)KMH21G0267
19211	KMH21G0268	(01)04026704347995(17)260628(10)KMH21G0268
19211	KMH21G0271	(01)04026704347995(17)260628(10)KMH21G0271
19211	KMH21G0279	(01)04026704347995(17)260628(10)KMH21G0279
19211	KMH21G0282	(01)04026704347995(17)260628(10)KMH21G0282
19211	KMH21G0285	(01)04026704347995(17)260628(10)KMH21G0285
19211	KMH21G0287	(01)04026704347995(17)260628(10)KMH21G0287
19211	KMH21G0290	(01)04026704347995(17)260628(10)KMH21G0290
19211	KMH21G0291	(01)04026704347995(17)260628(10)KMH21G0291
19211	KMH21G0296	(01)04026704347995(17)260628(10)KMH21G0296
19211	KMH21G0300	(01)04026704347995(17)260628(10)KMH21G0300
19211	KMH21G0307	(01)04026704347995(17)260628(10)KMH21G0307
19211	KMH21G0310	(01)04026704347995(17)260628(10)KMH21G0310
19211	KMH21G0314	(01)04026704347995(17)260628(10)KMH21G0314
19211	KMH21G0319	(01)04026704347995(17)260628(10)KMH21G0319
19211	KMH21G0320	(01)04026704347995(17)260628(10)KMH21G0320
19211	KMH21G0323	(01)04026704347995(17)260628(10)KMH21G0323
19211	KMH21G0332	(01)04026704347995(17)260628(10)KMH21G0332
19211	KMH21G0335	(01)04026704347995(17)260628(10)KMH21G0335
19211	KMH21G0338	(01)04026704347995(17)260628(10)KMH21G0338
19211	KMH21G0345	(01)04026704347995(17)260628(10)KMH21G0345
19211	KMH21G0353	(01)04026704347995(17)260628(10)KMH21G0353
19211	KMH21G0356	(01)04026704347995(17)260628(10)KMH21G0356
19211	KMH21G0358	(01)04026704347995(17)260628(10)KMH21G0358
19211	KMH21G0359	(01)04026704347995(17)260628(10)KMH21G0359
19211	KMH21G0361	(01)04026704347995(17)260628(10)KMH21G0361
19211	KMH21G0365	(01)04026704347995(17)260628(10)KMH21G0365
19211	KMH21G0368	(01)04026704347995(17)260628(10)KMH21G0368
19211	KMH21G0369	(01)04026704347995(17)260628(10)KMH21G0369



**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21G0371	(01)04026704347995(17)260628(10)KMH21G0371
19211	KMH21G0373	(01)04026704347995(17)260628(10)KMH21G0373
19211	KMH21G0376	(01)04026704347995(17)260628(10)KMH21G0376
19211	KMH21G0378	(01)04026704347995(17)260628(10)KMH21G0378
19211	KMH21G0381	(01)04026704347995(17)260628(10)KMH21G0381
19211	KMH21G0388	(01)04026704347995(17)260628(10)KMH21G0388
19211	KMH21G0390	(01)04026704347995(17)260628(10)KMH21G0390
19211	KMH21G0391	(01)04026704347995(17)260628(10)KMH21G0391
19211	KMH21G0425	(01)04026704347995(17)260628(10)KMH21G0425
19211	KMH21G0427	(01)04026704347995(17)260628(10)KMH21G0427
19211	KMH21G0433	(01)04026704347995(17)260628(10)KMH21G0433
19211	KMH21G0438	(01)04026704347995(17)260628(10)KMH21G0438
19211	KMH21H0002	(01)04026704347995(17)260728(10)KMH21H0002
19211	KMH21H0008	(01)04026704347995(17)260728(10)KMH21H0008
19211	KMH21H0014	(01)04026704347995(17)260728(10)KMH21H0014
19211	KMH21H0016	(01)04026704347995(17)260728(10)KMH21H0016
19211	KMH21H0019	(01)04026704347995(17)260728(10)KMH21H0019
19211	KMH21H0020	(01)04026704347995(17)260728(10)KMH21H0020
19211	KMH21H0021	(01)04026704347995(17)260728(10)KMH21H0021
19211	KMH21H0027	(01)04026704347995(17)260728(10)KMH21H0027
19211	KMH21H0029	(01)04026704347995(17)260728(10)KMH21H0029
19211	KMH21H0031	(01)04026704347995(17)260728(10)KMH21H0031
19211	KMH21H0033	(01)04026704347995(17)260728(10)KMH21H0033
19211	KMH21H0037	(01)04026704347995(17)260728(10)KMH21H0037
19211	KMH21H0038	(01)04026704347995(17)260728(10)KMH21H0038
19211	KMH21H0044	(01)04026704347995(17)260728(10)KMH21H0044
19211	KMH21H0045	(01)04026704347995(17)260728(10)KMH21H0045
19211	KMH21H0046	(01)04026704347995(17)260728(10)KMH21H0046
19211	KMH21H0052	(01)04026704347995(17)260728(10)KMH21H0052
19211	KMH21H0053	(01)04026704347995(17)260728(10)KMH21H0053
19211	KMH21H0055	(01)04026704347995(17)260728(10)KMH21H0055
19211	KMH21H0063	(01)04026704347995(17)260728(10)KMH21H0063
19211	KMH21H0071	(01)04026704347995(17)260728(10)KMH21H0071
19211	KMH21H0075	(01)04026704347995(17)260728(10)KMH21H0075
19211	KMH21H0076	(01)04026704347995(17)260728(10)KMH21H0076
19211	KMH21H0077	(01)04026704347995(17)260728(10)KMH21H0077
19211	KMH21H0082	(01)04026704347995(17)260728(10)KMH21H0082
19211	KMH21H0084	(01)04026704347995(17)260728(10)KMH21H0084
19211	KMH21H0090	(01)04026704347995(17)260728(10)KMH21H0090
19211	KMH21H0097	(01)04026704347995(17)260728(10)KMH21H0097
19211	KMH21H0159	(01)04026704347995(17)260728(10)KMH21H0159
19211	KMH21H0161	(01)04026704347995(17)260728(10)KMH21H0161
19211	KMH21H0167	(01)04026704347995(17)260728(10)KMH21H0167
19211	KMH21H0170	(01)04026704347995(17)260728(10)KMH21H0170
19211	KMH21H0188	(01)04026704347995(17)260728(10)KMH21H0188
19211	KMH21H0200	(01)04026704347995(17)260728(10)KMH21H0200
19211	KMH21H0208	(01)04026704347995(17)260728(10)KMH21H0208
19211	KMH21H0210	(01)04026704347995(17)260728(10)KMH21H0210
19211	KMH21H0214	(01)04026704347995(17)260728(10)KMH21H0214
19211	KMH21H0221	(01)04026704347995(17)260728(10)KMH21H0221
19211	KMH21H0223	(01)04026704347995(17)260728(10)KMH21H0223
19211	KMH21H0228	(01)04026704347995(17)260728(10)KMH21H0228

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21H0261	(01)04026704347995(17)260728(10)KMH21H0261
19211	KMH21H0263	(01)04026704347995(17)260728(10)KMH21H0263
19211	KMH21H0265	(01)04026704347995(17)260728(10)KMH21H0265
19211	KMH21J0007	(01)04026704347995(17)260828(10)KMH21J0007
19211	KMH21J0012	(01)04026704347995(17)260828(10)KMH21J0012
19211	KMH21J0083	(01)04026704347995(17)260828(10)KMH21J0083
19211	KMH21J0088	(01)04026704347995(17)260828(10)KMH21J0088
19211	KMH21J0089	(01)04026704347995(17)260828(10)KMH21J0089
19211	KMH21J0104	(01)04026704347995(17)260828(10)KMH21J0104
19211	KMH21J0154	(01)04026704347995(17)260828(10)KMH21J0154
19211	KMH21J0157	(01)04026704347995(17)260828(10)KMH21J0157
19211	KMH21J0158	(01)04026704347995(17)260828(10)KMH21J0158
19211	KMH21J0159	(01)04026704347995(17)260828(10)KMH21J0159
19211	KMH21J0164	(01)04026704347995(17)260828(10)KMH21J0164
19211	KMH21J0165	(01)04026704347995(17)260828(10)KMH21J0165
19211	KMH21J0167	(01)04026704347995(17)260828(10)KMH21J0167
19211	KMH21J0168	(01)04026704347995(17)260828(10)KMH21J0168
19211	KMH21J0170	(01)04026704347995(17)260828(10)KMH21J0170
19211	KMH21J0171	(01)04026704347995(17)260828(10)KMH21J0171
19211	KMH21J0175	(01)04026704347995(17)260828(10)KMH21J0175
19211	KMH21J0177	(01)04026704347995(17)260828(10)KMH21J0177
19211	KMH21J0181	(01)04026704347995(17)260828(10)KMH21J0181
19211	KMH21J0186	(01)04026704347995(17)260828(10)KMH21J0186
19211	KMH21J0188	(01)04026704347995(17)260828(10)KMH21J0188
19211	KMH21J0191	(01)04026704347995(17)260828(10)KMH21J0191
19211	KMH21J0192	(01)04026704347995(17)260828(10)KMH21J0192
19211	KMH21J0195	(01)04026704347995(17)260828(10)KMH21J0195
19211	KMH21J0199	(01)04026704347995(17)260828(10)KMH21J0199
19211	KMH21J0201	(01)04026704347995(17)260828(10)KMH21J0201
19211	KMH21J0204	(01)04026704347995(17)260828(10)KMH21J0204
19211	KMH21J0207	(01)04026704347995(17)260828(10)KMH21J0207
19211	KMH21J0212	(01)04026704347995(17)260828(10)KMH21J0212
19211	KMH21J0214	(01)04026704347995(17)260828(10)KMH21J0214
19211	KMH21J0218	(01)04026704347995(17)260828(10)KMH21J0218
19211	KMH21J0222	(01)04026704347995(17)260828(10)KMH21J0222
19211	KMH21J0226	(01)04026704347995(17)260828(10)KMH21J0226
19211	KMH21J0229	(01)04026704347995(17)260828(10)KMH21J0229
19211	KMH21J0232	(01)04026704347995(17)260828(10)KMH21J0232
19211	KMH21J0239	(01)04026704347995(17)260828(10)KMH21J0239
19211	KMH21J0240	(01)04026704347995(17)260828(10)KMH21J0240
19211	KMH21J0248	(01)04026704347995(17)260828(10)KMH21J0248
19211	KMH21J0255	(01)04026704347995(17)260828(10)KMH21J0255
19211	KMH21J0256	(01)04026704347995(17)260828(10)KMH21J0256
19211	KMH21J0262	(01)04026704347995(17)260828(10)KMH21J0262
19211	KMH21J0265	(01)04026704347995(17)260828(10)KMH21J0265
19211	KMH21J0268	(01)04026704347995(17)260828(10)KMH21J0268
19211	KMH21J0278	(01)04026704347995(17)260828(10)KMH21J0278
19211	KMH21J0281	(01)04026704347995(17)260828(10)KMH21J0281
19211	KMH21J0286	(01)04026704347995(17)260828(10)KMH21J0286
19211	KMH21J0287	(01)04026704347995(17)260828(10)KMH21J0287
19211	KMH21J0289	(01)04026704347995(17)260828(10)KMH21J0289
19211	KMH21J0291	(01)04026704347995(17)260828(10)KMH21J0291

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21J0294	(01)04026704347995(17)260828(10)KMH21J0294
19211	KMH21J0305	(01)04026704347995(17)260828(10)KMH21J0305
19211	KMH21J0306	(01)04026704347995(17)260828(10)KMH21J0306
19211	KMH21J0312	(01)04026704347995(17)260828(10)KMH21J0312
19211	KMH21J0314	(01)04026704347995(17)260828(10)KMH21J0314
19211	KMH21J0317	(01)04026704347995(17)260828(10)KMH21J0317
19211	KMH21J0328	(01)04026704347995(17)260828(10)KMH21J0328
19211	KMH21J0333	(01)04026704347995(17)260828(10)KMH21J0333
19211	KMH21J0334	(01)04026704347995(17)260828(10)KMH21J0334
19211	KMH21K0001	(01)04026704347995(17)260928(10)KMH21K0001
19211	KMH21K0003	(01)04026704347995(17)260928(10)KMH21K0003
19211	KMH21K0007	(01)04026704347995(17)260928(10)KMH21K0007
19211	KMH21K0011	(01)04026704347995(17)260928(10)KMH21K0011
19211	KMH21K0015	(01)04026704347995(17)260928(10)KMH21K0015
19211	KMH21K0018	(01)04026704347995(17)260928(10)KMH21K0018
19211	KMH21K0019	(01)04026704347995(17)260928(10)KMH21K0019
19211	KMH21K0029	(01)04026704347995(17)260928(10)KMH21K0029
19211	KMH21K0030	(01)04026704347995(17)260928(10)KMH21K0030
19211	KMH21K0037	(01)04026704347995(17)260928(10)KMH21K0037
19211	KMH21K0078	(01)04026704347995(17)260928(10)KMH21K0078
19211	KMH21K0080	(01)04026704347995(17)260928(10)KMH21K0080
19211	KMH21K0083	(01)04026704347995(17)260928(10)KMH21K0083
19211	KMH21K0092	(01)04026704347995(17)260928(10)KMH21K0092
19211	KMH21K0100	(01)04026704347995(17)260928(10)KMH21K0100
19211	KMH21K0112	(01)04026704347995(17)260928(10)KMH21K0112
19211	KMH21K0182	(01)04026704347995(17)260928(10)KMH21K0182
19211	KMH21K0186	(01)04026704347995(17)260928(10)KMH21K0186
19211	KMH21K0237	(01)04026704347995(17)260928(10)KMH21K0237
19211	KMH21K0401	(01)04026704347995(17)260928(10)KMH21K0401
19211	KMH21K0403	(01)04026704347995(17)260928(10)KMH21K0403
19211	KMH21K0406	(01)04026704347995(17)260928(10)KMH21K0406
19211	KMH21K0407	(01)04026704347995(17)260928(10)KMH21K0407
19211	KMH21K0408	(01)04026704347995(17)260928(10)KMH21K0408
19211	KMH21K0409	(01)04026704347995(17)260928(10)KMH21K0409
19211	KMH21L0012	(01)04026704347995(17)261028(10)KMH21L0012
19211	KMH21L0021	(01)04026704347995(17)261028(10)KMH21L0021
19211	KMH21L0023	(01)04026704347995(17)261028(10)KMH21L0023
19211	KMH21L0024	(01)04026704347995(17)261028(10)KMH21L0024
19211	KMH21L0061	(01)04026704347995(17)261028(10)KMH21L0061
19211	KMH21L0169	(01)04026704347995(17)261028(10)KMH21L0169
19211	KMH21L0178	(01)04026704347995(17)261028(10)KMH21L0178
19211	KMH21L0182	(01)04026704347995(17)261028(10)KMH21L0182
19211	KMH21L0190	(01)04026704347995(17)261028(10)KMH21L0190
19211	KMH21L0196	(01)04026704347995(17)261028(10)KMH21L0196
19211	KMH21L0211	(01)04026704347995(17)261028(10)KMH21L0211
19211	KMH21L0213	(01)04026704347995(17)261028(10)KMH21L0213
19211	KMH21L0225	(01)04026704347995(17)261028(10)KMH21L0225
19211	KMH21L0370	(01)04026704347995(17)261028(10)KMH21L0370
19211	KMH21M0002	(01)04026704347995(17)261128(10)KMH21M0002
19211	KMH21M0008	(01)04026704347995(17)261128(10)KMH21M0008
19211	KMH21M0014	(01)04026704347995(17)261128(10)KMH21M0014
19211	KMH21M0015	(01)04026704347995(17)261128(10)KMH21M0015

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH21M0016	(01)04026704347995(17)261128(10)KMH21M0016
19211	KMH21M0018	(01)04026704347995(17)261128(10)KMH21M0018
19211	KMH21M0019	(01)04026704347995(17)261128(10)KMH21M0019
19211	KMH21M0020	(01)04026704347995(17)261128(10)KMH21M0020
19211	KMH21M0031	(01)04026704347995(17)261128(10)KMH21M0031
19211	KMH21M0054	(01)04026704347995(17)261128(10)KMH21M0054
19211	KMH21M0055	(01)04026704347995(17)261128(10)KMH21M0055
19211	KMH21M0056	(01)04026704347995(17)261128(10)KMH21M0056
19211	KMH21M0057	(01)04026704347995(17)261128(10)KMH21M0057
19211	KMH21M0065	(01)04026704347995(17)261128(10)KMH21M0065
19211	KMH21M0069	(01)04026704347995(17)261128(10)KMH21M0069
19211	KMH21M0071	(01)04026704347995(17)261128(10)KMH21M0071
19211	KMH21M0082	(01)04026704347995(17)261128(10)KMH21M0082
19211	KMH21M0103	(01)04026704347995(17)261128(10)KMH21M0103
19211	KMH21M0106	(01)04026704347995(17)261128(10)KMH21M0106
19211	KMH21M0111	(01)04026704347995(17)261128(10)KMH21M0111
19211	KMH21M0135	(01)04026704347995(17)261128(10)KMH21M0135
19211	KMH21M0144	(01)04026704347995(17)261128(10)KMH21M0144
19211	KMH21M0145	(01)04026704347995(17)261128(10)KMH21M0145
19211	KMH21M0150	(01)04026704347995(17)261128(10)KMH21M0150
19211	KMH21M0157	(01)04026704347995(17)261128(10)KMH21M0157
19211	KMH21M0160	(01)04026704347995(17)261128(10)KMH21M0160
19211	KMH21M0162	(01)04026704347995(17)261128(10)KMH21M0162
19211	KMH21M0183	(01)04026704347995(17)261128(10)KMH21M0183
19211	KMH21M0185	(01)04026704347995(17)261128(10)KMH21M0185
19211	KMH21M0198	(01)04026704347995(17)261128(10)KMH21M0198
19211	KMH21M0202	(01)04026704347995(17)261128(10)KMH21M0202
19211	KMH21M0204	(01)04026704347995(17)261128(10)KMH21M0204
19211	KMH21M0206	(01)04026704347995(17)261128(10)KMH21M0206
19211	KMH21M0208	(01)04026704347995(17)261128(10)KMH21M0208
19211	KMH21M0215	(01)04026704347995(17)261128(10)KMH21M0215
19211	KMH21M0216	(01)04026704347995(17)261128(10)KMH21M0216
19211	KMH21M0226	(01)04026704347995(17)261128(10)KMH21M0226
19211	KMH21M0228	(01)04026704347995(17)261128(10)KMH21M0228
19211	KMH21M0231	(01)04026704347995(17)261128(10)KMH21M0231
19211	KMH21M0232	(01)04026704347995(17)261128(10)KMH21M0232
19211	KMH21M0242	(01)04026704347995(17)261128(10)KMH21M0242
19211	KMH21M0243	(01)04026704347995(17)261128(10)KMH21M0243
19211	KMH21M0244	(01)04026704347995(17)261128(10)KMH21M0244
19211	KMH21M0249	(01)04026704347995(17)261128(10)KMH21M0249
19211	KMH21M0250	(01)04026704347995(17)261128(10)KMH21M0250
19211	KMH22A0005	(01)04026704347995(17)261228(10)KMH22A0005
19211	KMH22A0013	(01)04026704347995(17)261228(10)KMH22A0013
19211	KMH22A0043	(01)04026704347995(17)261228(10)KMH22A0043
19211	KMH22A0046	(01)04026704347995(17)261228(10)KMH22A0046
19211	KMH22A0047	(01)04026704347995(17)261228(10)KMH22A0047
19211	KMH22A0051	(01)04026704347995(17)261228(10)KMH22A0051
19211	KMH22A0070	(01)04026704347995(17)261228(10)KMH22A0070
19211	KMH22A0078	(01)04026704347995(17)261228(10)KMH22A0078
19211	KMH22A0081	(01)04026704347995(17)261228(10)KMH22A0081
19211	KMH22A0082	(01)04026704347995(17)261228(10)KMH22A0082
19211	KMH22A0095	(01)04026704347995(17)261228(10)KMH22A0095

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH22A0096	(01)04026704347995(17)261228(10)KMH22A0096
19211	KMH22A0097	(01)04026704347995(17)261228(10)KMH22A0097
19211	KMH22A0113	(01)04026704347995(17)261228(10)KMH22A0113
19211	KMH22A0124	(01)04026704347995(17)261228(10)KMH22A0124
19211	KMH22A0126	(01)04026704347995(17)261228(10)KMH22A0126
19211	KMH22A0129	(01)04026704347995(17)261228(10)KMH22A0129
19211	KMH22A0142	(01)04026704347995(17)261228(10)KMH22A0142
19211	KMH22A0144	(01)04026704347995(17)261228(10)KMH22A0144
19211	KMH22A0146	(01)04026704347995(17)261228(10)KMH22A0146
19211	KMH22A0151	(01)04026704347995(17)261228(10)KMH22A0151
19211	KMH22A0155	(01)04026704347995(17)261228(10)KMH22A0155
19211	KMH22A0270	(01)04026704347995(17)261228(10)KMH22A0270
19211	KMH22A0281	(01)04026704347995(17)261228(10)KMH22A0281
19211	KMH22A0283	(01)04026704347995(17)261228(10)KMH22A0283
19211	KMH22A0286	(01)04026704347995(17)261228(10)KMH22A0286
19211	KMH22A0297	(01)04026704347995(17)261228(10)KMH22A0297
19211	KMH22A0300	(01)04026704347995(17)261228(10)KMH22A0300
19211	KMH22A0319	(01)04026704347995(17)261228(10)KMH22A0319
19211	KMH22A0326	(01)04026704347995(17)261228(10)KMH22A0326
19211	KMH22A0327	(01)04026704347995(17)261228(10)KMH22A0327
19211	KMH22A0328	(01)04026704347995(17)261228(10)KMH22A0328
19211	KMH22A0332	(01)04026704347995(17)261228(10)KMH22A0332
19211	KMH22B0001	(01)04026704347995(17)270128(10)KMH22B0001
19211	KMH22B0003	(01)04026704347995(17)270128(10)KMH22B0003
19211	KMH22B0006	(01)04026704347995(17)270128(10)KMH22B0006
19211	KMH22B0092	(01)04026704347995(17)270128(10)KMH22B0092
19211	KMH22B0094	(01)04026704347995(17)270128(10)KMH22B0094
19211	KMH22B0095	(01)04026704347995(17)270128(10)KMH22B0095
19211	KMH22B0096	(01)04026704347995(17)270128(10)KMH22B0096
19211	KMH22B0112	(01)04026704347995(17)270128(10)KMH22B0112
19211	KMH22B0113	(01)04026704347995(17)270128(10)KMH22B0113
19211	KMH22B0120	(01)04026704347995(17)270128(10)KMH22B0120
19211	KMH22B0122	(01)04026704347995(17)270128(10)KMH22B0122
19211	KMH22B0123	(01)04026704347995(17)270128(10)KMH22B0123
19211	KMH22B0124	(01)04026704347995(17)270128(10)KMH22B0124
19211	KMH22B0125	(01)04026704347995(17)270128(10)KMH22B0125
19211	KMH22B0127	(01)04026704347995(17)270128(10)KMH22B0127
19211	KMH22B0128	(01)04026704347995(17)270128(10)KMH22B0128
19211	KMH22B0129	(01)04026704347995(17)270128(10)KMH22B0129
19211	KMH22B0130	(01)04026704347995(17)270128(10)KMH22B0130
19211	KMH22B0133	(01)04026704347995(17)270128(10)KMH22B0133
19211	KMH22B0170	(01)04026704347995(17)270128(10)KMH22B0170
19211	KMH22B0171	(01)04026704347995(17)270128(10)KMH22B0171
19211	KMH22B0176	(01)04026704347995(17)270128(10)KMH22B0176
19211	KMH22B0184	(01)04026704347995(17)270128(10)KMH22B0184
19211	KMH22B0186	(01)04026704347995(17)270128(10)KMH22B0186
19211	KMH22B0194	(01)04026704347995(17)270128(10)KMH22B0194
19211	KMH22B0199	(01)04026704347995(17)270128(10)KMH22B0199
19211	KMH22B0206	(01)04026704347995(17)270128(10)KMH22B0206
19211	KMH22B0210	(01)04026704347995(17)270128(10)KMH22B0210
19211	KMH22B0214	(01)04026704347995(17)270128(10)KMH22B0214
19211	KMH22B0236	(01)04026704347995(17)270128(10)KMH22B0236

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19211	KMH22C0058	(01)04026704347995(17)270228(10)KMH22C0058
19211	KMH22C0063	(01)04026704347995(17)270228(10)KMH22C0063
19211	KMH22C0064	(01)04026704347995(17)270228(10)KMH22C0064
19211	KMH22C0079	(01)04026704347995(17)270228(10)KMH22C0079
19211	KMH22C0080	(01)04026704347995(17)270228(10)KMH22C0080
19211	KMH22C0083	(01)04026704347995(17)270228(10)KMH22C0083
19211	KMH22C0087	(01)04026704347995(17)270228(10)KMH22C0087
19211	KMH22C0095	(01)04026704347995(17)270228(10)KMH22C0095
19211	KMH22C0102	(01)04026704347995(17)270228(10)KMH22C0102
19211	KMH22C0103	(01)04026704347995(17)270228(10)KMH22C0103
19211	KMH22C0134	(01)04026704347995(17)270228(10)KMH22C0134
19211	KMH22C0147	(01)04026704347995(17)270228(10)KMH22C0147
19211	KMH22C0150	(01)04026704347995(17)270228(10)KMH22C0150
19211T	20FT38	(01)04026704400478(17)250528(10)20FT38
19211T	20FT63	(01)04026704400478(17)250528(10)20FT63
19211T	20GT03	(01)04026704400478(17)250628(10)20GT03
19211T	KMH20H0303	(01)04026704400478(17)250728(10)KMH20H0303
19211T	KMH20J0148	(01)04026704400478(17)250828(10)KMH20J0148
19211T	KMH20L0135	(01)04026704400478(17)251028(10)KMH20L0135
19211T	KMH20L0137	(01)04026704400478(17)251028(10)KMH20L0137
19211T	KMH20L0362	(01)04026704400478(17)251028(10)KMH20L0362
19211T	KMH20M0236	(01)04026704400478(17)251128(10)KMH20M0236
19211T	KMH20M0276	(01)04026704400478(17)251128(10)KMH20M0276
19211T	KMH20M0477	(01)04026704400478(17)251128(10)KMH20M0477
19211T	KMH21A0397	(01)04026704400478(17)251228(10)KMH21A0397
19211T	KMH21C0228	(01)04026704400478(17)260228(10)KMH21C0228
19211T	KMH21J0146	(01)04026704400478(17)260828(10)KMH21J0146
19211T	KMH21K0122	(01)04026704400478(17)260928(10)KMH21K0122
19211T	KMH21K0130	(01)04026704400478(17)260928(10)KMH21K0130
19211T	KMH21K0133	(01)04026704400478(17)260928(10)KMH21K0133
19211T	KMH21K0137	(01)04026704400478(17)260928(10)KMH21K0137
19211T	KMH21K0155	(01)04026704400478(17)260928(10)KMH21K0155
19212	202023	(01)04026704348008(17)250528(10)202023
19212	202024	(01)04026704348008(17)250528(10)202024
19212	202025	(01)04026704348008(17)250528(10)202025
19212	202026	(01)04026704348008(17)250528(10)202026
19212	202027	(01)04026704348008(17)250628(10)202027
19212	202028	(01)04026704348008(17)250628(10)202028
19212	202029	(01)04026704348008(17)250628(10)202029
19212	202030	(01)04026704348008(17)250628(10)202030
19212	KMZ20G0098	(01)04026704348008(17)250628(10)KMZ20G0098
19212	KMZ20G0120	(01)04026704348008(17)250628(10)KMZ20G0120
19212	KMZ20G0121	(01)04026704348008(17)250628(10)KMZ20G0121
19212	KMZ20G0133	(01)04026704348008(17)250628(10)KMZ20G0133
19212	KMZ20G0147	(01)04026704348008(17)250628(10)KMZ20G0147
19212	KMZ20G0149	(01)04026704348008(17)250628(10)KMZ20G0149
19212	KMZ20G0151	(01)04026704348008(17)250628(10)KMZ20G0151
19212	KMZ20G0181	(01)04026704348008(17)250628(10)KMZ20G0181
19212	KMZ20H0138	(01)04026704348008(17)250728(10)KMZ20H0138
19212	KMZ20H0148	(01)04026704348008(17)250728(10)KMZ20H0148
19212	KMZ20H0213	(01)04026704348008(17)250728(10)KMZ20H0213
19212	KMZ20H0238	(01)04026704348008(17)250728(10)KMZ20H0238

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212	KMZ20H0281	(01)04026704348008(17)250728(10)KMZ20H0281
19212	KMZ20H0282	(01)04026704348008(17)250728(10)KMZ20H0282
19212	KMZ20H0308	(01)04026704348008(17)250728(10)KMZ20H0308
19212	KMZ20H0384	(01)04026704348008(17)250728(10)KMZ20H0384
19212	KMZ20H0387	(01)04026704348008(17)250728(10)KMZ20H0387
19212	KMZ20H0401	(01)04026704348008(17)250728(10)KMZ20H0401
19212	KMZ20H0411	(01)04026704348008(17)250728(10)KMZ20H0411
19212	KMZ20H0426	(01)04026704348008(17)250728(10)KMZ20H0426
19212	KMZ20H0428	(01)04026704348008(17)250728(10)KMZ20H0428
19212	KMZ20H0436	(01)04026704348008(17)250728(10)KMZ20H0436
19212	KMZ20H0490	(01)04026704348008(17)250728(10)KMZ20H0490
19212	KMZ20H0495	(01)04026704348008(17)250728(10)KMZ20H0495
19212	KMZ20H0497	(01)04026704348008(17)250728(10)KMZ20H0497
19212	KMZ20H0503	(01)04026704348008(17)250728(10)KMZ20H0503
19212	KMZ20H0517	(01)04026704348008(17)250728(10)KMZ20H0517
19212	KMZ20H0562	(01)04026704348008(17)250728(10)KMZ20H0562
19212	KMZ20H0595	(01)04026704348008(17)250728(10)KMZ20H0595
19212	KMZ20H0606	(01)04026704348008(17)250728(10)KMZ20H0606
19212	KMZ20H0618	(01)04026704348008(17)250728(10)KMZ20H0618
19212	KMZ20H0666	(01)04026704348008(17)250728(10)KMZ20H0666
19212	KMZ20H0677	(01)04026704348008(17)250728(10)KMZ20H0677
19212	KMZ20H0678	(01)04026704348008(17)250728(10)KMZ20H0678
19212	KMZ20H0687	(01)04026704348008(17)250728(10)KMZ20H0687
19212	KMZ20H0691	(01)04026704348008(17)250728(10)KMZ20H0691
19212	KMZ20H0750	(01)04026704348008(17)250728(10)KMZ20H0750
19212	KMZ20H0751	(01)04026704348008(17)250728(10)KMZ20H0751
19212	KMZ20H0767	(01)04026704348008(17)250728(10)KMZ20H0767
19212	KMZ20H0771	(01)04026704348008(17)250728(10)KMZ20H0771
19212	KMZ20H0850	(01)04026704348008(17)250728(10)KMZ20H0850
19212	KMZ20H0852	(01)04026704348008(17)250728(10)KMZ20H0852
19212	KMZ20H0864	(01)04026704348008(17)250728(10)KMZ20H0864
19212	KMZ20H0867	(01)04026704348008(17)250728(10)KMZ20H0867
19212	KMZ20H0882	(01)04026704348008(17)250728(10)KMZ20H0882
19212	KMZ20H0946	(01)04026704348008(17)250728(10)KMZ20H0946
19212	KMZ20H0955	(01)04026704348008(17)250728(10)KMZ20H0955
19212	KMZ20H0972	(01)04026704348008(17)250728(10)KMZ20H0972
19212	KMZ20H0976	(01)04026704348008(17)250728(10)KMZ20H0976
19212	KMZ20H0991	(01)04026704348008(17)250728(10)KMZ20H0991
19212	KMZ20H0992	(01)04026704348008(17)250728(10)KMZ20H0992
19212	KMZ20H1042	(01)04026704348008(17)250728(10)KMZ20H1042
19212	KMZ20H1043	(01)04026704348008(17)250728(10)KMZ20H1043
19212	KMZ20J0167	(01)04026704348008(17)250828(10)KMZ20J0167
19212	KMZ20J0215	(01)04026704348008(17)250828(10)KMZ20J0215
19212	KMZ20J0601	(01)04026704348008(17)250828(10)KMZ20J0601
19212	KMZ20J0944	(01)04026704348008(17)250828(10)KMZ20J0944
19212	KMZ20K0322	(01)04026704348008(17)250928(10)KMZ20K0322
19212	KMZ20K0415	(01)04026704348008(17)250928(10)KMZ20K0415
19212	KMZ20K0444	(01)04026704348008(17)250928(10)KMZ20K0444
19212	KMZ20K0459	(01)04026704348008(17)250928(10)KMZ20K0459
19212	KMZ20K0465	(01)04026704348008(17)250928(10)KMZ20K0465
19212	KMZ20K0468	(01)04026704348008(17)250928(10)KMZ20K0468
19212	KMZ20K0517	(01)04026704348008(17)250928(10)KMZ20K0517

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212	KMZ20K0531	(01)04026704348008(17)250928(10)KMZ20K0531
19212	KMZ20K0538	(01)04026704348008(17)250928(10)KMZ20K0538
19212	KMZ20K0591	(01)04026704348008(17)250928(10)KMZ20K0591
19212	KMZ20K0594	(01)04026704348008(17)250928(10)KMZ20K0594
19212	KMZ20K0595	(01)04026704348008(17)250928(10)KMZ20K0595
19212	KMZ20K0614	(01)04026704348008(17)250928(10)KMZ20K0614
19212	KMZ20K0632	(01)04026704348008(17)250928(10)KMZ20K0632
19212	KMZ20K0638	(01)04026704348008(17)250928(10)KMZ20K0638
19212	KMZ20K0641	(01)04026704348008(17)250928(10)KMZ20K0641
19212	KMZ20K0653	(01)04026704348008(17)250928(10)KMZ20K0653
19212	KMZ20K0674	(01)04026704348008(17)250928(10)KMZ20K0674
19212	KMZ20K0679	(01)04026704348008(17)250928(10)KMZ20K0679
19212	KMZ20K0686	(01)04026704348008(17)250928(10)KMZ20K0686
19212	KMZ20K0725	(01)04026704348008(17)250928(10)KMZ20K0725
19212	KMZ20K0805	(01)04026704348008(17)250928(10)KMZ20K0805
19212	KMZ20K1020	(01)04026704348008(17)250928(10)KMZ20K1020
19212	KMZ20K1045	(01)04026704348008(17)250928(10)KMZ20K1045
19212	KMZ20K1054	(01)04026704348008(17)250928(10)KMZ20K1054
19212	KMZ20K1057	(01)04026704348008(17)250928(10)KMZ20K1057
19212	KMZ20K1065	(01)04026704348008(17)250928(10)KMZ20K1065
19212	KMZ20K1179	(01)04026704348008(17)250928(10)KMZ20K1179
19212	KMZ20K1320	(01)04026704348008(17)250928(10)KMZ20K1320
19212	KMZ20L0019	(01)04026704348008(17)251028(10)KMZ20L0019
19212	KMZ20L0034	(01)04026704348008(17)251028(10)KMZ20L0034
19212	KMZ20L0035	(01)04026704348008(17)251028(10)KMZ20L0035
19212	KMZ20L0106	(01)04026704348008(17)251028(10)KMZ20L0106
19212	KMZ20L0119	(01)04026704348008(17)251028(10)KMZ20L0119
19212	KMZ20L0166	(01)04026704348008(17)251028(10)KMZ20L0166
19212	KMZ20L0167	(01)04026704348008(17)251028(10)KMZ20L0167
19212	KMZ20L0188	(01)04026704348008(17)251028(10)KMZ20L0188
19212	KMZ20L0208	(01)04026704348008(17)251028(10)KMZ20L0208
19212	KMZ20L0209	(01)04026704348008(17)251028(10)KMZ20L0209
19212	KMZ20L0226	(01)04026704348008(17)251028(10)KMZ20L0226
19212	KMZ20L0241	(01)04026704348008(17)251028(10)KMZ20L0241
19212	KMZ20L0266	(01)04026704348008(17)251028(10)KMZ20L0266
19212	KMZ20L0309	(01)04026704348008(17)251028(10)KMZ20L0309
19212	KMZ20L0319	(01)04026704348008(17)251028(10)KMZ20L0319
19212	KMZ20L0361	(01)04026704348008(17)251028(10)KMZ20L0361
19212	KMZ20L0377	(01)04026704348008(17)251028(10)KMZ20L0377
19212	KMZ20L0392	(01)04026704348008(17)251028(10)KMZ20L0392
19212	KMZ20L0393	(01)04026704348008(17)251028(10)KMZ20L0393
19212	KMZ20L0404	(01)04026704348008(17)251028(10)KMZ20L0404
19212	KMZ20L0434	(01)04026704348008(17)251028(10)KMZ20L0434
19212	KMZ20L0441	(01)04026704348008(17)251028(10)KMZ20L0441
19212	KMZ20L0461	(01)04026704348008(17)251028(10)KMZ20L0461
19212	KMZ20L0475	(01)04026704348008(17)251028(10)KMZ20L0475
19212	KMZ20L0517	(01)04026704348008(17)251028(10)KMZ20L0517
19212	KMZ20L0522	(01)04026704348008(17)251028(10)KMZ20L0522
19212	KMZ20L0550	(01)04026704348008(17)251028(10)KMZ20L0550
19212	KMZ20L0618	(01)04026704348008(17)251028(10)KMZ20L0618
19212	KMZ20L0622	(01)04026704348008(17)251028(10)KMZ20L0622
19212	KMZ20L0637	(01)04026704348008(17)251028(10)KMZ20L0637



**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212	KMZ20L1178	(01)04026704348008(17)251028(10)KMZ20L1178
19212	KMZ20M0251	(01)04026704348008(17)251128(10)KMZ20M0251
19212	KMZ20M0432	(01)04026704348008(17)251128(10)KMZ20M0432
19212	KMZ20M0456	(01)04026704348008(17)251128(10)KMZ20M0456
19212	KMZ20M0558	(01)04026704348008(17)251128(10)KMZ20M0558
19212	KMZ20M0849	(01)04026704348008(17)251128(10)KMZ20M0849
19212	KMZ20M0864	(01)04026704348008(17)251128(10)KMZ20M0864
19212	KMZ20M0878	(01)04026704348008(17)251128(10)KMZ20M0878
19212	KMZ20M0895	(01)04026704348008(17)251128(10)KMZ20M0895
19212	KMZ20M0904	(01)04026704348008(17)251128(10)KMZ20M0904
19212	KMZ20M0946	(01)04026704348008(17)251128(10)KMZ20M0946
19212	KMZ20M1159	(01)04026704348008(17)251128(10)KMZ20M1159
19212	KMZ20M1195	(01)04026704348008(17)251128(10)KMZ20M1195
19212	KMZ20M1204	(01)04026704348008(17)251128(10)KMZ20M1204
19212	KMZ20M1243	(01)04026704348008(17)251128(10)KMZ20M1243
19212	KMZ20M1255	(01)04026704348008(17)251128(10)KMZ20M1255
19212	KMZ20M1271	(01)04026704348008(17)251128(10)KMZ20M1271
19212	KMZ21A0489	(01)04026704348008(17)251228(10)KMZ21A0489
19212	KMZ21A0538	(01)04026704348008(17)251228(10)KMZ21A0538
19212	KMZ21A0587	(01)04026704348008(17)251228(10)KMZ21A0587
19212	KMZ21A0622	(01)04026704348008(17)251228(10)KMZ21A0622
19212	KMZ21A0633	(01)04026704348008(17)251228(10)KMZ21A0633
19212	KMZ21A0644	(01)04026704348008(17)251228(10)KMZ21A0644
19212	KMZ21A0652	(01)04026704348008(17)251228(10)KMZ21A0652
19212	KMZ21A0698	(01)04026704348008(17)251228(10)KMZ21A0698
19212	KMZ21A0707	(01)04026704348008(17)251228(10)KMZ21A0707
19212	KMZ21A0709	(01)04026704348008(17)251228(10)KMZ21A0709
19212	KMZ21A0726	(01)04026704348008(17)251228(10)KMZ21A0726
19212	KMZ21A0778	(01)04026704348008(17)251228(10)KMZ21A0778
19212	KMZ21A1007	(01)04026704348008(17)251228(10)KMZ21A1007
19212	KMZ21A1057	(01)04026704348008(17)251228(10)KMZ21A1057
19212	KMZ21A1086	(01)04026704348008(17)251228(10)KMZ21A1086
19212	KMZ21A1100	(01)04026704348008(17)251228(10)KMZ21A1100
19212	KMZ21A1104	(01)04026704348008(17)251228(10)KMZ21A1104
19212	KMZ21A1121	(01)04026704348008(17)251228(10)KMZ21A1121
19212	KMZ21A1126	(01)04026704348008(17)251228(10)KMZ21A1126
19212	KMZ21A1127	(01)04026704348008(17)251228(10)KMZ21A1127
19212	KMZ21B0054	(01)04026704348008(17)260128(10)KMZ21B0054
19212	KMZ21B0338	(01)04026704348008(17)260128(10)KMZ21B0338
19212	KMZ21B0519	(01)04026704348008(17)260128(10)KMZ21B0519
19212	KMZ21B0543	(01)04026704348008(17)260128(10)KMZ21B0543
19212	KMZ21B0548	(01)04026704348008(17)260128(10)KMZ21B0548
19212	KMZ21B0565	(01)04026704348008(17)260128(10)KMZ21B0565
19212	KMZ21B0612	(01)04026704348008(17)260128(10)KMZ21B0612
19212	KMZ21B0780	(01)04026704348008(17)260128(10)KMZ21B0780
19212	KMZ21B0816	(01)04026704348008(17)260128(10)KMZ21B0816
19212	KMZ21B1062	(01)04026704348008(17)260128(10)KMZ21B1062
19212	KMZ21B1071	(01)04026704348008(17)260128(10)KMZ21B1071
19212	KMZ21B1092	(01)04026704348008(17)260128(10)KMZ21B1092
19212	KMZ21B1098	(01)04026704348008(17)260128(10)KMZ21B1098
19212	KMZ21B1131	(01)04026704348008(17)260128(10)KMZ21B1131
19212	KMZ21B1183	(01)04026704348008(17)260128(10)KMZ21B1183

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212	KMZ21B1196	(01)04026704348008(17)260128(10)KMZ21B1196
19212	KMZ21C0027	(01)04026704348008(17)260228(10)KMZ21C0027
19212	KMZ21C0059	(01)04026704348008(17)260228(10)KMZ21C0059
19212	KMZ21C0564	(01)04026704348008(17)260228(10)KMZ21C0564
19212	KMZ21C0870	(01)04026704348008(17)260228(10)KMZ21C0870
19212	KMZ21C0908	(01)04026704348008(17)260228(10)KMZ21C0908
19212	KMZ21C1208	(01)04026704348008(17)260228(10)KMZ21C1208
19212	KMZ21D0018	(01)04026704348008(17)260328(10)KMZ21D0018
19212	KMZ21D0179	(01)04026704348008(17)260328(10)KMZ21D0179
19212	KMZ21D0195	(01)04026704348008(17)260328(10)KMZ21D0195
19212	KMZ21D0462	(01)04026704348008(17)260328(10)KMZ21D0462
19212	KMZ21D0937	(01)04026704348008(17)260328(10)KMZ21D0937
19212	KMZ21D0942	(01)04026704348008(17)260328(10)KMZ21D0942
19212	KMZ21D0970	(01)04026704348008(17)260328(10)KMZ21D0970
19212	KMZ21D0980	(01)04026704348008(17)260328(10)KMZ21D0980
19212	KMZ21D0994	(01)04026704348008(17)260328(10)KMZ21D0994
19212	KMZ21D0995	(01)04026704348008(17)260328(10)KMZ21D0995
19212	KMZ21D1027	(01)04026704348008(17)260328(10)KMZ21D1027
19212	KMZ21D1059	(01)04026704348008(17)260328(10)KMZ21D1059
19212	KMZ21D1066	(01)04026704348008(17)260328(10)KMZ21D1066
19212	KMZ21D1078	(01)04026704348008(17)260328(10)KMZ21D1078
19212	KMZ21D1092	(01)04026704348008(17)260328(10)KMZ21D1092
19212	KMZ21D1114	(01)04026704348008(17)260328(10)KMZ21D1114
19212	KMZ21E0016	(01)04026704348008(17)260428(10)KMZ21E0016
19212	KMZ21E0021	(01)04026704348008(17)260428(10)KMZ21E0021
19212	KMZ21E0047	(01)04026704348008(17)260428(10)KMZ21E0047
19212	KMZ21E0054	(01)04026704348008(17)260428(10)KMZ21E0054
19212	KMZ21E0075	(01)04026704348008(17)260428(10)KMZ21E0075
19212	KMZ21E0083	(01)04026704348008(17)260428(10)KMZ21E0083
19212	KMZ21E0099	(01)04026704348008(17)260428(10)KMZ21E0099
19212	KMZ21E0118	(01)04026704348008(17)260428(10)KMZ21E0118
19212	KMZ21E0120	(01)04026704348008(17)260428(10)KMZ21E0120
19212	KMZ21E0617	(01)04026704348008(17)260428(10)KMZ21E0617
19212	KMZ21E0629	(01)04026704348008(17)260428(10)KMZ21E0629
19212	KMZ21E0635	(01)04026704348008(17)260428(10)KMZ21E0635
19212	KMZ21E0647	(01)04026704348008(17)260428(10)KMZ21E0647
19212	KMZ21E0656	(01)04026704348008(17)260428(10)KMZ21E0656
19212	KMZ21E0662	(01)04026704348008(17)260428(10)KMZ21E0662
19212	KMZ21E0666	(01)04026704348008(17)260428(10)KMZ21E0666
19212	KMZ21F0005	(01)04026704348008(17)260528(10)KMZ21F0005
19212	KMZ21F0088	(01)04026704348008(17)260528(10)KMZ21F0088
19212	KMZ21F0102	(01)04026704348008(17)260528(10)KMZ21F0102
19212	KMZ21F0411	(01)04026704348008(17)260528(10)KMZ21F0411
19212	KMZ21F0420	(01)04026704348008(17)260528(10)KMZ21F0420
19212	KMZ21F0965	(01)04026704348008(17)260528(10)KMZ21F0965
19212	KMZ21G0003	(01)04026704348008(17)260628(10)KMZ21G0003
19212	KMZ21H0247	(01)04026704348008(17)260728(10)KMZ21H0247
19212	KMZ21H0291	(01)04026704348008(17)260728(10)KMZ21H0291
19212	KMZ21H0297	(01)04026704348008(17)260728(10)KMZ21H0297
19212	KMZ21H0304	(01)04026704348008(17)260728(10)KMZ21H0304
19212	KMZ21H0308	(01)04026704348008(17)260728(10)KMZ21H0308
19212	KMZ21H0312	(01)04026704348008(17)260728(10)KMZ21H0312

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212	KMZ21H0333	(01)04026704348008(17)260728(10)KMZ21H0333
19212	KMZ21H0338	(01)04026704348008(17)260728(10)KMZ21H0338
19212	KMZ21H0378	(01)04026704348008(17)260728(10)KMZ21H0378
19212	KMZ21H0384	(01)04026704348008(17)260728(10)KMZ21H0384
19212	KMZ21H0465	(01)04026704348008(17)260728(10)KMZ21H0465
19212	KMZ21H0468	(01)04026704348008(17)260728(10)KMZ21H0468
19212	KMZ21H0553	(01)04026704348008(17)260728(10)KMZ21H0553
19212	KMZ21H0555	(01)04026704348008(17)260728(10)KMZ21H0555
19212	KMZ21H0556	(01)04026704348008(17)260728(10)KMZ21H0556
19212	KMZ21J0003	(01)04026704348008(17)260828(10)KMZ21J0003
19212	KMZ21J0562	(01)04026704348008(17)260828(10)KMZ21J0562
19212	KMZ21J0563	(01)04026704348008(17)260828(10)KMZ21J0563
19212	KMZ21J0586	(01)04026704348008(17)260828(10)KMZ21J0586
19212	KMZ21J0587	(01)04026704348008(17)260828(10)KMZ21J0587
19212	KMZ21J0662	(01)04026704348008(17)260828(10)KMZ21J0662
19212	KMZ21J0701	(01)04026704348008(17)260828(10)KMZ21J0701
19212	KMZ21J0703	(01)04026704348008(17)260828(10)KMZ21J0703
19212	KMZ21J0708	(01)04026704348008(17)260828(10)KMZ21J0708
19212	KMZ21J0713	(01)04026704348008(17)260828(10)KMZ21J0713
19212	KMZ21J0725	(01)04026704348008(17)260828(10)KMZ21J0725
19212	KMZ21J0728	(01)04026704348008(17)260828(10)KMZ21J0728
19212	KMZ21K0046	(01)04026704348008(17)260928(10)KMZ21K0046
19212	KMZ21K0064	(01)04026704348008(17)260928(10)KMZ21K0064
19212	KMZ21K0075	(01)04026704348008(17)260928(10)KMZ21K0075
19212	KMZ21K0076	(01)04026704348008(17)260928(10)KMZ21K0076
19212	KMZ21K0084	(01)04026704348008(17)260928(10)KMZ21K0084
19212	KMZ21K0098	(01)04026704348008(17)260928(10)KMZ21K0098
19212	KMZ21K0106	(01)04026704348008(17)260928(10)KMZ21K0106
19212	KMZ21K0161	(01)04026704348008(17)260928(10)KMZ21K0161
19212	KMZ21K0182	(01)04026704348008(17)260928(10)KMZ21K0182
19212	KMZ21K0183	(01)04026704348008(17)260928(10)KMZ21K0183
19212	KMZ21K0523	(01)04026704348008(17)260928(10)KMZ21K0523
19212	KMZ21K0528	(01)04026704348008(17)260928(10)KMZ21K0528
19212	KMZ21K0530	(01)04026704348008(17)260928(10)KMZ21K0530
19212	KMZ21K0597	(01)04026704348008(17)260928(10)KMZ21K0597
19212	KMZ21K0598	(01)04026704348008(17)260928(10)KMZ21K0598
19212	KMZ21K0611	(01)04026704348008(17)260928(10)KMZ21K0611
19212	KMZ21K0613	(01)04026704348008(17)260928(10)KMZ21K0613
19212	KMZ21K0943	(01)04026704348008(17)260928(10)KMZ21K0943
19212	KMZ21K0944	(01)04026704348008(17)260928(10)KMZ21K0944
19212	KMZ21K0945	(01)04026704348008(17)260928(10)KMZ21K0945
19212	KMZ21K0946	(01)04026704348008(17)260928(10)KMZ21K0946
19212	KMZ21L0888	(01)04026704348008(17)261028(10)KMZ21L0888
19212	KMZ21L0892	(01)04026704348008(17)261028(10)KMZ21L0892
19212	KMZ21L0900	(01)04026704348008(17)261028(10)KMZ21L0900
19212	KMZ21M0255	(01)04026704348008(17)261128(10)KMZ21M0255
19212	KMZ21M0261	(01)04026704348008(17)261128(10)KMZ21M0261
19212	KMZ21M0311	(01)04026704348008(17)261128(10)KMZ21M0311
19212	KMZ21M0318	(01)04026704348008(17)261128(10)KMZ21M0318
19212	KMZ21M0332	(01)04026704348008(17)261128(10)KMZ21M0332
19212	KMZ21M0333	(01)04026704348008(17)261128(10)KMZ21M0333
19212	KMZ21M0338	(01)04026704348008(17)261128(10)KMZ21M0338

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212	KMZ21M0346	(01)04026704348008(17)261128(10)KMZ21M0346
19212	KMZ21M0355	(01)04026704348008(17)261128(10)KMZ21M0355
19212	KMZ21M0363	(01)04026704348008(17)261128(10)KMZ21M0363
19212	KMZ21M0380	(01)04026704348008(17)261128(10)KMZ21M0380
19212	KMZ21M0434	(01)04026704348008(17)261128(10)KMZ21M0434
19212	KMZ21M0448	(01)04026704348008(17)261128(10)KMZ21M0448
19212	KMZ21M0470	(01)04026704348008(17)261128(10)KMZ21M0470
19212	KMZ21M0512	(01)04026704348008(17)261128(10)KMZ21M0512
19212	KMZ21M0522	(01)04026704348008(17)261128(10)KMZ21M0522
19212	KMZ21M0530	(01)04026704348008(17)261128(10)KMZ21M0530
19212	KMZ21M0591	(01)04026704348008(17)261128(10)KMZ21M0591
19212	KMZ21M0639	(01)04026704348008(17)261128(10)KMZ21M0639
19212	KMZ21M0707	(01)04026704348008(17)261128(10)KMZ21M0707
19212	KMZ21M0724	(01)04026704348008(17)261128(10)KMZ21M0724
19212	KMZ21M0888	(01)04026704348008(17)261128(10)KMZ21M0888
19212	KMZ21M0896	(01)04026704348008(17)261128(10)KMZ21M0896
19212	KMZ21M0906	(01)04026704348008(17)261128(10)KMZ21M0906
19212	KMZ21M0914	(01)04026704348008(17)261128(10)KMZ21M0914
19212	KMZ21M0943	(01)04026704348008(17)261128(10)KMZ21M0943
19212	KMZ21M0954	(01)04026704348008(17)261128(10)KMZ21M0954
19212	KMZ21M0989	(01)04026704348008(17)261128(10)KMZ21M0989
19212	KMZ21M1006	(01)04026704348008(17)261128(10)KMZ21M1006
19212	KMZ21M1037	(01)04026704348008(17)261128(10)KMZ21M1037
19212	KMZ21M1053	(01)04026704348008(17)261128(10)KMZ21M1053
19212	KMZ22A0013	(01)04026704348008(17)261228(10)KMZ22A0013
19212	KMZ22A0016	(01)04026704348008(17)261228(10)KMZ22A0016
19212	KMZ22A0021	(01)04026704348008(17)261228(10)KMZ22A0021
19212	KMZ22A0025	(01)04026704348008(17)261228(10)KMZ22A0025
19212	KMZ22A0068	(01)04026704348008(17)261228(10)KMZ22A0068
19212	KMZ22A0087	(01)04026704348008(17)261228(10)KMZ22A0087
19212	KMZ22A0124	(01)04026704348008(17)261228(10)KMZ22A0124
19212	KMZ22A0132	(01)04026704348008(17)261228(10)KMZ22A0132
19212	KMZ22A0264	(01)04026704348008(17)261228(10)KMZ22A0264
19212	KMZ22A0274	(01)04026704348008(17)261228(10)KMZ22A0274
19212	KMZ22A0280	(01)04026704348008(17)261228(10)KMZ22A0280
19212	KMZ22A0439	(01)04026704348008(17)261228(10)KMZ22A0439
19212	KMZ22A0661	(01)04026704348008(17)261228(10)KMZ22A0661
19212	KMZ22A0662	(01)04026704348008(17)261228(10)KMZ22A0662
19212	KMZ22A0702	(01)04026704348008(17)261228(10)KMZ22A0702
19212	KMZ22A0713	(01)04026704348008(17)261228(10)KMZ22A0713
19212	KMZ22A0730	(01)04026704348008(17)261228(10)KMZ22A0730
19212	KMZ22A0744	(01)04026704348008(17)261228(10)KMZ22A0744
19212	KMZ22A0755	(01)04026704348008(17)261228(10)KMZ22A0755
19212	KMZ22A0767	(01)04026704348008(17)261228(10)KMZ22A0767
19212	KMZ22A0777	(01)04026704348008(17)261228(10)KMZ22A0777
19212	KMZ22A0904	(01)04026704348008(17)261228(10)KMZ22A0904
19212	KMZ22A1002	(01)04026704348008(17)261228(10)KMZ22A1002
19212	KMZ22A1017	(01)04026704348008(17)261228(10)KMZ22A1017
19212	KMZ22B0001	(01)04026704348008(17)270128(10)KMZ22B0001
19212	KMZ22B0016	(01)04026704348008(17)270128(10)KMZ22B0016
19212	KMZ22B0022	(01)04026704348008(17)270128(10)KMZ22B0022
19212	KMZ22B0064	(01)04026704348008(17)270128(10)KMZ22B0064

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212	KMZ22B0071	(01)04026704348008(17)270128(10)KMZ22B0071
19212	KMZ22B0080	(01)04026704348008(17)270128(10)KMZ22B0080
19212	KMZ22B0090	(01)04026704348008(17)270128(10)KMZ22B0090
19212	KMZ22B0097	(01)04026704348008(17)270128(10)KMZ22B0097
19212	KMZ22B0106	(01)04026704348008(17)270128(10)KMZ22B0106
19212	KMZ22B0263	(01)04026704348008(17)270128(10)KMZ22B0263
19212	KMZ22B0596	(01)04026704348008(17)270128(10)KMZ22B0596
19212	KMZ22B0631	(01)04026704348008(17)270128(10)KMZ22B0631
19212	KMZ22B0636	(01)04026704348008(17)270128(10)KMZ22B0636
19212	KMZ22B0651	(01)04026704348008(17)270128(10)KMZ22B0651
19212	KMZ22B0743	(01)04026704348008(17)270128(10)KMZ22B0743
19212	KMZ22B0859	(01)04026704348008(17)270128(10)KMZ22B0859
19212	KMZ22C0005	(01)04026704348008(17)270228(10)KMZ22C0005
19212	KMZ22C0117	(01)04026704348008(17)270228(10)KMZ22C0117
19212	KMZ22C0226	(01)04026704348008(17)270228(10)KMZ22C0226
19212	KMZ22C0238	(01)04026704348008(17)270228(10)KMZ22C0238
19212	KMZ22C0241	(01)04026704348008(17)270228(10)KMZ22C0241
19212	KMZ22C0328	(01)04026704348008(17)270228(10)KMZ22C0328
19212	KMZ22C0476	(01)04026704348008(17)270228(10)KMZ22C0476
19212	KMZ22C0817	(01)04026704348008(17)270228(10)KMZ22C0817
19212T	202025	(01)04026704400485(17)250528(10)202025
19212T	202026	(01)04026704400485(17)250528(10)202026
19212T	202027	(01)04026704400485(17)250628(10)202027
19212T	202029	(01)04026704400485(17)250628(10)202029
19212T	202030	(01)04026704400485(17)250628(10)202030
19212T	KMZ20H0321	(01)04026704400485(17)250728(10)KMZ20H0321
19212T	KMZ20H0341	(01)04026704400485(17)250728(10)KMZ20H0341
19212T	KMZ20H0407	(01)04026704400485(17)250728(10)KMZ20H0407
19212T	KMZ20H0519	(01)04026704400485(17)250728(10)KMZ20H0519
19212T	KMZ20H0602	(01)04026704400485(17)250728(10)KMZ20H0602
19212T	KMZ20H0619	(01)04026704400485(17)250728(10)KMZ20H0619
19212T	KMZ20H0785	(01)04026704400485(17)250728(10)KMZ20H0785
19212T	KMZ20H0792	(01)04026704400485(17)250728(10)KMZ20H0792
19212T	KMZ20H0804	(01)04026704400485(17)250728(10)KMZ20H0804
19212T	KMZ20H1025	(01)04026704400485(17)250728(10)KMZ20H1025
19212T	KMZ20H1026	(01)04026704400485(17)250728(10)KMZ20H1026
19212T	KMZ20J0386	(01)04026704400485(17)250828(10)KMZ20J0386
19212T	KMZ20J0428	(01)04026704400485(17)250828(10)KMZ20J0428
19212T	KMZ20J0466	(01)04026704400485(17)250828(10)KMZ20J0466
19212T	KMZ20J0475	(01)04026704400485(17)250828(10)KMZ20J0475
19212T	KMZ20K0161	(01)04026704400485(17)250928(10)KMZ20K0161
19212T	KMZ20K0436	(01)04026704400485(17)250928(10)KMZ20K0436
19212T	KMZ20K0478	(01)04026704400485(17)250928(10)KMZ20K0478
19212T	KMZ20K0526	(01)04026704400485(17)250928(10)KMZ20K0526
19212T	KMZ20K0527	(01)04026704400485(17)250928(10)KMZ20K0527
19212T	KMZ20K0799	(01)04026704400485(17)250928(10)KMZ20K0799
19212T	KMZ20K0917	(01)04026704400485(17)250928(10)KMZ20K0917
19212T	KMZ20K1295	(01)04026704400485(17)250928(10)KMZ20K1295
19212T	KMZ20K1311	(01)04026704400485(17)250928(10)KMZ20K1311
19212T	KMZ20L0144	(01)04026704400485(17)251028(10)KMZ20L0144
19212T	KMZ20L0160	(01)04026704400485(17)251028(10)KMZ20L0160
19212T	KMZ20L0338	(01)04026704400485(17)251028(10)KMZ20L0338

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212T	KMZ20L0350	(01)04026704400485(17)251028(10)KMZ20L0350
19212T	KMZ20L0912	(01)04026704400485(17)251028(10)KMZ20L0912
19212T	KMZ20L0926	(01)04026704400485(17)251028(10)KMZ20L0926
19212T	KMZ20L0927	(01)04026704400485(17)251028(10)KMZ20L0927
19212T	KMZ20L0928	(01)04026704400485(17)251028(10)KMZ20L0928
19212T	KMZ20L0944	(01)04026704400485(17)251028(10)KMZ20L0944
19212T	KMZ20M0411	(01)04026704400485(17)251128(10)KMZ20M0411
19212T	KMZ20M1224	(01)04026704400485(17)251128(10)KMZ20M1224
19212T	KMZ21A0011	(01)04026704400485(17)251228(10)KMZ21A0011
19212T	KMZ21A0214	(01)04026704400485(17)251228(10)KMZ21A0214
19212T	KMZ21A0301	(01)04026704400485(17)251228(10)KMZ21A0301
19212T	KMZ21A0369	(01)04026704400485(17)251228(10)KMZ21A0369
19212T	KMZ21A0419	(01)04026704400485(17)251228(10)KMZ21A0419
19212T	KMZ21A0533	(01)04026704400485(17)251228(10)KMZ21A0533
19212T	KMZ21A0807	(01)04026704400485(17)251228(10)KMZ21A0807
19212T	KMZ21A0867	(01)04026704400485(17)251228(10)KMZ21A0867
19212T	KMZ21A1017	(01)04026704400485(17)251228(10)KMZ21A1017
19212T	KMZ21B0688	(01)04026704400485(17)260128(10)KMZ21B0688
19212T	KMZ21B1032	(01)04026704400485(17)260128(10)KMZ21B1032
19212T	KMZ21B1156	(01)04026704400485(17)260128(10)KMZ21B1156
19212T	KMZ21B1188	(01)04026704400485(17)260128(10)KMZ21B1188
19212T	KMZ21C0712	(01)04026704400485(17)260228(10)KMZ21C0712
19212T	KMZ21C0760	(01)04026704400485(17)260228(10)KMZ21C0760
19212T	KMZ21C0882	(01)04026704400485(17)260228(10)KMZ21C0882
19212T	KMZ21C1163	(01)04026704400485(17)260228(10)KMZ21C1163
19212T	KMZ21C1329	(01)04026704400485(17)260228(10)KMZ21C1329
19212T	KMZ21D0390	(01)04026704400485(17)260328(10)KMZ21D0390
19212T	KMZ21D0394	(01)04026704400485(17)260328(10)KMZ21D0394
19212T	KMZ21D0400	(01)04026704400485(17)260328(10)KMZ21D0400
19212T	KMZ21D0518	(01)04026704400485(17)260328(10)KMZ21D0518
19212T	KMZ21D0534	(01)04026704400485(17)260328(10)KMZ21D0534
19212T	KMZ21D1020	(01)04026704400485(17)260328(10)KMZ21D1020
19212T	KMZ21D1051	(01)04026704400485(17)260328(10)KMZ21D1051
19212T	KMZ21E0245	(01)04026704400485(17)260428(10)KMZ21E0245
19212T	KMZ21E0253	(01)04026704400485(17)260428(10)KMZ21E0253
19212T	KMZ21E0420	(01)04026704400485(17)260428(10)KMZ21E0420
19212T	KMZ21F0202	(01)04026704400485(17)260528(10)KMZ21F0202
19212T	KMZ21F0206	(01)04026704400485(17)260528(10)KMZ21F0206
19212T	KMZ21F0217	(01)04026704400485(17)260528(10)KMZ21F0217
19212T	KMZ21F0339	(01)04026704400485(17)260528(10)KMZ21F0339
19212T	KMZ21F0438	(01)04026704400485(17)260528(10)KMZ21F0438
19212T	KMZ21F0773	(01)04026704400485(17)260528(10)KMZ21F0773
19212T	KMZ21F0786	(01)04026704400485(17)260528(10)KMZ21F0786
19212T	KMZ21G0060	(01)04026704400485(17)260628(10)KMZ21G0060
19212T	KMZ21G0353	(01)04026704400485(17)260628(10)KMZ21G0353
19212T	KMZ21G0360	(01)04026704400485(17)260628(10)KMZ21G0360
19212T	KMZ21G0392	(01)04026704400485(17)260628(10)KMZ21G0392
19212T	KMZ21H0410	(01)04026704400485(17)260728(10)KMZ21H0410
19212T	KMZ21H0424	(01)04026704400485(17)260728(10)KMZ21H0424
19212T	KMZ21H0534	(01)04026704400485(17)260728(10)KMZ21H0534
19212T	KMZ21H0538	(01)04026704400485(17)260728(10)KMZ21H0538
19212T	KMZ21J0036	(01)04026704400485(17)260828(10)KMZ21J0036

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212T	KMZ21J0321	(01)04026704400485(17)260828(10)KMZ21J0321
19212T	KMZ21J0338	(01)04026704400485(17)260828(10)KMZ21J0338
19212T	KMZ21J0403	(01)04026704400485(17)260828(10)KMZ21J0403
19212T	KMZ21J0407	(01)04026704400485(17)260828(10)KMZ21J0407
19212T	KMZ21J0428	(01)04026704400485(17)260828(10)KMZ21J0428
19212T	KMZ21J0435	(01)04026704400485(17)260828(10)KMZ21J0435
19212T	KMZ21K0028	(01)04026704400485(17)260928(10)KMZ21K0028
19212T	KMZ21K0113	(01)04026704400485(17)260928(10)KMZ21K0113
19212T	KMZ21K0191	(01)04026704400485(17)260928(10)KMZ21K0191
19212T	KMZ21K0223	(01)04026704400485(17)260928(10)KMZ21K0223
19212T	KMZ21K0230	(01)04026704400485(17)260928(10)KMZ21K0230
19212T	KMZ21K0239	(01)04026704400485(17)260928(10)KMZ21K0239
19212T	KMZ21K0260	(01)04026704400485(17)260928(10)KMZ21K0260
19212T	KMZ21K0329	(01)04026704400485(17)260928(10)KMZ21K0329
19212T	KMZ21K0381	(01)04026704400485(17)260928(10)KMZ21K0381
19212T	KMZ21K0391	(01)04026704400485(17)260928(10)KMZ21K0391
19212T	KMZ21K0403	(01)04026704400485(17)260928(10)KMZ21K0403
19212T	KMZ21K0438	(01)04026704400485(17)260928(10)KMZ21K0438
19212T	KMZ21K0457	(01)04026704400485(17)260928(10)KMZ21K0457
19212T	KMZ21K0498	(01)04026704400485(17)260928(10)KMZ21K0498
19212T	KMZ21K0552	(01)04026704400485(17)260928(10)KMZ21K0552
19212T	KMZ21K0562	(01)04026704400485(17)260928(10)KMZ21K0562
19212T	KMZ21K0579	(01)04026704400485(17)260928(10)KMZ21K0579
19212T	KMZ21K0610	(01)04026704400485(17)260928(10)KMZ21K0610
19212T	KMZ21L0114	(01)04026704400485(17)261028(10)KMZ21L0114
19212T	KMZ21L0209	(01)04026704400485(17)261028(10)KMZ21L0209
19212T	KMZ21L0321	(01)04026704400485(17)261028(10)KMZ21L0321
19212T	KMZ21L0536	(01)04026704400485(17)261028(10)KMZ21L0536
19212T	KMZ21L0665	(01)04026704400485(17)261028(10)KMZ21L0665
19212T	KMZ21M0066	(01)04026704400485(17)261128(10)KMZ21M0066
19212T	KMZ21M0476	(01)04026704400485(17)261128(10)KMZ21M0476
19212T	KMZ21M0486	(01)04026704400485(17)261128(10)KMZ21M0486
19212T	KMZ21M0526	(01)04026704400485(17)261128(10)KMZ21M0526
19212T	KMZ21M0572	(01)04026704400485(17)261128(10)KMZ21M0572
19212T	KMZ21M0586	(01)04026704400485(17)261128(10)KMZ21M0586
19212T	KMZ21M0601	(01)04026704400485(17)261128(10)KMZ21M0601
19212T	KMZ21M0606	(01)04026704400485(17)261128(10)KMZ21M0606
19212T	KMZ21M0658	(01)04026704400485(17)261128(10)KMZ21M0658
19212T	KMZ21M0742	(01)04026704400485(17)261128(10)KMZ21M0742
19212T	KMZ21M0775	(01)04026704400485(17)261128(10)KMZ21M0775
19212T	KMZ22A0035	(01)04026704400485(17)261228(10)KMZ22A0035
19212T	KMZ22A0059	(01)04026704400485(17)261228(10)KMZ22A0059
19212T	KMZ22A0076	(01)04026704400485(17)261228(10)KMZ22A0076
19212T	KMZ22A0085	(01)04026704400485(17)261228(10)KMZ22A0085
19212T	KMZ22A0471	(01)04026704400485(17)261228(10)KMZ22A0471
19212T	KMZ22A0482	(01)04026704400485(17)261228(10)KMZ22A0482
19212T	KMZ22A0483	(01)04026704400485(17)261228(10)KMZ22A0483
19212T	KMZ22B0043	(01)04026704400485(17)270128(10)KMZ22B0043
19212T	KMZ22B0044	(01)04026704400485(17)270128(10)KMZ22B0044
19212T	KMZ22B0123	(01)04026704400485(17)270128(10)KMZ22B0123
19212T	KMZ22B0439	(01)04026704400485(17)270128(10)KMZ22B0439
19212T	KMZ22B0480	(01)04026704400485(17)270128(10)KMZ22B0480

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19212T	KMZ22B0783	(01)04026704400485(17)270128(10)KMZ22B0783
19212T	KMZ22B0871	(01)04026704400485(17)270128(10)KMZ22B0871
19212T	KMZ22C0231	(01)04026704400485(17)270228(10)KMZ22C0231
19212T	KMZ22C0398	(01)04026704400485(17)270228(10)KMZ22C0398
19212T	KMZ22C0444	(01)04026704400485(17)270228(10)KMZ22C0444
19212T	KMZ22C0606	(01)04026704400485(17)270228(10)KMZ22C0606
19212T	KMZ22C0611	(01)04026704400485(17)270228(10)KMZ22C0611
19212T	KMZ22C0616	(01)04026704400485(17)270228(10)KMZ22C0616
19212T	KMZ22C0624	(01)04026704400485(17)270228(10)KMZ22C0624
19261	KMH20L0234	(01)04026704388486(17)251028(10)KMH20L0234
19261	KMH20L0283	(01)04026704388486(17)251028(10)KMH20L0283
19261	KMH20L0288	(01)04026704388486(17)251028(10)KMH20L0288
19261	KMH20L0333	(01)04026704388486(17)251028(10)KMH20L0333
19261	KMH21A0035	(01)04026704388486(17)251228(10)KMH21A0035
19261	KMH21B0093	(01)04026704388486(17)260128(10)KMH21B0093
19261	KMH21B0193	(01)04026704388486(17)260128(10)KMH21B0193
19261	KMH21J0246	(01)04026704388486(17)260828(10)KMH21J0246
19261	KMH21L0294	(01)04026704388486(17)261028(10)KMH21L0294
19261	KMH22A0230	(01)04026704388486(17)261228(10)KMH22A0230
19261T	20FT51	(01)04026704400492(17)250528(10)20FT51
19261T	KMH20L0409	(01)04026704400492(17)251028(10)KMH20L0409
19261T	KMH21A0086	(01)04026704400492(17)251228(10)KMH21A0086
19261T	KMH21A0190	(01)04026704400492(17)251228(10)KMH21A0190
19261T	KMH21B0163	(01)04026704400492(17)260128(10)KMH21B0163
19261T	KMH21B0343	(01)04026704400492(17)260128(10)KMH21B0343
19261T	KMH21B0355	(01)04026704400492(17)260128(10)KMH21B0355
19261T	KMH21C0177	(01)04026704400492(17)260228(10)KMH21C0177
19261T	KMH21C0387	(01)04026704400492(17)260228(10)KMH21C0387
19261T	KMH21F0156	(01)04026704400492(17)260528(10)KMH21F0156
19261T	KMH21F0337	(01)04026704400492(17)260528(10)KMH21F0337
19261T	KMH21G0426	(01)04026704400492(17)260628(10)KMH21G0426
19261T	KMH21J0020	(01)04026704400492(17)260828(10)KMH21J0020
19261T	KMH21J0245	(01)04026704400492(17)260828(10)KMH21J0245
19261T	KMH21K0142	(01)04026704400492(17)260928(10)KMH21K0142
19261T	KMH21M0182	(01)04026704400492(17)261128(10)KMH21M0182
19261T	KMH22B0151	(01)04026704400492(17)270128(10)KMH22B0151
19261T	KMH22B0248	(01)04026704400492(17)270128(10)KMH22B0248
19262T	202025	(01)04026704400508(17)250528(10)202025
19262T	202026	(01)04026704400508(17)250528(10)202026
19262T	202030	(01)04026704400508(17)250628(10)202030
19262T	KMZ20J0163	(01)04026704400508(17)250828(10)KMZ20J0163
19262T	KMZ20J0171	(01)04026704400508(17)250828(10)KMZ20J0171
19262T	KMZ20J0582	(01)04026704400508(17)250828(10)KMZ20J0582
19262T	KMZ20J0586	(01)04026704400508(17)250828(10)KMZ20J0586
19262T	KMZ20J0587	(01)04026704400508(17)250828(10)KMZ20J0587
19262T	KMZ20K0681	(01)04026704400508(17)250928(10)KMZ20K0681
19262T	KMZ20L0968	(01)04026704400508(17)251028(10)KMZ20L0968
19262T	KMZ20L0977	(01)04026704400508(17)251028(10)KMZ20L0977
19262T	KMZ20L1275	(01)04026704400508(17)251028(10)KMZ20L1275
19262T	KMZ20M0544	(01)04026704400508(17)251128(10)KMZ20M0544
19262T	KMZ20M0661	(01)04026704400508(17)251128(10)KMZ20M0661
19262T	KMZ21A0175	(01)04026704400508(17)251228(10)KMZ21A0175



**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19262T	KMZ21A0186	(01)04026704400508(17)251228(10)KMZ21A0186
19262T	KMZ21A0187	(01)04026704400508(17)251228(10)KMZ21A0187
19262T	KMZ21A0208	(01)04026704400508(17)251228(10)KMZ21A0208
19262T	KMZ21B0483	(01)04026704400508(17)260128(10)KMZ21B0483
19262T	KMZ21B0485	(01)04026704400508(17)260128(10)KMZ21B0485
19262T	KMZ21B0495	(01)04026704400508(17)260128(10)KMZ21B0495
19262T	KMZ21B0811	(01)04026704400508(17)260128(10)KMZ21B0811
19262T	KMZ21B0918	(01)04026704400508(17)260128(10)KMZ21B0918
19262T	KMZ21B1172	(01)04026704400508(17)260128(10)KMZ21B1172
19262T	KMZ21C0884	(01)04026704400508(17)260228(10)KMZ21C0884
19262T	KMZ21C0910	(01)04026704400508(17)260228(10)KMZ21C0910
19262T	KMZ21C0929	(01)04026704400508(17)260228(10)KMZ21C0929
19262T	KMZ21C0931	(01)04026704400508(17)260228(10)KMZ21C0931
19262T	KMZ21D0419	(01)04026704400508(17)260328(10)KMZ21D0419
19262T	KMZ21F0223	(01)04026704400508(17)260528(10)KMZ21F0223
19262T	KMZ21F0260	(01)04026704400508(17)260528(10)KMZ21F0260
19262T	KMZ21F0528	(01)04026704400508(17)260528(10)KMZ21F0528
19262T	KMZ21G0728	(01)04026704400508(17)260628(10)KMZ21G0728
19262T	KMZ21G0735	(01)04026704400508(17)260628(10)KMZ21G0735
19262T	KMZ21G0739	(01)04026704400508(17)260628(10)KMZ21G0739
19262T	KMZ21H0055	(01)04026704400508(17)260728(10)KMZ21H0055
19262T	KMZ21J0006	(01)04026704400508(17)260828(10)KMZ21J0006
19262T	KMZ21J0007	(01)04026704400508(17)260828(10)KMZ21J0007
19262T	KMZ21J0008	(01)04026704400508(17)260828(10)KMZ21J0008
19262T	KMZ21K0510	(01)04026704400508(17)260928(10)KMZ21K0510
19262T	KMZ21K0524	(01)04026704400508(17)260928(10)KMZ21K0524
19262T	KMZ21K0539	(01)04026704400508(17)260928(10)KMZ21K0539
19262T	KMZ21K0572	(01)04026704400508(17)260928(10)KMZ21K0572
19262T	KMZ21K0578	(01)04026704400508(17)260928(10)KMZ21K0578
19262T	KMZ21K0592	(01)04026704400508(17)260928(10)KMZ21K0592
19262T	KMZ21K0952	(01)04026704400508(17)260928(10)KMZ21K0952
19262T	KMZ21L0218	(01)04026704400508(17)261028(10)KMZ21L0218
19262T	KMZ21M0674	(01)04026704400508(17)261128(10)KMZ21M0674
19262T	KMZ22B0125	(01)04026704400508(17)270128(10)KMZ22B0125
19262T	KMZ22B0345	(01)04026704400508(17)270128(10)KMZ22B0345
19262T	KMZ22B0809	(01)04026704400508(17)270128(10)KMZ22B0809
19262T	KMZ22B0828	(01)04026704400508(17)270128(10)KMZ22B0828
19262T	KMZ22B0838	(01)04026704400508(17)270128(10)KMZ22B0838
19262T	KMZ22B0869	(01)04026704400508(17)270128(10)KMZ22B0869
19262T	KMZ22C0659	(01)04026704400508(17)270228(10)KMZ22C0659
19272	KMZ21E0613	(01)04026704388509(17)260428(10)KMZ21E0613
19272	KMZ21F0939	(01)04026704388509(17)260528(10)KMZ21F0939
19272T	202025	(01)04026704400515(17)250528(10)202025
19272T	202026	(01)04026704400515(17)250528(10)202026
19272T	202027	(01)04026704400515(17)250628(10)202027
19272T	202028	(01)04026704400515(17)250628(10)202028
19272T	202029	(01)04026704400515(17)250628(10)202029
19272T	202030	(01)04026704400515(17)250628(10)202030
19272T	202031	(01)04026704400515(17)250628(10)202031
19272T	KMZ20G0005	(01)04026704400515(17)250628(10)KMZ20G0005
19272T	KMZ20G0026	(01)04026704400515(17)250628(10)KMZ20G0026
19272T	KMZ20G0105	(01)04026704400515(17)250628(10)KMZ20G0105

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19272T	KMZ20G0112	(01)04026704400515(17)250628(10)KMZ20G0112
19272T	KMZ20H0081	(01)04026704400515(17)250728(10)KMZ20H0081
19272T	KMZ20H0096	(01)04026704400515(17)250728(10)KMZ20H0096
19272T	KMZ20H0158	(01)04026704400515(17)250728(10)KMZ20H0158
19272T	KMZ20H0379	(01)04026704400515(17)250728(10)KMZ20H0379
19272T	KMZ20H0530	(01)04026704400515(17)250728(10)KMZ20H0530
19272T	KMZ20H0540	(01)04026704400515(17)250728(10)KMZ20H0540
19272T	KMZ20H0564	(01)04026704400515(17)250728(10)KMZ20H0564
19272T	KMZ20H0812	(01)04026704400515(17)250728(10)KMZ20H0812
19272T	KMZ20H0841	(01)04026704400515(17)250728(10)KMZ20H0841
19272T	KMZ20H1021	(01)04026704400515(17)250728(10)KMZ20H1021
19272T	KMZ20H1024	(01)04026704400515(17)250728(10)KMZ20H1024
19272T	KMZ20H1031	(01)04026704400515(17)250728(10)KMZ20H1031
19272T	KMZ20J0092	(01)04026704400515(17)250828(10)KMZ20J0092
19272T	KMZ20J0102	(01)04026704400515(17)250828(10)KMZ20J0102
19272T	KMZ20J0227	(01)04026704400515(17)250828(10)KMZ20J0227
19272T	KMZ20J0441	(01)04026704400515(17)250828(10)KMZ20J0441
19272T	KMZ20J0535	(01)04026704400515(17)250828(10)KMZ20J0535
19272T	KMZ20J0551	(01)04026704400515(17)250828(10)KMZ20J0551
19272T	KMZ20J0566	(01)04026704400515(17)250828(10)KMZ20J0566
19272T	KMZ20J0640	(01)04026704400515(17)250828(10)KMZ20J0640
19272T	KMZ20J0654	(01)04026704400515(17)250828(10)KMZ20J0654
19272T	KMZ20J0663	(01)04026704400515(17)250828(10)KMZ20J0663
19272T	KMZ20J0893	(01)04026704400515(17)250828(10)KMZ20J0893
19272T	KMZ20J0958	(01)04026704400515(17)250828(10)KMZ20J0958
19272T	KMZ20K0034	(01)04026704400515(17)250928(10)KMZ20K0034
19272T	KMZ20K0196	(01)04026704400515(17)250928(10)KMZ20K0196
19272T	KMZ20K0236	(01)04026704400515(17)250928(10)KMZ20K0236
19272T	KMZ20K0430	(01)04026704400515(17)250928(10)KMZ20K0430
19272T	KMZ20K0659	(01)04026704400515(17)250928(10)KMZ20K0659
19272T	KMZ20K0662	(01)04026704400515(17)250928(10)KMZ20K0662
19272T	KMZ20K0838	(01)04026704400515(17)250928(10)KMZ20K0838
19272T	KMZ20K0946	(01)04026704400515(17)250928(10)KMZ20K0946
19272T	KMZ20K1260	(01)04026704400515(17)250928(10)KMZ20K1260
19272T	KMZ20K1289	(01)04026704400515(17)250928(10)KMZ20K1289
19272T	KMZ20K1294	(01)04026704400515(17)250928(10)KMZ20K1294
19272T	KMZ20K1310	(01)04026704400515(17)250928(10)KMZ20K1310
19272T	KMZ20K1339	(01)04026704400515(17)250928(10)KMZ20K1339
19272T	KMZ20K1340	(01)04026704400515(17)250928(10)KMZ20K1340
19272T	KMZ20L0003	(01)04026704400515(17)251028(10)KMZ20L0003
19272T	KMZ20L0030	(01)04026704400515(17)251028(10)KMZ20L0030
19272T	KMZ20L0031	(01)04026704400515(17)251028(10)KMZ20L0031
19272T	KMZ20L0117	(01)04026704400515(17)251028(10)KMZ20L0117
19272T	KMZ20L0214	(01)04026704400515(17)251028(10)KMZ20L0214
19272T	KMZ20L0351	(01)04026704400515(17)251028(10)KMZ20L0351
19272T	KMZ20L0359	(01)04026704400515(17)251028(10)KMZ20L0359
19272T	KMZ20L0499	(01)04026704400515(17)251028(10)KMZ20L0499
19272T	KMZ20L0537	(01)04026704400515(17)251028(10)KMZ20L0537
19272T	KMZ20L0706	(01)04026704400515(17)251028(10)KMZ20L0706
19272T	KMZ20L0724	(01)04026704400515(17)251028(10)KMZ20L0724
19272T	KMZ20L0885	(01)04026704400515(17)251028(10)KMZ20L0885
19272T	KMZ20L0886	(01)04026704400515(17)251028(10)KMZ20L0886

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19272T	KMZ20L0933	(01)04026704400515(17)251028(10)KMZ20L0933
19272T	KMZ20L0945	(01)04026704400515(17)251028(10)KMZ20L0945
19272T	KMZ20L1056	(01)04026704400515(17)251028(10)KMZ20L1056
19272T	KMZ20L1091	(01)04026704400515(17)251028(10)KMZ20L1091
19272T	KMZ20L1103	(01)04026704400515(17)251028(10)KMZ20L1103
19272T	KMZ20L1104	(01)04026704400515(17)251028(10)KMZ20L1104
19272T	KMZ20M0387	(01)04026704400515(17)251128(10)KMZ20M0387
19272T	KMZ20M0777	(01)04026704400515(17)251128(10)KMZ20M0777
19272T	KMZ20M0778	(01)04026704400515(17)251128(10)KMZ20M0778
19272T	KMZ20M0787	(01)04026704400515(17)251128(10)KMZ20M0787
19272T	KMZ20M0822	(01)04026704400515(17)251128(10)KMZ20M0822
19272T	KMZ20M0842	(01)04026704400515(17)251128(10)KMZ20M0842
19272T	KMZ20M0938	(01)04026704400515(17)251128(10)KMZ20M0938
19272T	KMZ20M0967	(01)04026704400515(17)251128(10)KMZ20M0967
19272T	KMZ20M0990	(01)04026704400515(17)251128(10)KMZ20M0990
19272T	KMZ20M1249	(01)04026704400515(17)251128(10)KMZ20M1249
19272T	KMZ20M1250	(01)04026704400515(17)251128(10)KMZ20M1250
19272T	KMZ20M1257	(01)04026704400515(17)251128(10)KMZ20M1257
19272T	KMZ20M1267	(01)04026704400515(17)251128(10)KMZ20M1267
19272T	KMZ21A0012	(01)04026704400515(17)251228(10)KMZ21A0012
19272T	KMZ21A0024	(01)04026704400515(17)251228(10)KMZ21A0024
19272T	KMZ21A0117	(01)04026704400515(17)251228(10)KMZ21A0117
19272T	KMZ21A0126	(01)04026704400515(17)251228(10)KMZ21A0126
19272T	KMZ21A0155	(01)04026704400515(17)251228(10)KMZ21A0155
19272T	KMZ21A0162	(01)04026704400515(17)251228(10)KMZ21A0162
19272T	KMZ21A0164	(01)04026704400515(17)251228(10)KMZ21A0164
19272T	KMZ21A0203	(01)04026704400515(17)251228(10)KMZ21A0203
19272T	KMZ21A0218	(01)04026704400515(17)251228(10)KMZ21A0218
19272T	KMZ21A0226	(01)04026704400515(17)251228(10)KMZ21A0226
19272T	KMZ21A0368	(01)04026704400515(17)251228(10)KMZ21A0368
19272T	KMZ21A0618	(01)04026704400515(17)251228(10)KMZ21A0618
19272T	KMZ21A0635	(01)04026704400515(17)251228(10)KMZ21A0635
19272T	KMZ21A0693	(01)04026704400515(17)251228(10)KMZ21A0693
19272T	KMZ21A0706	(01)04026704400515(17)251228(10)KMZ21A0706
19272T	KMZ21A0714	(01)04026704400515(17)251228(10)KMZ21A0714
19272T	KMZ21A0725	(01)04026704400515(17)251228(10)KMZ21A0725
19272T	KMZ21A0734	(01)04026704400515(17)251228(10)KMZ21A0734
19272T	KMZ21A0950	(01)04026704400515(17)251228(10)KMZ21A0950
19272T	KMZ21B0101	(01)04026704400515(17)260128(10)KMZ21B0101
19272T	KMZ21B0163	(01)04026704400515(17)260128(10)KMZ21B0163
19272T	KMZ21B0181	(01)04026704400515(17)260128(10)KMZ21B0181
19272T	KMZ21B0183	(01)04026704400515(17)260128(10)KMZ21B0183
19272T	KMZ21B0185	(01)04026704400515(17)260128(10)KMZ21B0185
19272T	KMZ21B0189	(01)04026704400515(17)260128(10)KMZ21B0189
19272T	KMZ21B0227	(01)04026704400515(17)260128(10)KMZ21B0227
19272T	KMZ21B0234	(01)04026704400515(17)260128(10)KMZ21B0234
19272T	KMZ21B0305	(01)04026704400515(17)260128(10)KMZ21B0305
19272T	KMZ21B0320	(01)04026704400515(17)260128(10)KMZ21B0320
19272T	KMZ21B0417	(01)04026704400515(17)260128(10)KMZ21B0417
19272T	KMZ21B0421	(01)04026704400515(17)260128(10)KMZ21B0421
19272T	KMZ21B0435	(01)04026704400515(17)260128(10)KMZ21B0435
19272T	KMZ21B0465	(01)04026704400515(17)260128(10)KMZ21B0465

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19272T	KMZ21B0471	(01)04026704400515(17)260128(10)KMZ21B0471
19272T	KMZ21B0479	(01)04026704400515(17)260128(10)KMZ21B0479
19272T	KMZ21B0533	(01)04026704400515(17)260128(10)KMZ21B0533
19272T	KMZ21B0689	(01)04026704400515(17)260128(10)KMZ21B0689
19272T	KMZ21B1171	(01)04026704400515(17)260128(10)KMZ21B1171
19272T	KMZ21C0270	(01)04026704400515(17)260228(10)KMZ21C0270
19272T	KMZ21C0388	(01)04026704400515(17)260228(10)KMZ21C0388
19272T	KMZ21C0398	(01)04026704400515(17)260228(10)KMZ21C0398
19272T	KMZ21C0883	(01)04026704400515(17)260228(10)KMZ21C0883
19272T	KMZ21C1159	(01)04026704400515(17)260228(10)KMZ21C1159
19272T	KMZ21C1160	(01)04026704400515(17)260228(10)KMZ21C1160
19272T	KMZ21C1358	(01)04026704400515(17)260228(10)KMZ21C1358
19272T	KMZ21D0420	(01)04026704400515(17)260328(10)KMZ21D0420
19272T	KMZ21F0227	(01)04026704400515(17)260528(10)KMZ21F0227
19272T	KMZ21F0265	(01)04026704400515(17)260528(10)KMZ21F0265
19272T	KMZ21F0296	(01)04026704400515(17)260528(10)KMZ21F0296
19272T	KMZ21F0335	(01)04026704400515(17)260528(10)KMZ21F0335
19272T	KMZ21F0375	(01)04026704400515(17)260528(10)KMZ21F0375
19272T	KMZ21F0430	(01)04026704400515(17)260528(10)KMZ21F0430
19272T	KMZ21F0467	(01)04026704400515(17)260528(10)KMZ21F0467
19272T	KMZ21F0468	(01)04026704400515(17)260528(10)KMZ21F0468
19272T	KMZ21F0478	(01)04026704400515(17)260528(10)KMZ21F0478
19272T	KMZ21F0501	(01)04026704400515(17)260528(10)KMZ21F0501
19272T	KMZ21F0517	(01)04026704400515(17)260528(10)KMZ21F0517
19272T	KMZ21F0534	(01)04026704400515(17)260528(10)KMZ21F0534
19272T	KMZ21F0550	(01)04026704400515(17)260528(10)KMZ21F0550
19272T	KMZ21F0551	(01)04026704400515(17)260528(10)KMZ21F0551
19272T	KMZ21F0552	(01)04026704400515(17)260528(10)KMZ21F0552
19272T	KMZ21F0557	(01)04026704400515(17)260528(10)KMZ21F0557
19272T	KMZ21F0567	(01)04026704400515(17)260528(10)KMZ21F0567
19272T	KMZ21F0568	(01)04026704400515(17)260528(10)KMZ21F0568
19272T	KMZ21F0588	(01)04026704400515(17)260528(10)KMZ21F0588
19272T	KMZ21F0594	(01)04026704400515(17)260528(10)KMZ21F0594
19272T	KMZ21F0742	(01)04026704400515(17)260528(10)KMZ21F0742
19272T	KMZ21F0777	(01)04026704400515(17)260528(10)KMZ21F0777
19272T	KMZ21G0399	(01)04026704400515(17)260628(10)KMZ21G0399
19272T	KMZ21G0403	(01)04026704400515(17)260628(10)KMZ21G0403
19272T	KMZ21G0404	(01)04026704400515(17)260628(10)KMZ21G0404
19272T	KMZ21G0751	(01)04026704400515(17)260628(10)KMZ21G0751
19272T	KMZ21G0765	(01)04026704400515(17)260628(10)KMZ21G0765
19272T	KMZ21G0770	(01)04026704400515(17)260628(10)KMZ21G0770
19272T	KMZ21G0778	(01)04026704400515(17)260628(10)KMZ21G0778
19272T	KMZ21G0786	(01)04026704400515(17)260628(10)KMZ21G0786
19272T	KMZ21G0800	(01)04026704400515(17)260628(10)KMZ21G0800
19272T	KMZ21G0806	(01)04026704400515(17)260628(10)KMZ21G0806
19272T	KMZ21G0838	(01)04026704400515(17)260628(10)KMZ21G0838
19272T	KMZ21G0911	(01)04026704400515(17)260628(10)KMZ21G0911
19272T	KMZ21G0912	(01)04026704400515(17)260628(10)KMZ21G0912
19272T	KMZ21G0917	(01)04026704400515(17)260628(10)KMZ21G0917
19272T	KMZ21G0954	(01)04026704400515(17)260628(10)KMZ21G0954
19272T	KMZ21G0955	(01)04026704400515(17)260628(10)KMZ21G0955
19272T	KMZ21G0956	(01)04026704400515(17)260628(10)KMZ21G0956

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19272T	KMZ21G0964	(01)04026704400515(17)260628(10)KMZ21G0964
19272T	KMZ21G0971	(01)04026704400515(17)260628(10)KMZ21G0971
19272T	KMZ21G0982	(01)04026704400515(17)260628(10)KMZ21G0982
19272T	KMZ21G0986	(01)04026704400515(17)260628(10)KMZ21G0986
19272T	KMZ21G0992	(01)04026704400515(17)260628(10)KMZ21G0992
19272T	KMZ21G0993	(01)04026704400515(17)260628(10)KMZ21G0993
19272T	KMZ21G0995	(01)04026704400515(17)260628(10)KMZ21G0995
19272T	KMZ21G1001	(01)04026704400515(17)260628(10)KMZ21G1001
19272T	KMZ21G1008	(01)04026704400515(17)260628(10)KMZ21G1008
19272T	KMZ21G1017	(01)04026704400515(17)260628(10)KMZ21G1017
19272T	KMZ21G1037	(01)04026704400515(17)260628(10)KMZ21G1037
19272T	KMZ21G1038	(01)04026704400515(17)260628(10)KMZ21G1038
19272T	KMZ21H0115	(01)04026704400515(17)260728(10)KMZ21H0115
19272T	KMZ21H0257	(01)04026704400515(17)260728(10)KMZ21H0257
19272T	KMZ21H0299	(01)04026704400515(17)260728(10)KMZ21H0299
19272T	KMZ21H0405	(01)04026704400515(17)260728(10)KMZ21H0405
19272T	KMZ21H0429	(01)04026704400515(17)260728(10)KMZ21H0429
19272T	KMZ21H0431	(01)04026704400515(17)260728(10)KMZ21H0431
19272T	KMZ21H0433	(01)04026704400515(17)260728(10)KMZ21H0433
19272T	KMZ21H0435	(01)04026704400515(17)260728(10)KMZ21H0435
19272T	KMZ21H0504	(01)04026704400515(17)260728(10)KMZ21H0504
19272T	KMZ21H0508	(01)04026704400515(17)260728(10)KMZ21H0508
19272T	KMZ21H0526	(01)04026704400515(17)260728(10)KMZ21H0526
19272T	KMZ21H0560	(01)04026704400515(17)260728(10)KMZ21H0560
19272T	KMZ21H0561	(01)04026704400515(17)260728(10)KMZ21H0561
19272T	KMZ21J0050	(01)04026704400515(17)260828(10)KMZ21J0050
19272T	KMZ21J0077	(01)04026704400515(17)260828(10)KMZ21J0077
19272T	KMZ21J0088	(01)04026704400515(17)260828(10)KMZ21J0088
19272T	KMZ21J0089	(01)04026704400515(17)260828(10)KMZ21J0089
19272T	KMZ21J0446	(01)04026704400515(17)260828(10)KMZ21J0446
19272T	KMZ21J0447	(01)04026704400515(17)260828(10)KMZ21J0447
19272T	KMZ21J0448	(01)04026704400515(17)260828(10)KMZ21J0448
19272T	KMZ21J0484	(01)04026704400515(17)260828(10)KMZ21J0484
19272T	KMZ21J0485	(01)04026704400515(17)260828(10)KMZ21J0485
19272T	KMZ21J0534	(01)04026704400515(17)260828(10)KMZ21J0534
19272T	KMZ21J0546	(01)04026704400515(17)260828(10)KMZ21J0546
19272T	KMZ21J0554	(01)04026704400515(17)260828(10)KMZ21J0554
19272T	KMZ21J0555	(01)04026704400515(17)260828(10)KMZ21J0555
19272T	KMZ21J0570	(01)04026704400515(17)260828(10)KMZ21J0570
19272T	KMZ21J0571	(01)04026704400515(17)260828(10)KMZ21J0571
19272T	KMZ21J0579	(01)04026704400515(17)260828(10)KMZ21J0579
19272T	KMZ21J0580	(01)04026704400515(17)260828(10)KMZ21J0580
19272T	KMZ21K0045	(01)04026704400515(17)260928(10)KMZ21K0045
19272T	KMZ21K0114	(01)04026704400515(17)260928(10)KMZ21K0114
19272T	KMZ21K0163	(01)04026704400515(17)260928(10)KMZ21K0163
19272T	KMZ21K0164	(01)04026704400515(17)260928(10)KMZ21K0164
19272T	KMZ21K0266	(01)04026704400515(17)260928(10)KMZ21K0266
19272T	KMZ21K0274	(01)04026704400515(17)260928(10)KMZ21K0274
19272T	KMZ21K0295	(01)04026704400515(17)260928(10)KMZ21K0295
19272T	KMZ21K0299	(01)04026704400515(17)260928(10)KMZ21K0299
19272T	KMZ21K0308	(01)04026704400515(17)260928(10)KMZ21K0308
19272T	KMZ21K0314	(01)04026704400515(17)260928(10)KMZ21K0314

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19272T	KMZ21K0318	(01)04026704400515(17)260928(10)KMZ21K0318
19272T	KMZ21K0323	(01)04026704400515(17)260928(10)KMZ21K0323
19272T	KMZ21K0452	(01)04026704400515(17)260928(10)KMZ21K0452
19272T	KMZ21K0453	(01)04026704400515(17)260928(10)KMZ21K0453
19272T	KMZ21K0509	(01)04026704400515(17)260928(10)KMZ21K0509
19272T	KMZ21K0926	(01)04026704400515(17)260928(10)KMZ21K0926
19272T	KMZ21K0953	(01)04026704400515(17)260928(10)KMZ21K0953
19272T	KMZ21K0954	(01)04026704400515(17)260928(10)KMZ21K0954
19272T	KMZ21L0047	(01)04026704400515(17)261028(10)KMZ21L0047
19272T	KMZ21L0070	(01)04026704400515(17)261028(10)KMZ21L0070
19272T	KMZ21L0074	(01)04026704400515(17)261028(10)KMZ21L0074
19272T	KMZ21L0077	(01)04026704400515(17)261028(10)KMZ21L0077
19272T	KMZ21L0079	(01)04026704400515(17)261028(10)KMZ21L0079
19272T	KMZ21L0219	(01)04026704400515(17)261028(10)KMZ21L0219
19272T	KMZ21L0220	(01)04026704400515(17)261028(10)KMZ21L0220
19272T	KMZ21L0238	(01)04026704400515(17)261028(10)KMZ21L0238
19272T	KMZ21L0239	(01)04026704400515(17)261028(10)KMZ21L0239
19272T	KMZ21L0244	(01)04026704400515(17)261028(10)KMZ21L0244
19272T	KMZ21L0295	(01)04026704400515(17)261028(10)KMZ21L0295
19272T	KMZ21L0296	(01)04026704400515(17)261028(10)KMZ21L0296
19272T	KMZ21L0307	(01)04026704400515(17)261028(10)KMZ21L0307
19272T	KMZ21L0308	(01)04026704400515(17)261028(10)KMZ21L0308
19272T	KMZ21L0316	(01)04026704400515(17)261028(10)KMZ21L0316
19272T	KMZ21L0333	(01)04026704400515(17)261028(10)KMZ21L0333
19272T	KMZ21L0334	(01)04026704400515(17)261028(10)KMZ21L0334
19272T	KMZ21L0337	(01)04026704400515(17)261028(10)KMZ21L0337
19272T	KMZ21L0344	(01)04026704400515(17)261028(10)KMZ21L0344
19272T	KMZ21L0367	(01)04026704400515(17)261028(10)KMZ21L0367
19272T	KMZ21L0368	(01)04026704400515(17)261028(10)KMZ21L0368
19272T	KMZ21L0382	(01)04026704400515(17)261028(10)KMZ21L0382
19272T	KMZ21L0389	(01)04026704400515(17)261028(10)KMZ21L0389
19272T	KMZ21L0541	(01)04026704400515(17)261028(10)KMZ21L0541
19272T	KMZ21L0798	(01)04026704400515(17)261028(10)KMZ21L0798
19272T	KMZ21L0878	(01)04026704400515(17)261028(10)KMZ21L0878
19272T	KMZ21L0880	(01)04026704400515(17)261028(10)KMZ21L0880
19272T	KMZ21L0901	(01)04026704400515(17)261028(10)KMZ21L0901
19272T	KMZ21L0924	(01)04026704400515(17)261028(10)KMZ21L0924
19272T	KMZ21M0050	(01)04026704400515(17)261128(10)KMZ21M0050
19272T	KMZ21M0690	(01)04026704400515(17)261128(10)KMZ21M0690
19272T	KMZ21M0692	(01)04026704400515(17)261128(10)KMZ21M0692
19272T	KMZ21M0719	(01)04026704400515(17)261128(10)KMZ21M0719
19272T	KMZ21M0756	(01)04026704400515(17)261128(10)KMZ21M0756
19272T	KMZ21M0781	(01)04026704400515(17)261128(10)KMZ21M0781
19272T	KMZ21M0786	(01)04026704400515(17)261128(10)KMZ21M0786
19272T	KMZ21M0790	(01)04026704400515(17)261128(10)KMZ21M0790
19272T	KMZ21M0805	(01)04026704400515(17)261128(10)KMZ21M0805
19272T	KMZ22A0635	(01)04026704400515(17)261228(10)KMZ22A0635
19272T	KMZ22A0636	(01)04026704400515(17)261228(10)KMZ22A0636
19272T	KMZ22A0650	(01)04026704400515(17)261228(10)KMZ22A0650
19272T	KMZ22A0871	(01)04026704400515(17)261228(10)KMZ22A0871
19272T	KMZ22B0130	(01)04026704400515(17)270128(10)KMZ22B0130
19272T	KMZ22B0131	(01)04026704400515(17)270128(10)KMZ22B0131

**Appendix 2: Product in scope of EIF-000513 and EIF-000513-01**

Product Code	Batch	Unique Device Identifier
19272T	KMZ22B0192	(01)04026704400515(17)270128(10)KMZ22B0192
19272T	KMZ22B0509	(01)04026704400515(17)270128(10)KMZ22B0509
19272T	KMZ22B0524	(01)04026704400515(17)270128(10)KMZ22B0524
19272T	KMZ22B0525	(01)04026704400515(17)270128(10)KMZ22B0525
19272T	KMZ22B0556	(01)04026704400515(17)270128(10)KMZ22B0556
19272T	KMZ22B0557	(01)04026704400515(17)270128(10)KMZ22B0557
19272T	KMZ22B0784	(01)04026704400515(17)270128(10)KMZ22B0784
19272T	KMZ22B0785	(01)04026704400515(17)270128(10)KMZ22B0785
19272T	KMZ22C0395	(01)04026704400515(17)270228(10)KMZ22C0395
19272T	KMZ22C0620	(01)04026704400515(17)270228(10)KMZ22C0620
19272T	KMZ22C0802	(01)04026704400515(17)270228(10)KMZ22C0802
19272T	KMZ22C0803	(01)04026704400515(17)270228(10)KMZ22C0803

- end of document -