



Dear Customer

17 December 2021

This letter contains important information about your **Laerdal Suction Unit with Reusable Canister (LSU) and Laerdal Suction Unit Reusable Canister**. Please read this letter and the attached Urgent Field Safety Notice.

Laerdal Medical has become aware of a problem that may cause leakage in the canister or canister lid. This will lead to lack of suctioning or reduced suctioning.

Please note that this Urgent Field Safety Notice is only applies to Laerdal Suction Unit Reusable Canisters **manufactured** in the period 17 April 2021 and 2 July 2021.

The intention of this letter and the attachments is to provide you with a description of the problem and information about:

- Which units are affected;
- Actions to be taken by you as customer
- Actions to be taken by Laerdal Medical.

If you should have any questions, please contact your local Laerdal Medical Company/distributor customer service organization.

This Urgent Field Safety Notice will be reported to the appropriate authorities.

Laerdal Medical apologises for the inconvenience caused by this problem. Releasing products and services of the highest quality is at the root of our mission - helping save lives - and your continued trust in us will always be our top priority.

Kindest regards,

