

Urgent – Product Recall 1st Notification – Urgent

Date: May 6, 2011

Recalled Product Identification:

Smith & Nephew, DYONICS[®] POWER Shaver Blades, Full Radius, Long Curved Convex 4.5 mm, and
Smith & Nephew, DYONICS POWER Shaver Blades, Full Radius, Long Straight 4.5 mm

This letter is to inform you of a voluntary, Class II recall from Smith & Nephew Endoscopy for two lots of our DYONICS POWER Shaver Full Radius Blades. We have recently identified an issue with insufficient weld penetration, which could result in the tip separating from the shaft tube during use. Product details are indicated below:

Product #	Description	Lot #	Individual Units to be Returned	No Product to Return
7205336	DYONICS POWER Shaver Blades, Full Radius, Long Curved Convex 4.5 mm, package of 3	50540107		
7205341	DYONICS POWER Shaver Blades, Full Radius, Long Straight 4.5 mm, package of 3	50561144		

Please complete the last two columns in the table above, indicating the quantities to be returned, and include your phone and fax number. If you do not have product to return, please place an "x" in the column "no product to return".

To return product and obtain credit or order replacement blades, please contact the Smith & Nephew Endoscopy Returns at 800-343-5717 (option 3) for a Return Authorization (RA) or send an email to endo.andreturns@smith-nephew.com.

Completed by Smith & Nephew Sales Rep:

Name (Print): _____ Signature: _____ Date: _____

Contact Phone #: _____ Fax: _____

Completed by Hospital Rep:

Name (Print): _____ Signature: _____ Date: _____

Contact Phone #: _____ Fax: _____

Please return this form to Hoangthi Le by email at [redacted]@smith-nephew.com or by fax at 508-261-3620. If you have questions, you may also contact [redacted] by phone at 508-337-3731.

We apologize for any inconvenience this may cause.

Respectfully,

[redacted]

[redacted]

Director, Quality

[®]Trademark of Smith & Nephew. Registered US Patent and Trademark Office.