



Urgent Field Safety Notice

Dolomite Soprano rollator

Date: 25th February 2009
Internal ref: 2009-02-04:1
Item no 12030-44

Safety Information:

In the interests of the safety of our customers we would like to bring to your attention to known incidents with product Dolomite Soprano where the brake function can be insufficient due to an early wear of the components and the rear wheel can become loose, the rollator is then unstable, after some time of use in this state the wheel can fall off.

The affected products serial numbers are within range of G065-000001 to G065-000800 with manufacturing dates 2008-05 to 2008-08.

During you wait for the replacement use the product with caution.

Dolomite will replace products on the market from above mentioned time period.

Your local sales representative contacts will contact you to set up time schedule for exchange of your product.

Authorities have been informed regarding this case.

If you have any questions relating to this bulletin, or should you require any additional information, please contact your local Invacare sales (address below).

Contacts

ADDRESS AND PHONE NO OF LOCAL SALES

XXXXXX

YYYYYYY

Invacare Rea AB

Box 200, Växjövägen 303, SE 343 75 Diö Sweden Phone: +46 476 535 00 Fax: +46 476 535 99