

Breeze C Scooter Field Safety Notice

October 28, 2013

Kibbutz Afikim, 15148, Israel
+972-4-6754180 / www.afiscooters.com

On October 28, 2013, Afikim Electric Vehicles initiated a recall of 13 Breeze C Scooters. The Breeze C Scooter may become unresponsive during operation due to an electrical contact which may fault due to extensive wear, which potentially could result in the driver not being able to stop the scooter during operation except by turning off the ignition key.

Recalled Breeze C Scooters were sold between 15/08/2012 and 02/04/2013
The following styles have been recalled:

Breeze C, 4 Wheel Scooters, blue and silver color

Serial Numbers:

S/N	INV. DATE
12490071	15/08/2012
13490501	02/04/2013
13490502	02/04/2013
13490503	02/04/2013
13490509	02/04/2013
13490510	02/04/2013
13490511	02/04/2013
13490512	02/04/2013
13490513	02/04/2013
13490520	02/04/2013
13490521	02/04/2013
13490522	02/04/2013
13490524	02/04/2013

Total Number Recalled is 13.

Product(s) can be identified by their color; the recalled products are either metallic blue or metallic silver. The products have two rear wheels and have two wheels in the front. The scooters have an adjustable swivel seat and a front storage basket; some scooters may be fitted with an optional larger rear basket.

Afikim Electric Vehicles voluntarily recalled Breeze C Scooters after becoming aware of a potential for malfunction.

During product testing we became aware of a potential scenario where an internal potentiometer wiper could become disconnected due to extensive wear. If this were to occur the device could continue to drive forward at a slow speed of approximately 3-5 Km/h and the only way to stop the device would be by turning off the main key switch. The number of scooters affected is 13; no injuries have been reported to date.

Afikim Electric Vehicles is notifying its distributor by email and is arranging for correction of all recalled Breeze C Scooters.

All above mentioned 13 scooters were sold by Afikim Electric Vehicles to OxyCare GmbH in Bremen.

Any customer/distributor who is in possession of one of the above mentioned scooters should contact the manufacturer at +972-46754180 or via Email: ortal@afiscooters.com

Authorized Representative in EU for Afikim Electric Vehicles:

5 Beaumont, Gate, Shenley Hill
Radlett, Hertfordshire
WD7 7AR UK

 arazygroup.com