



FSN Ref: 2023-10(02)
Date: 13 Nov 2023

FSCA Ref: 2023-10(02)

Urgent Field Safety Notice
Mölnlycke® Procedure Trays

For Attention of: Theatre Manager

Contact details of local representative (name, e-mail, telephone, address etc.)

Name: Local Customer Care contact will be added for each specific market

Email: XXX.XXX@mölnlycke.com

Telephone: +XXXXXXXXXXXXXXXXXX

Urgent Field Safety Notice (FSN)
Mölnlycke® Procedure Trays
MEDECO Gauze Swabs within Mölnlycke® Procedure Trays

1. Information on Affected Devices	
1.	<p style="text-align: center;">1. Device Type(s)</p> <p>Component from MEDECO: Gauze lap 20x25cm 20Th 16P Loop XRD Whit</p>  <p>Gauze swab 10x10cm 17Th 12P XRD White 5pc</p>  <p>Gauze swab 10x20cm 17Th 12P XRD White 5pc</p>  <p>Included in various Mölnlycke® Procedure Trays. Mölnlycke® Procedure Trays consist of customized configurations of several sterilized components, which are assembled and delivered sterile within one procedure Tray.</p>

FSN Ref: 2023-10(02)

FSCA Ref: 2023-10(02)

Date: 13 Nov 2023

1.	2. Commercial name(s) See Appendix I Product Table
1.	3. Primary clinical purpose of device(s) The Gauze swab XRD products are used to absorb fluids during surgery in OR. The clinical purpose of Mölnlycke® Procedure Trays is to provide a customized sterile co-packing of components for different clinical interventions.
1.	4. Device Model/Catalogue/part number(s) See Appendix I Product Table
1.	5. Affected serial or lot number range See Appendix I Product Table

2 Reason for Field Safety Corrective Action (FSCA)	
2	1. Description of the product problem* Mölnlycke has recently been informed by the legal manufacturer MEDECO that they are performing a recall on components Klinion Kliniray Gauze Compress X-Ray and Klinion Kliniray Abdominal Gauze Compress X-Ray, which Mölnlycke includes in some of the Mölnlycke® Procedure trays. The reason for this recall is that after assessment of these products, it became known that x-ray threads in gauze compresses can disintegrate. It has been confirmed that these products are not meeting the required criteria. Theoretically, small pieces of thread could be left behind if the x-ray thread breaks or frays. No patient harm has been reported. Mölnlycke has decided to follow the legal manufacturer FSN and perform a Field Safety Corrective Action. At the point of use of these Mölnlycke® Procedure trays, the user is required to identify the affected component and discard it.
2	2. Hazard giving rise to the FSCA* Information from MEDECO: Theoretically, small pieces of thread could be left behind if the x-ray thread breaks or frays. This may, for example, lead to inflammation and/or granuloma formation when remaining in the body.



3. Type of Action to mitigate the risk	
3.	1. Action To Be Taken by the User <input checked="" type="checkbox"/> Identify Device <input checked="" type="checkbox"/> Destroy Device We need your help in ensuring that all affected products are located and that below actions are performed. Please follow below instructions: 1. Identify and isolate the unused Mölnlycke® Procedure Trays at your facility, please see Appendix I for affected product information. 2. Attach the Appendix II to all unused Mölnlycke® Procedure trays. 3. At the point of use of the tray , the user is required to remove the affected component from the Mölnlycke® Procedure tray and destroy them .

FSN Ref: 2023-10(02)

FSCA Ref: 2023-10(02)

Date: 13 Nov 2023

	<p>4. Fill out the Customer Reply Form or Distributor Reply Form with the quantity of identified affected products. Please sign and email/fax the Customer Reply Form or Distributor Reply Form per its instructions within 10 business days.</p> <p>5. Even if you no longer have any concerned Mölnlycke® Procedure trays, fill out the Customer Reply Form or Distributor Reply Form and return it back within 10 business days. Mölnlycke needs to be sure all customers are aware of the situation.</p> <p>6. Mölnlycke will contact you as soon as you return the Customer Reply Form or Distributor Reply Form. Mölnlycke will issue a credit for the goods destroyed.</p> <p>7. If you have forwarded any affected products to other healthcare institutions, please send them a copy of this Field Safety Notice. Make sure they act accordingly.</p> <p>8. If you are a distributor, please inform your customers by sending them a copy of this Field Safety Notice. Make sure they act accordingly and return the Distributor Reply Form with information collected from your end-users.</p> <p>We apologize for any inconvenience this will cause you, and rest assured it is our utmost intent to make this process as easy for you as possible. In addition, Mölnlycke appreciates your help in collecting data on product complaints and/or incidents related to the concerned product. Please follow the reporting procedures established by your facility.</p>	
3.	2. Is customer Reply Required?	Yes (Within 10 business days)

4. General Information		
4.	1. FSN Type	New
4.	2. Further advice or information already expected in follow-up FSN?	No
4.	3. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	Mölnlycke Health Care AB
	b. Address	Box 130 80, SE-402 52 Gothenburg, Sweden
	c. Website address	www.molnlycke.com
4.	4. The Competent (Regulatory) Authority of your country has been informed about this communication to customers.	
4.	5. List of attachments/appendices:	Appendix I Product table Appendix II Tag to attach on affected Trays Customer Reply Form Distributor Reply Form
4.	6. Name/Signature	
		

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>

FSN Ref: 2023-10(02)

Date: 13 Nov 2023

Appendix I

FSCA Ref: 2023-10(02)

Product table

To be added for each market

FSN Ref: 2023-10(02)
Date: 13 Nov 2023

FSCA Ref: 2023-10(02)

Appendix II

Tag to be attached to affected Mölnlycke® Procedure Trays(unused)

Action To Be Taken by the User

At the point of use, the user is required to **remove affected component** from the Mölnlycke® Procedure tray and **destroy them**.

Gauze lap 20x25cm 20Th 16P Loop XRD Whit – Mölnlycke Component code: 2310347-00



Gauze swab 10x10cm 17Th 12P XRD White 5pc – Mölnlycke Component code: 2318508-00



Gauze swab 10x20cm 17Th 12P XRD White 5pc – Mölnlycke Component code: 2318534-00



Mölnlycke Ref: 2023-10(02) - Product Table							
Field Safety Notice Ref: 2023-10(02) - Mölnlycke® Procedure Trays				Affected Products numbers and class globally			
Product	LOT	Product Name	Material group description	GMDN code	GMDN term	GMDN definition	MDD Class
97006182-13	21015983	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006182-13	21032739	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006182-13	21047918	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006182-13	21052878	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006182-13	21069530	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-11	21069286	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006184-11	21082469	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006186-08	19451600	Aneurysma pakket	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97006186-08	21077910	Aneurysma pakket	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97006187-21	20534457	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-21	20534461	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-21	21044558	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III

97006187-21	21050844	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-21	21051879	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006287-12	21031627	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	21048619	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	21052036	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	21060961	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006287-12	21234748	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	21257651	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006290-08	21049475	Carotis Pakket	Carotis	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97006387-09	21077984	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006391-10	21025596	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006391-10	21049314	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97007231-08	21069272	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	21069350	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	21255823	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007233-05	21057147	Vaginaal Pakket	Vaginal hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97007234-08	21016117	Perifeer vaat pakket	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97007234-08	21026784	Perifeer vaat pakket	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA

97007465-06	21073307	Tulband Pakket	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007477-09	21017591	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007477-09	21032765	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007477-09	21052879	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007477-09	21058055	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007479-11	19452572	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97007479-11	21015988	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	21058020	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	21026811	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	21049333	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	21069508	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	21186087	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97007482-06	21209475	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	21057112	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	21073592	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-18	21037443	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-18	21038021	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-18	21049236	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97008039-18	21226654	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008652-07	21074754	Neuro Rug Nek Pakket	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008654-05	21058131	Gatlagen Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008656-04	19451393	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008656-04	21017770	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008656-04	21025587	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97008656-04	21039277	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008656-04	21069507	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97011376-10	21061956	VKB Pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97019003-10	21019206	Sectio Tray	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97019003-10	21019207	Sectio Tray	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97019003-11	21026824	Sectio Tray	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97035304-13	21037998	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97035304-13	21038024	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97035304-13	21221682	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97035304-13	21231579	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036257-12	21056373	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97036257-12	21056518	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97038473-07	21062879	ARTHROSCOPIE TRAY	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97046828-05	21062890	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	21063062	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	21213388	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97052281-05	21018231	DIEP Pakket	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	21016039	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97058036-07	21047986	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97058036-07	21052842	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97058036-07	21054857	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97061798-04	21016048	FESS navigatie tray	MIS Nose/Sinus	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97062064-07	21021712	Laparoscopie abdo-perineaal tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97062077-04	21061697	Laparoscopie algemeen tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97064671-05	21019609	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	21019610	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	21019611	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-06	21044159	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-06	21044238	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-06	21242931	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97068089-07	21037137	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97068089-07	21037187	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97068089-07	21038238	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97069959-04	21016041	Bekkenbodem Tray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97071539-07	21015647	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97071539-07	21015731	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA

97071539-07	21032650	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97071539-07	21032727	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97071539-07	21063001	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97071539-07	21063166	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97074259-01	21057869	Ortho Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074553-02	21024252	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97074553-02	21047269	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074606-01	21034788	Ortho Hand/Voet tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074635-04	21026709	Extractie Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074870-01	21069271	Basis groot tray SKZ	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074874-04	21032887	Splitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074875-01	21015244	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97074875-01	21017520	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074875-01	21022695	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074875-01	21045597	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074875-01	21052783	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074881-03	21057220	Hand- ok Tray EMC Sophia	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074882-04	21015625	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III

97074882-04	21015719	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074882-04	21046707	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074882-04	21046743	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074882-04	21053073	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074882-04	21069499	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074987-02	21077117	TUR Tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97074989-02	21032467	Gynaecologie Laparoscopie Tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97074990-04	21015720	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074990-04	21015544	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074990-04	21022704	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074990-04	21032805	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074990-04	21034865	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97075012-04	21014290	Vrije Lap Tray Erasmus MC	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97075012-04	21077062	Vrije Lap Tray Erasmus MC	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97075013-02	21016428	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075016-02	21016411	VKB Tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	III
97075019-02	21015565	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075019-02	21022623	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97075019-02	21024123	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075019-02	21034758	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075019-02	21052898	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075020-04	20117897	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075020-05	20224381	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075020-05	20237166	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075020-05	20216658	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075020-05	20257466	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075020-05	20257502	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075020-05	20469210	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075020-05	21015590	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075020-05	21032810	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075020-05	21054182	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	21015530	EMC Extremitäten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	21047995	EMC Extremitäten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	21052789	EMC Extremitäten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075025-04	21032469	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97075025-04	21064386	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III

97075042-02	21015527	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97075042-02	21032606	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97075042-02	21052814	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076706-01	21014166	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21014219	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21014257	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076706-01	21014261	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non- medicated, single- use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general- purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21014285	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non- medicated, single- use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general- purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21051295	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non- medicated, single- use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general- purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21051845	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non- medicated, single- use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general- purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21051846	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non- medicated, single- use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general- purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21058339	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non- medicated, single- use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general- purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076706-01	21058342	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21077871	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	21077959	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076735-00	21041131	Aangezichtstray	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	21025072	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	21025073	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076737-00	21045860	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non- medicated, single- use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	21045860	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non- medicated, single- use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	21058436	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non- medicated, single- use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	21058586	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non- medicated, single- use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	21077850	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non- medicated, single- use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	21249304	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non- medicated, single- use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076737-00	21258738	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076754-01	21069061	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076754-01	21236606	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	21048277	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	21097680	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	21117621	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076755-02	21179054	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	21214043	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	21256155	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	21344418	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076773-01	21045904	Craniosynostose Open/Endo tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	20534610	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076774-01	21038993	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	21060461	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	21191261	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076776-00	21024178	CVL tray	CVC/impl. chbr plac	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076781-02	21228557	Diagn. Laparoscopie tray Gynaecologie	MIS Extraterine (Myomectomy/Salpingo- oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076815-01	21051847	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97076815-01	21074736	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076816-01	21184243	DIEP tray Onder	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076877-02	21021925	Drain Volwassene en Kind tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-02	21069513	Drain Volwassene en Kind tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-02	21199177	Drain Volwassene en Kind tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-02	21215254	Drain Volwassene en Kind tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076879-01	21019174	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076879-01	21019205	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076882-00	21213387	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	21213394	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076883-01	21034876	Fem-Pop tray	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076885-01	21015963	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076885-01	21021735	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	21024170	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	21031721	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	21047924	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	21054603	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	21197331	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076885-01	21216730	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	21069038	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	21069158	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076888-02	21021957	Incontinentietray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076892-01	21069465	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	21213505	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076894-01	21191043	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076894-01	21214697	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076895-01	19462219	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	21037121	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	21037208	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	21074621	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97076901-03	21025764	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076901-03	21213419	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076902-01	21024336	Pediatric Laparoscopie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076902-01	21047283	Pediatric Laparoscopie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	21016172	Pediatric tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076909-01	21054806	Rad. Vulvectomy tray	Abdom. hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB

97076932-01	21024270	Re-Sternotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076987-01	21010211	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	21032599	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	21045589	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	21052077	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	21191264	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076987-01	21213308	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077001-01	21016171	Stent tray	Endovascular Others	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97077001-01	21192903	Stent tray	Endovascular Others	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97077006-01	21019261	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97077006-01	21026607	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97077006-01	21034776	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB

97077006-01	21208362	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97077007-01	21022960	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97077007-01	21197338	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97077008-00	21038156	Verticaal isolatielaken tray	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077234-02	21021764	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97077234-02	21032777	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97077235-02	21058428	KNO oor tray	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97082276-01	21032658	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97082276-01	21056271	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97082276-01	21058025	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97082386-04	21023572	AVR tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97082386-04	21051896	AVR tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97082386-04	21095109	AVR tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97084957-02	20058639	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97084957-02	20089511	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97084957-02	20277559	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97084957-02	20363993	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97084957-02	21022908	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97084957-02	21191044	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97085666-00	21025492	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97085666-00	21047275	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97087964-02	21044165	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97087964-02	21089326	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97087964-02	21255559	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB

97087965-02	21044524	Laparoscopische uterus tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB
97087965-02	21047183	Laparoscopische uterus tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB
97091291-00	21034710	Mamma Tray chirurgie	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97091612-02	21032722	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97091612-02	21053760	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97091738-01	21055297	Angio tray	Angioplasty (Coronary - PTCA, Carotid or Peripheral) - dilatation with or w/o stent)	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97092416-02	21056613	KNO Bijholte Tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97092568-01	21069129	Hysteroscopie tray	MIS Extraterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97096718-01	21048106	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97096853-01	21058225	Neuro nervus ulnaris	Peripheral nerves	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97100613-03	21103546	TAAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97103564-00	21025565	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97103564-00	21047282	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97103564-00	21054771	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97113929-00	21065547	Lipofilling tray	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	19445906	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076893-02	19455990	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076737-00	19457277	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076706-01	19458742	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	19458750	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97071539-06	19459710	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97006187-18	19459915	Knienpakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-17	19450017	Knienpakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006287-11	19463082	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006287-11	19463083	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	19463277	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006391-10	19463365	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	19463516	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-03	19464249	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97074875-01	19464510	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075042-02	19464531	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97074874-03	19464623	Spitlaken Tray SKZ	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97085665-00	19465451	Universeel klein	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97008654-05	19465474	Gatlaken Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	19466844	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97050152-09	19467543	Knief tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB

97058036-07	19467654	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076877-01	19477872	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074875-01	19479734	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075019-02	19479779	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97074874-03	19479875	Splitlaken Tray SKZ	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97091291-00	19470703	Mamma Tray chirurgie	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97008656-04	19470723	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006391-10	19470745	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	19470783	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97096718-00	19470799	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97058036-07	19470814	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97007230-10	19473544	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA

97036257-10	19484510	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97007230-10	19484615	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97036257-10	19484663	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97075019-02	19485411	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97071539-06	19485420	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97085665-00	19485428	Universeel klein	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075042-02	19485429	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97074882-03	19485450	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97071539-06	19487038	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97006391-10	19487054	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	19487565	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97091291-00	19488238	Mamma Tray chirurgie	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97006391-10	19488317	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	19488795	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97068089-04	19499073	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97007234-08	19489127	Perifeer vaat pakket	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97008654-05	19489304	Gatlaken Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008656-04	19489357	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97075015-02	19499472	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-03	19499503	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97087964-02	19499514	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97064671-03	19499615	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97091612-01	19489779	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97082276-01	19490232	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97050152-09	19490235	Knie tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97011376-10	19490308	VKB Pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97007230-10	19490337	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97006187-18	19490605	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97078860-00	19490653	Gyn Hysteroscopie tray	MIS Hysterectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97007477-09	19490670	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97062064-07	19491376	Laparoscopie abdo-perineaal tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97046828-05	19491388	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076754-01	19491481	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074820-03	19491536	PBC Algemeen Erasmus MC Chir.Dagbehandeling	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075012-04	19491688	Vrije Lap Tray Erasmus MC	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97092126-00	19491692	Hoofd/Hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97074987-02	19491720	TUR Tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007234-08	19492503	Perifeer vaat pakket	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076893-02	19492535	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97008039-17	19492607	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	19492662	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97077008-00	19492941	Verticaal isolatielaken tray	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97035304-12	19492960	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077006-12	19494325	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	19494345	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97102963-01	19494445	Tonsillectomie Tray	Tonsillectomy/Adenoïdectomy	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97019003-10	19494461	Sectio Tray	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97077234-02	19494486	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97077235-02	19494489	KNO oor tray	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97092416-02	19494507	KNO Bijholte Tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97090390-01	19494624	Laparoscopische Urologie Tray	MIS Nephrectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB
97046828-05	19494688	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	19494789	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97082859-00	19494819	Kinder tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006182-11	19494865	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075025-01	19495045	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97075042-02	19495274	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97075018-01	19495290	Algemene Hand Tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	19495334	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97076882-00	19495362	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076987-01	19505364	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	19495368	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075015-02	19495492	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074990-03	19495522	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97076987-01	19506114	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	19506476	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076877-01	19506481	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076893-02	19506827	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076886-01	19507052	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	19507053	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076815-01	19507065	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076816-01	19507092	DIEP tray Onder	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97076886-01	19507167	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-01	19507348	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-01	19507350	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075016-02	19507450	VKB Tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	III
97064671-03	19507487	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97091738-00	19507587	Angio tray	Angioplasty (Coronary - PTCA, Carotid or Peripheral) - dilatation with or w/o stent	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97065522-05	19508280	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97075021-04	19508565	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076815-01	19508611	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97006184-10	19508663	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	19500168	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97079569-00	19511538	Plastisch DIEP tray	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97064671-03	19511855	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076932-01	19512795	Re-Sternotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076877-01	19514905	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	19515104	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	19515123	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077005-01	19516154	Thoracotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076893-02	19526288	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97046828-05	19516437	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076887-01	19516474	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	19516510	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-01	19516528	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077007-01	19516729	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076773-01	19526862	Craniosynostose Open/Endo tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	19527354	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076885-01	19527402	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	19527744	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076706-01	19527783	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	19527785	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076895-01	19527787	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	19527791	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97087965-02	19528054	Laparoscopische uterus tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB
97076987-01	20011038	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20011046	Extremiteten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20022166	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006184-10	20022601	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076754-01	20023363	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	20023537	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077006-01	20025583	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076706-01	20025719	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20026243	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97008039-17	20026756	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20026789	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20026819	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20026835	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077001-01	20026994	Stent tray	Endovascular Others	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076902-01	20027038	Pediatric Laparoscopie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076886-01	20027213	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006187-18	20027273	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076886-01	20027289	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20027354	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97100613-02	20028114	TAAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97068089-04	20038437	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97076816-01	20038463	DIEP tray Onder	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076774-01	20039864	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076900-02	20020032	Oncologie tray KNO	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20032945	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20033027	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20033194	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97076706-01	20033195	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20033317	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076888-01	20034682	Incontinentietray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076706-01	20034848	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-03	20045459	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97077005-01	20035814	Thoracotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076987-01	20035912	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076754-01	20046158	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-03	20046537	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97007482-06	20047927	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20041842	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20041858	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97087964-02	20042074	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076877-01	20043077	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076781-01	20053537	Diagn. Laparoscopie tray Gynaecologie	MIS Extruterine (Myomectomy/Salpingo- oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076987-01	20054233	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20055178	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20056408	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076933-01	20056519	Rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97078039-17	20058258	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20069579	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	20069580	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97065522-05	20050059	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076882-00	20064972	Extremitäten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076885-01	20064973	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076776-00	20064983	CVL tray	CVC/impl. chbr plac	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20065027	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076987-01	20065047	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20065138	Extremiteten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077006-01	20065440	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB

97077007-01	20066355	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076904-02	20066490	Pediatric tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076902-01	20066571	Pediatric Laparoscopy tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076933-02	20066584	Rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20066619	Pediatric tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20076698	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97076706-01	20076700	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20076797	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076893-02	20076820	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076737-00	20076825	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20076885	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006187-18	20077440	Knepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III

97046828-05	20077660	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20077668	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-01	20077686	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076751-00	20070644	Arthroscopie tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97076892-01	20071605	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	20083014	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97076895-01	20083063	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076781-01	20086001	Diagn. Laparoscopie tray Gynaecologie	MIS Extraterine (Myomectomy/Salpingo- oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076882-00	20086959	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20086960	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076932-01	20087299	Re-Sternotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076754-01	20088270	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97077006-01	20088359	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076774-01	20089275	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	20089303	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20089320	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006187-18	20091280	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076893-02	20091305	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076706-01	20091593	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20091603	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20091647	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006187-18	20092710	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076894-01	20092991	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076892-01	20093096	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006287-12	20093524	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	20093561	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036246-09	20094417	Laparotomie tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20094880	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-10	20105477	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20095834	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076756-01	20106992	CAD tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007233-05	20107436	Vaginaal Pakket	Vaginal hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20108546	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074875-01	20108553	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075018-01	20108563	Algemene Hand Tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074882-03	20108578	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III

97074882-03	20108644	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074874-03	20108660	Splitlaken Tray SKZ	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074990-03	20108699	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97085666-00	20108732	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075013-02	20108865	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074989-02	20108887	Gynaecologie Laparoscopie Tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97075021-04	20108922	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075021-04	20109043	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075021-04	20109044	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97092126-00	20109054	Hoofd/Hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075021-04	20109107	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075021-04	20109108	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075021-04	20109109	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075021-04	20109110	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075015-02	20109144	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97087965-02	20109148	Laparoscopische uterus tray	MIS Extruterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB
97092354-02	20109176	Abdominaal Gynaecologie	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97038473-06	20109590	ARTHROSCOPIE TRAY	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA

97068089-06	20109630	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20109633	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97075017-01	20109634	Wervel Tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97076877-01	20109692	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007477-09	20109723	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97006391-10	20109730	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97076892-01	20109735	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97099157-00	20109884	Kaak tray SKZ	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97007465-06	20109927	Tulband Pakket	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076879-01	20109932	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97103564-00	20109962	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97076706-01	20109971	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97062077-04	20100044	Laparoscopie algemeen tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97076737-00	20100051	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20100060	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076706-01	20100140	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20100281	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	20100308	Extremiteiten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97007488-06	20100311	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20100320	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20100321	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	20100385	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-05	20100542	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20100543	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97064671-05	20100544	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076987-01	20101198	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076883-01	20101217	Fem-Pop tray	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076885-01	20101219	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076879-01	20101632	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97062077-04	20101709	Laparoscopie algemeen tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97074259-01	20101719	Ortho Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077003-00	20101771	Subtray Otologie	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077235-02	20101815	KNO oor tray	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074635-02	20101928	Extractie Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20102587	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	20102598	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97008654-05	20102721	Gataken Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007465-06	20112814	Tuband Pakket	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97103564-00	20112824	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97076895-01	20112858	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076751-00	20102869	Arthroscopie tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97076895-01	20112896	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97081333-01	20112939	Gyn kleine verrichtingen tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	20113716	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97006290-08	20114051	Carotis Pakket	Carotis	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97074874-03	20114920	Splittaken Tray SKZ	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074990-03	20115029	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97091612-01	20115047	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97071539-06	20115087	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97071539-06	20115088	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97006187-18	20115327	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-18	20115347	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97075014-01	20115493	Vorbereidings Tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076900-02	20116865	Oncologie tray KNO	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97085491-00	20117038	Arthroscopie tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97091290-00	20117589	Hoofd/ hals chirurgie Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076737-00	20117651	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-10	20117769	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97096573-01	20117780	Hoofd/hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	20117807	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076877-01	20118094	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20118104	Extremiteten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97082276-01	20118120	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007477-09	20118148	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007479-11	20118170	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20118216	Extremiteten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97058036-07	20118232	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97096718-00	20118249	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97008656-04	20118320	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008654-05	20118324	Gatlagen Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075015-02	20118554	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20128927	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076774-01	20128928	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20128934	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20129037	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006187-18	20129041	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006387-09	20129131	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20129134	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97008652-07	20129157	Neuro Rug Nek Pakket	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20129162	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20129164	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20129195	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97082276-01	20119389	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97046828-05	20119399	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006391-10	20119496	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075025-02	20129782	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076887-01	20129870	Hypofyse tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077006-01	20120460	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076751-00	20120558	Arthroscopie tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97076816-01	20121098	DIEP tray Onder	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97085666-00	20121285	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97085666-00	20121290	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
9708039-17	20123235	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97065522-05	20123273	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
9708039-17	20123330	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97085666-00	20134061	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97085666-00	20134062	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076815-01	20134072	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076877-01	20134824	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077004-01	20135752	TAVI tray	Valve replacement (TAVI)	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076751-00	20137413	Arthroscopie tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97076987-01	20148014	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97007482-06	20138103	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076776-00	20138113	CVL tray	CVC/impl. chbr plac	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20138603	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	20138676	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97074987-02	20141373	TUR Tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97075024-03	20141402	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075021-04	20141407	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076751-00	20141944	Arthroscopie tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97074870-01	20142066	Basis groot tray SKZ	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075019-02	20153096	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97007482-06	20154903	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20166389	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97087964-02	20168289	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076751-00	20168296	Arthroscopie tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97100613-02	20168875	TAAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076877-01	20179951	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006186-08	20170192	Aneurysma pakket	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076774-01	20171748	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97077005-01	20172309	Thoracotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97079569-00	20173293	Plastisch DIEP tray	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076885-01	20183605	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20183695	Pediatric tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076751-00	20183802	Arthroscopie tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97076901-03	20183911	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97074877-03	20184831	Laparoscopie Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97061798-04	20185173	FESS navigatie tray	MIS Nose/Sinus	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006290-08	20185369	Carotis Pakket	Carotis	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97064671-05	20195874	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20195905	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20195906	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97077001-01	20196938	Stent tray	Endovascular Others	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076987-01	20198956	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20199736	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076909-01	20209955	Rad. Vulvectomie tray	Abdom. hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97077008-00	20190174	Verticaal isolatielaken tray	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20202562	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97087964-02	20202587	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076706-01	20202600	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20202603	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076893-02	20202633	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076706-01	20202655	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20202660	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076893-02	20202722	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97038615-07	20203252	Mammareductie/Abdominoplastiek tray	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-17	20203268	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97091291-00	20203282	Mamma Tray chirurgie	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97062064-07	20203283	Laparoscopie abdo-perineaal tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97036257-10	20203319	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97036257-10	20203388	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97019003-10	20203412	Sectio Tray	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97096718-00	20203430	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97006391-10	20203437	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97007234-08	20203441	Perifeer vaat pakket	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97077234-02	20203477	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97065522-05	20203509	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97074606-01	20203721	Ortho Hand/Voet tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97097266-04	20203761	Procedure tray basis neuro	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97092416-02	20203829	KNO Bijholte Tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	20203897	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97087965-02	20203901	Laparoscopische uterus tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB

97006287-12	20203904	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97099157-01	20203933	Kaak tray SKZ	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97071539-07	20204293	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97071539-07	20204432	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97096854-02	20204646	Neuro rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97068089-06	20204828	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97068089-06	20204829	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076879-01	20215300	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97097363-03	20215392	TUR TRAY	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97096573-01	20215394	Hoofd/hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20205404	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20205499	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076886-01	20205503	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076933-02	20205794	Rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20205907	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076735-00	20205926	Aangezichtstray	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076902-01	20206034	Pediatric Laparoscopie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20206096	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97076774-01	20206155	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	20216772	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077004-01	20216775	TAVI tray	Valve replacement (TAVI)	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076877-01	20216834	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	20216839	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-05	20216869	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97064671-05	20216870	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97071539-07	20207013	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97077234-02	20217076	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97097266-04	20218272	Procedure tray basis neuro	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97006187-18	20219177	Knepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-18	20219198	Knepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III

97096854-02	20219298	Neuro rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97011376-10	20219509	VKB Pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97074259-01	20210096	Ortho Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20210137	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20210147	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20210208	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97007477-09	20210354	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97006391-10	20210930	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	20210954	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	20211010	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97058036-07	20211058	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076987-01	20221862	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97077006-01	20222140	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97074259-01	20222215	Ortho Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97103564-00	20222734	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97076754-01	20222766	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077007-01	20222970	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97071539-07	20212978	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA

97058036-07	20223006	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97007482-06	20223057	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97091738-01	20224336	Angio tray	Angioplasty (Coronary - PTCA, Carotid or Peripheral) - dilatation with or w/o stent)	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97006184-10	20224395	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-10	20224398	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20225053	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97008656-04	20225068	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97096718-00	20225123	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97008656-04	20225190	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077235-02	20225231	KNO oor tray	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97092126-01	20225291	Hoofd/Hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97079569-00	20225354	Plastisch DIEP tray	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97038473-06	20225830	ARTHROSCOPIE TRAY	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97041288-06	20225915	CTS tray POLI POK	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97069959-02	20236679	Bekkenbodem Tray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97076987-01	20236681	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-17	20236756	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20236808	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97006391-10	20236869	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97074877-03	20237296	Laparoscopie Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97075021-04	20237482	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	20237483	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074987-02	20237650	TUR Tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076737-00	20237772	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006182-13	20237844	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97007477-09	20237916	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97019003-10	20237936	Sectio Tray	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97008654-05	20237937	Gatlaken Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076879-01	20238090	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97075019-02	20238093	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97076932-01	20238113	Re-Sternotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076892-01	20238117	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20238118	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074882-04	20238136	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97075021-04	20238814	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97061798-04	20239309	FESS navigatie tray	MIS Nose/Sinus	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97075021-04	20239322	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97007488-06	20239419	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	20239539	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	20239582	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007226-07	20239688	KNO oor pakket	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	20230037	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97007230-10	20230069	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97085666-00	20230316	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074606-01	20230376	Ortho Hand/Voet tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97085666-00	20230411	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97077005-01	20230894	Thoracotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97035304-12	20230977	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97091290-00	20241222	Hoofd/ hals chirurgie Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075015-02	20241549	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075042-02	20241800	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97074553-01	20242007	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97061797-06	20242030	Neus tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074989-02	20242118	Gynaecologie Laparoscopie Tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97076885-01	20243058	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	20243074	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97074875-01	20243076	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97008656-04	20243158	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075024-03	20243288	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074820-03	20243291	PBC Algemeen Erasmus MC Chir.Dagbehandeling	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97077006-01	20243320	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076781-01	20243324	Diagn. Laparoscopie tray Gynaecologie	MIS Extruterine (Myomectomy/Salpingo- oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97092697-00	20243365	Oor tray	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97076901-03	20244081	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97050152-10	20244123	Knie tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97078860-00	20245198	Gyn Hysteroscopie tray	MIS Hysterectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97006391-10	20245242	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97008652-07	20245674	Neuro Rug Nek Pakket	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008652-07	20245766	Neuro Rug Nek Pakket	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97062064-07	20246034	Laparoscopie abdo-perineaal tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97058036-07	20246183	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074635-02	20246303	Extractie Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006387-09	20256654	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20256689	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20256692	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074877-03	20256698	Laparoscopie Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97008039-17	20257027	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	20257063	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97007479-11	20257116	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20257129	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007233-05	20257134	Vaginaal Pakket	Vaginal hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20257352	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076755-02	20257381	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97074987-02	20257437	TUR Tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97075016-02	20258073	VKB Tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	III
97076895-01	20258083	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075021-04	20258106	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075021-04	20258141	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076895-01	20258145	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075021-04	20258156	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97074881-03	20258311	Hand- ok Tray EMC Sophia	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076776-00	20263409	CVL tray	CVC/impl. chbr plac	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	20277331	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20277338	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20277339	Extremiteten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075042-02	20277340	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97046828-05	20277347	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20277353	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076892-01	20277357	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97085666-00	20277366	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075019-02	20277369	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97076987-01	20277370	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076879-01	20277376	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97007477-09	20277383	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97091290-00	20277387	Hoofd/ hals chirurgie Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076877-01	20277395	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	20277426	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97007230-10	20277454	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA

97096718-00	20277459	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97007482-06	20277463	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	20277481	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97011376-10	20277482	VKB Pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97006391-10	20277485	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97008656-04	20277494	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97085666-00	20277509	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076882-00	20277515	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97091612-02	20277531	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074882-04	20277583	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97006290-08	20277585	Carotis Pakket	Carotis	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97074874-04	20277608	Splittaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97008656-04	20277612	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97096853-01	20277615	Neuro nervus ulnaris	Peripheral nerves	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97076888-01	20277626	Incontinentietray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	20277653	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97097362-05	20277700	Uretrplastiek perineaal	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIB
97092697-00	20277706	Oor tray	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97075024-03	20277720	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074872-02	20277776	Oor Tray SKZ	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97104554-01	20277784	Maxilla Tray	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97075012-04	20277805	Vrije Lap Tray Erasmus MC	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97076815-01	20277809	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97068089-06	20277900	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97075025-03	20277934	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97074881-03	20277939	Hand- ok Tray EMC Sophia	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97099157-01	20277957	Kaak tray SKZ	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076893-03	20277959	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97107024-01	20277969	Bilateraal Splittaken tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97068089-06	20277978	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97046828-05	20278742	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007477-09	20278826	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97058036-07	20278840	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074990-04	20278844	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97075019-02	20278864	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97007233-05	20278872	Vaginaal Pakket	Vaginal hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97006391-10	20278879	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97074553-01	20278960	Basis Kaak Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97091612-02	20278968	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074990-04	20278992	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97076902-01	20279003	Pediatric Laparoscopy tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074882-04	20279028	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III

97074874-04	20279038	Splitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076706-01	20279142	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20279161	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97092126-01	20279232	Hoofd/Hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075015-03	20279270	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076900-02	20279328	Oncologie tray KNO	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97075015-03	20279371	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97092354-02	20279494	Abdominaal Gynaecologie	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97006182-13	20270236	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97082859-00	20270257	Kinder tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97082619-01	20270341	CAD laparoscopie tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97075021-04	20280589	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075021-04	20280590	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074875-01	20280646	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074875-01	20280647	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97006287-12	20281387	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075025-03	20281388	Heup/Knie Laminaar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006287-12	20281393	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97075014-01	20281399	Vorbereidings Tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075025-03	20281923	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97075015-03	20281948	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	20282808	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075018-01	20282810	Algemene Hand Tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97071539-07	20282859	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA

97074990-04	20282864	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074882-04	20282880	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97077234-02	20282897	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074882-04	20282921	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97071539-07	20282965	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97074987-02	20283023	TUR Tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97006187-20	20283059	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-20	20283063	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076892-01	20283542	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20283544	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20283591	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076877-01	20283598	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076879-01	20283610	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97091291-00	20283613	Mamma Tray chirurgie	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074874-04	20283739	Splitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97046150-08	20283804	Vaginale uterus set	Vaginal hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97007465-06	20293815	Tulband Pakket	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046150-08	20283895	Vaginale uterus set	Vaginal hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97038655-04	20293897	Handset Plastische tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97038655-04	20283904	Handset Plastische tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97038655-04	20283945	Handset Plastische tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007477-09	20294143	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97074606-01	20294191	Ortho Hand/Voet tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97062077-04	20294303	Laparoscopie algemeen tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97064671-05	20294339	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97074820-03	20294368	PBC Algemeen Erasmus MC Chir.Dagbehandeling	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97007226-07	20294403	KNO oor pakket	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20294457	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077003-00	20294685	Subtray Otologie	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074874-04	20294831	Splittaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075019-02	20294855	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075024-03	20294946	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075017-04	20294980	Wervel Tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074990-04	20295831	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97075013-02	20296037	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97006387-09	20296084	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97074989-02	20296087	Gynaecologie Laparoscopie Tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075015-03	20296131	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20296506	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20296514	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	20296530	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97069959-02	20296546	Bekkenbodem Tray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III

97041288-06	20296591	CTS tray POLI POK	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	20296602	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97071539-07	20296625	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076987-01	20296627	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	20296635	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97007230-10	20296646	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA

97071539-07	20296679	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97075013-02	20296873	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076754-01	20296874	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	20296924	Extremiteiten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-05	20296987	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20296990	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97064671-05	20296991	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97075025-03	20297030	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97007488-06	20297035	Extremiteiten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	20297038	Extremiteiten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97087964-02	20307339	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97036244-09	20307425	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97036244-09	20307430	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97092126-01	20307437	Hoofd/Hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20307476	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20307477	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20307479	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20307480	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97007231-08	20307724	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97085666-00	20297959	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076883-01	20298004	Fem-Pop tray	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97006391-10	20298178	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006182-13	20308591	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075018-01	20308608	Algemene Hand Tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076879-01	20308615	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97074990-04	20308658	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97007482-06	20308663	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076776-00	20308667	CVL tray	CVC/impl. chbr plac	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	20308672	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075025-03	20309167	Heup/Knie Laminaar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III

97075042-02	20309836	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97074874-04	20309845	Splitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	20309903	EMC Extremitäten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075013-02	20309916	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075016-02	20309924	VKB Tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	III
97075013-02	20309988	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97075044-04	20300011	Craniotomie Tray EMC Sophia	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97085491-00	20300028	Arthroscopie tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075014-01	20300100	Vorbereidings Tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	20300544	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	20300628	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97100613-03	20300796	TAAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA

97008039-17	20300831	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20300879	Extremitäten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-17	20301028	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074882-04	20301077	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074259-01	20301821	Ortho Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074259-01	20301839	Ortho Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076706-01	20301845	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20301845	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075015-03	20301859	Craniotomie Tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20301915	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20301917	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076756-01	20313667	CAD tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076882-00	20314695	Extremitäten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97065522-05	20314769	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076987-01	20315822	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20315831	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077001-01	20316044	Stent tray	Endovascular Others	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97046828-05	20317206	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97035304-12	20317234	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97035304-12	20317396	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97052281-05	20328639	DIEP Pakket	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-10	20328758	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-05	20328759	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20328858	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97006184-10	20328859	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076816-01	20320795	DIEP tray Onder	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076737-00	20333068	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20333126	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20333136	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20323510	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076933-02	20323626	Rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	20333937	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077007-01	20336991	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076892-01	20337013	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20337024	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076885-01	20337036	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076877-01	20338087	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20348098	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076888-01	20338292	Incontinentietray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076706-01	20348431	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076933-02	20340428	Rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20341419	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076895-01	20353310	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20353384	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97006187-20	20353444	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97006187-20	20353505	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076886-01	20354085	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20354143	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97046828-05	20356131	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20356207	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97087965-02	20356724	Laparoscopische uterus tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB
97006287-12	20356769	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	20356770	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20357633	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97077006-01	20357741	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076706-01	20357786	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20357786	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-05	20358088	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20358090	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076883-01	20358341	Fem-Pop tray	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97076885-01	20358373	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20358487	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20358615	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20368726	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076893-03	20359074	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076893-03	20359119	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076774-01	20369559	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20360005	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	20360219	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	20360312	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20360658	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20360660	Extremiteten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076882-00	20360820	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20363255	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20364445	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077008-00	20375237	Verticaal isolatielaken tray	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-01	20365726	Drain Volwassene en Kind tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	20377848	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076706-01	20378699	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20378730	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-10	20378924	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-10	20379086	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20370499	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20370526	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97035304-12	20370544	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97035304-12	20370634	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077001-01	20370763	Stent tray	Endovascular Others	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076885-01	20371894	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	20372035	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077005-01	20383306	Thoracotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076904-02	20383347	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20383354	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20383586	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-17	20384726	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-17	20384846	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077006-01	20385025	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB

97076815-01	20385300	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076755-02	20386371	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97068089-06	20386691	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97068089-06	20386761	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076754-01	20387644	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97087964-02	20387700	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB

97006187-20	20389119	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97041535-11	20399967	AAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97065522-05	20391243	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076987-01	20408027	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20408116	Extremiteten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20408183	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076932-01	20408231	Re-Sternotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076882-00	20408286	Extremitäten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97041535-11	20408456	AAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076781-01	20408456	Diagn. Laparoscopie tray Gynaecologie	MIS Extraterine (Myomectomy/Salpingo- oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97006187-20	20408597	Knip pakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97064671-05	20408659	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III

97064671-05	20408718	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076776-00	20409034	CVL tray	CVC/impl. chbr plac	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97041535-11	20409101	AAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076933-02	20409125	Rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	20409412	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20409582	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076816-01	20409761	DIEP tray Onder	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076754-01	20419769	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20409808	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97035304-12	20409809	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97035304-12	20409898	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97052281-05	20409975	DIEP Pakket	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006187-20	20400037	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97064671-05	20400103	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97068089-06	20400113	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20400144	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97068089-06	20400196	Kleppen Tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076888-01	20400841	Incontinentietray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB

97076902-01	20400879	Pediatric Laparoscopie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077006-01	20410969	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97087965-02	20411343	Laparoscopische uterus tray	MIS Extruterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB
97046828-05	20412126	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20412222	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	20412601	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006287-12	20412625	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20413475	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076895-01	20413824	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20413857	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20413858	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20413860	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076737-00	20413861	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20414574	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97006186-08	20414691	Aneurysma pakket	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076900-02	20414761	Oncologie tray KNO	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20415145	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-02	20415279	Drain Volwassene en Kind tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076892-01	20417166	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076893-03	20428881	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076893-03	20428882	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076892-01	20429160	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-17	20420668	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008039-17	20420712	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97065522-05	20420744	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076774-01	20422048	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20422238	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20422298	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20422344	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076815-01	20422427	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97006187-20	20422550	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97100613-03	20422598	TAAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97006187-20	20422640	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076987-01	20423888	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	20434435	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97041535-11	20435858	AAA tray	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA

97007231-08	20435993	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007231-08	20436115	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20437849	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20437871	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076887-01	20438272	Hypofyse tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20439657	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076706-01	20439677	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20439684	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006187-20	20439792	Knienpakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97064671-05	20439871	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20439873	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97087965-02	20449880	Laparoscopische uterus tray	MIS Extrauterine (Myomectomy/Salpingo- oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB

97006287-12	20439891	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077007-01	20430493	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97007482-06	20430623	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20440905	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20440910	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006186-08	20441000	Aneurysma pakket	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA

97076932-01	20441809	Re-Sternotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076754-01	20442004	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	20442704	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20444080	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-02	20444758	Drain Volwassene en Kind tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20446293	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97084957-02	20446314	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97052281-05	20456684	DIEP Pakket	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076774-01	20457647	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20457656	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20457846	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076901-03	20458070	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97076895-01	20458071	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20458074	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	20458555	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076883-01	20459062	Fem-Pop tray	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97077001-01	20452052	Stent tray	Endovascular Others	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97008039-17	20452902	Heuptray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97077006-01	20452924	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076776-00	20453028	CVL tray	CVC/impl. chbr plac	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076933-02	20453203	Rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97064671-05	20463945	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97064671-05	20463946	CABG Tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97076987-01	20465750	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076774-01	20465758	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20466138	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20466215	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-11	20467472	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-11	20467502	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97046828-05	20467604	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076816-01	20468096	DIEP tray Onder	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076781-01	20468150	Diagn. Laparoscopie tray Gynaecologie	MIS Extraterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076894-01	20468153	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97046828-05	20468221	ASI Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076904-02	20468351	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076879-01	20469042	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97035304-12	20479071	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97041288-06	20469104	CTS tray POLI POK	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007482-06	20469234	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077234-02	20469285	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97007230-10	20479297	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97007230-10	20479340	Arthroscopie pakket	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA

97075025-03	20479421	Heup/Knie Laminar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076888-01	20469460	Incontinentietray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97077005-01	20479504	Thoracotomie tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97077008-00	20479630	Verticaal isolatielaken tray	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97052281-05	20479725	DIEP Pakket	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97082276-01	20470621	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97006182-13	20470678	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97065522-05	20470687	LAPAROSCOPISCHE LAR/SIGMOID TRAY	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076706-01	20470895	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076706-01	20470929	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20470973	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076737-00	20471094	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006287-12	20471291	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	20471513	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075042-02	20472229	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97074875-01	20472249	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97091612-02	20472314	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97058036-07	20472318	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97074882-04	20472351	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97096718-01	20472389	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97074882-04	20472401	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97085666-00	20472441	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97038473-06	20474125	ARTHROSCOPIE TRAY	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97058036-07	20474152	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97074875-01	20474165	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074606-01	20474189	Ortho Hand/Voet tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007479-11	20474279	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20474870	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	20475363	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076879-01	20475373	Endo Uro tray	Kidney procedures	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97076886-01	20486120	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20486219	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97006187-20	20476383	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076894-01	20486534	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97074875-01	20486546	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97006187-20	20476835	Kniepakket	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III

97076737-00	20487197	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20487531	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97075042-02	20487533	PM/ICD Pakket	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97075019-02	20487539	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075018-01	20487551	Algemene Hand Tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074989-02	20487557	Gynaecologie Laparoscopie Tray	MIS Extrauterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97076987-01	20487569	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074990-04	20487581	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97071539-07	20487653	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97074990-04	20487691	Neuro tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97071539-07	20487705	Pacemaker Tray	Pacemaker or ICD implantation	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97074882-04	20487767	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III

97075044-05	20488406	Craniotomie Tray EMC Sophia	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077007-01	20488540	Vats/ Nuss-Bar tray	Cardiovascular Other	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97058036-07	20489519	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074606-01	20489543	Ortho Hand/Voet tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074875-01	20489557	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074882-04	20489664	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III

97077006-01	20489771	Vaginaal tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97074989-02	20492274	Gynaecologie Laparoscopie Tray	MIS Extruterine (Myomectomy/Salpingo-oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97075019-02	20492304	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97087964-02	20492530	Debulking tray	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076885-01	20492607	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076894-01	20492713	Laparoscopische urologie tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97052281-05	20492929	DIEP Pakket	Breast procedures (esthetic)	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97069880-02	20492998	Bot in gnatho tray SKZ	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076901-03	20493904	Oncologie tray Urologie	Urology Others	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97103564-00	20494053	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97006182-13	20504982	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97082276-01	20494998	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97075024-03	20494999	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076774-01	20495008	Craniotomie tray	Craniotomy	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97008656-04	20505058	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006391-10	20505059	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97058036-07	20505080	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97082276-01	20495113	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA

97104554-02	20495192	Maxilla Tray	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97074874-04	20495267	Spitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076754-01	20495414	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20495424	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97097363-03	20495547	TUR TRAY	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97076815-01	20495639	DIEP tray Boven	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97007477-09	20495705	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97075014-01	20495817	Vorbereidings Tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97078860-00	20495947	Gyn Hysteroscopie tray	MIS Hysterectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97076886-01	20506263	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075013-02	20496286	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075013-02	20506318	Basis Chirurgie Klein Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97076909-01	20506640	Rad. Vulvectomy tray	Abdom. hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97091291-00	20506705	Mamma Tray chirurgie	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97061798-04	20497003	FESS navigatie tray	MIS Nose/Sinus	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076987-01	20507618	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97084957-02	20507716	Da Vinci robot tray Chirurgie	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97085666-00	20508443	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97058036-07	20508515	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97082859-02	20508539	Kinder tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074872-02	20508956	Oor Tray SKZ	Ear	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97006387-09	20508986	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20509079	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006184-11	20509086	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006184-11	20509092	Heuppakket	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97096854-02	20509781	Neuro rug tray	Laminectomy/Spinal	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076885-01	20500109	Hand tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076882-00	20500118	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076877-02	20500201	Drain Volwassene en Kind tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036257-11	20500232	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA

97076882-00	20500261	Extremitäten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97096718-01	20500264	Da Vinci robot tray	MIS Prostatectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97036257-11	20500351	Laparoscopie tray	Laparoscopy others	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97036244-09	20501578	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20501655	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97036244-09	20501659	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97036244-09	20501713	Klein chirurgische tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006387-09	20501755	Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97077234-02	20502241	KNO Septum tray	Nose	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97007234-08	20502346	Perifeer vaat pakket	Femoral	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97075025-03	20502360	Heup/Knie Laminaar Flow Tray	Hip replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076904-02	20512425	Pediatrie tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97058036-07	20513007	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97006182-13	20513013	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075018-01	20503033	Algemene Hand Tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075024-03	20503058	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97007479-11	20513062	Schouder Pakket	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97103564-00	20513077	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III

97074882-04	20503122	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97074874-04	20503126	Spitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074882-04	20503163	Sectio tray Erasmus MC Sophia	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III
97076706-01	20513183	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075021-04	20513242	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076706-01	20513252	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97076706-01	20513253	Algemene Basis tray Heelkunde	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075021-04	20513288	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076737-00	20513338	Algemene Basis tray KNO/MKA	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20503343	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97007231-08	20513414	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074870-01	20513448	Basis groot tray SKZ	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97074877-03	20513510	Laparoscopie Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97007231-08	20513568	Abdo-Perineaal Pakket	Abdo/perineal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97041288-06	20513641	CTS tray POLI POK	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97035304-12	20513647	Totale knie procedure tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074874-04	20514421	Splitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075021-04	20514521	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97076895-01	20514555	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97076895-01	20514563	Mamma tray	Breast procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97075019-02	20515101	Grote Chirurgie Laparoscopie Tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97103564-00	20515107	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97007477-09	20515112	TUR pakket	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97075011-03	20515126	Vaginaal Tray	Vaginal minor procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA

97008656-04	20515230	Hand Pakket	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006391-10	20515236	Klein Chir en Gyn pakket	Gynaecology Others	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97006186-08	20515511	Aneurysma pakket	Aortic bypass	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIA
97076887-01	20515529	Hypofyse tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97107024-01	20515592	Bilateraal Splittaken tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97099157-01	20515692	Kaak tray SKZ	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA

97006287-12	20515731	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006287-12	20515734	Groot Universeel Pakket	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97038659-08	20515999	VKB tray	MIS Knee arthroscopy	P44093	Arthroscopic surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical and endoscopic instruments, dressings, and additional materials intended to be used specifically to perform an arthroscopic surgical procedure. It does not contain any pharmaceuticals. This is a single-use device.	IIA
97085666-00	20516156	U-Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialties typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97091738-01	20516535	Angio tray	Angioplasty (Coronary - PTCA, Carotid or Peripheral) - dilatation with or w/o stent)	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	III
97069959-04	20517306	Bekkenbodem Tray	Incontinence care procedures	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	III

97076987-01	20517321	Schouder Heup tray Orthopedie	Shoulder replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074635-04	20517328	Extractie Tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97102963-01	20517387	Tonsillectomie Tray	Tonsillectomy/Adenoidectomy	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97104554-02	20517408	Maxilla Tray	ENT Others	P44069	ENT surgical procedure kit, non-medicated, single-use	A collection of various sterile ear/nose/throat (ENT) surgical instruments, dressings and the necessary materials used to perform an ENT surgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III
97007482-06	20517472	Extensie Pakket	DHS - Gamma nailing	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97090390-01	20517549	Laparoscopische Urologie Tray	MIS Nephrectomy	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIB

97007488-06	20528698	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97007488-06	20528730	Extremiteten Pakket	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074875-01	20518791	Basis klein tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075024-03	20518794	EMC Extremiteten Tray	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIB
97103564-00	20528839	Laparoscopie 3D tray	MIS Anti-reflux	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	III
97058036-07	20528857	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB

97074606-01	20528861	Ortho Hand/Voet tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97006182-13	20528985	Sectio pakket	C Section	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIA
97074874-04	20520116	Splitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076882-00	20520147	Extremiteiten tray Heelkunde	Foot/Ankle/Low limb	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076892-01	20521407	Kinderurologie basis tray	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97091612-02	20521429	Neuro rug tray	Neurosurgery Others	P44063	Neurosurgical procedure kit, non-medicated, single-use	A collection of various sterile neurosurgical instruments, dressings and the necessary materials used to perform a neurosurgical procedure. It does not contain pharmaceuticals. This is a single-use device.	III

97076909-01	20521571	Rad. Vulvectomy tray	Abdom. hysterectomy	P60644	Obstetrical/gynaecological surgical procedure kit, non-medicated, single-use	A collection of various sterile surgical instruments, dressings/drapes, and materials intended to be used during obstetrical/gynaecological surgery. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of obstetrical/gynaecological surgical procedure, but are often supplied together to perform a specific type of obstetrical/gynaecological procedure; it does not contain pharmaceuticals. This is a single-use device.	IIB
97076735-00	20521695	Aangezichtstray	Plastic surgery Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076781-01	20521722	Diagn. Laparoscopie tray Gynaecologie	MIS Extrauterine (Myomectomy/Salpingo- oophorectomy)	P44096	Laparoscopic surgical procedure kit, non-medicated, single-use	Laparoscopic surgical procedure kit, non-medicated, single-use A collection of various sterile laparoscopic and surgical instruments (e.g., access cannulae, forceps), dressings, and additional materials intended to be used specifically to perform a laparoscopic surgical procedure. It does not contain pharmaceuticals. It may be referred to as a procedure tray. This is a single-use device.	IIA
97058036-07	20521735	Universeel tray groot	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97074987-02	20521905	TUR Tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA
97038655-04	20521931	Handset Plastische tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97038655-04	20521987	Handset Plastische tray	Hand/Arm	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97074874-04	20522511	Spitlaken Tray SKZ	Paediatric surgery	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97075021-04	20522542	Basis Groot Chirurgische Tray	General/abdo Others	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIB
97076755-02	20522609	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97092126-01	20522627	Hoofd/Hals tray	Maxillo-facial	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97075014-01	20522728	Voorbereidings Tray	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA

97074869-04	20522729	Rug Tray 1 SKZ	Orthopaedic Others	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	III
97076754-01	20533311	Buik Beensteunen tray	Abdominal invasive	P33961	General surgical procedure kit, non-medicated, single-use	A collection of various surgical instruments, sometimes with dressings and/or other materials, intended to be used to: 1) assist a range of surgical procedures across multiple clinical specialities typically as a general-purpose tray; or 2) assist with a plastic surgery procedure. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076755-02	20533324	CABG tray	CABG	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97076886-01	20533345	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076886-01	20533353	Heup en Knieprothese tray	Knee replacement	P44052	Orthopaedic surgical procedure kit, non-medicated, single-use	A collection of various sterile orthopaedic surgical instruments, dressings, and other materials intended to be used to perform an orthopaedic surgical procedure, however the kit is not dedicated to orthopaedic implantation. It does not contain pharmaceuticals. This is a single-use device.	IIA
97076893-03	20533407	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB

97076893-03	20533477	Klep tray	Valve open heart surgery	P61371	Cardiothoracic surgical procedure kit, non-medicated, single-use	A collection of various sterile cardiothoracic surgical instruments, dressings, and materials intended to be used to assist a range of cardiothoracic surgical procedures; it does not contain pharmaceuticals. The kit is intended only for use during cardiothoracic surgery but is not dedicated to a specific type of cardiothoracic surgical procedure. This is a single-use device.	IIB
97082276-01	20533668	Uro transurethrale tray	TUR	P60755	Urological surgical procedure kit, non-medicated, single-use	A collection of various sterile urological surgical instruments and the necessary materials intended to be used to perform a urological surgical procedure. The kit may be referred to as a procedure pack and is typically comprised of devices that individually are not dedicated to a specific type of urological surgical procedure, but are often supplied together to perform a specific type of urological procedure; it does not contain any pharmaceuticals. This is a single-use device.	IIA